

ABI PRISM[®] GeneScan Analysis Software

Version 3.7 for the Windows NT[®] Platform

User Guide

ABI PRISM[®]
GeneScan[®] Analysis
Software

Version 3.7 for the Windows
NT[®] Platform

© Copyright 2001, Applied Biosystems. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

Notice to Purchaser: Limited License

For Software License and Warranty Information, please see Appendix G, "License and Warranty," of this manual.

Notice to Purchaser: License Disclaimer

Purchase of this software product alone does not imply any license under any process, instrument or other apparatus, system, composition, reagent or kit rights under patent claims owned or otherwise controlled by Applied Biosystems, either expressly, implied, or by estoppel.

ABI PRISM and its design, Applied Biosystems, BioLIMS, GeneScan, and Genotyper are registered trademarks and ABI, POP-5 and POP-6 are trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries.

Windows NT and Microsoft are registered trademarks of the Microsoft Corporation.

Oracle is a registered trademark and Oracle7 is a trademark of Oracle Corporation.

Sybase is a registered trademark and Sybase SQL Server and SyBooks are trademarks of Sybase, Inc.

Apple, Macintosh, and Power Macintosh are registered trademarks of Apple Computer, Inc.

All other trademarks are the sole property of their respective owners.

Applied Biosystems Corporation is committed to providing the world's leading technology and information for life scientists. Applied Biosystems Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Contents

1 GeneScan Analysis Software Overview

Overview	1-1
Introduction.....	1-1
In This Chapter	1-1
About the GeneScan Analysis Software	1-2
What Does the Software Do	1-2
Instruments	1-2
GeneScan Software on the 310 Instrument	1-3
Flowchart	1-3
GeneScan Analysis on the 377 Instrument	1-4
Flowchart	1-4
GeneScan Analysis on the 3100 Instrument	1-5
Flowchart	1-5
GeneScan Analysis on the 3700 Instrument	1-6
Flowchart	1-6
Sequence Collector Database Option	1-7
Related Manuals	1-7
Registering the Software	1-8
License and Warranty	1-8
Registering Your Software	1-8
Hardware and Software Requirements	1-9
Introduction.....	1-9
Computers Connected to Applied Biosystems Instruments	1-9
System Requirements	1-10
Minimum System Recommendations.....	1-10
Hard Drive Partitions	1-11

Installing the GeneScan Analysis Software	1–12
Preparation	1–12
Procedure	1–12
Starting the GeneScan Analysis Software for the First Time	1–13
Procedure	1–13
Removing the Software.	1–15
If Files or Folders Have Been Moved	1–15
Procedure	1–15

2 *Creating a Project*

Overview	2–1
In This Chapter.	2–1
Creating Projects and Auto-Analysis on 310 Instruments	2–2
Process Using the 310 Instrument	2–2
Creating Projects and Auto-Analysis on 377 Instruments	2–3
Process Using the 377 Instrument	2–3
Setting Up for Automatic Analysis on 3100 and 3700 Instruments	2–4
Introduction	2–4
Procedure	2–4
What Is a Project	2–7
Why Create a Project	2–7
Where to Store Projects	2–7
Working with Project Files	2–8
Opening an Existing Project	2–8
Creating a New Project	2–9
Unlocking Sample Files.	2–12
Removing Samples from a Project.	2–13
Finding Missing Sample Files	2–14
When Are Files Considered Lost	2–14
When an Alert Appears	2–14
Searching for Missing Sample Files	2–14
Re-Establishing Links	2–15
Procedure	2–16

3 Analyzing Project Files

Overview	3-1
In This Chapter	3-1
Analyzing Project Files: About the Analysis Control Window	3-2
Introduction	3-2
Analysis Control Window Example	3-2
Analysis Control Window Callouts Described	3-3
Customizing the Display	3-4
Using the Analysis Control Window	3-5
Analyzing Sample Files: Using the Analysis Control Window	3-6
Introduction	3-6
Accessing Sample Files	3-6
Analyzing Sample Files	3-7
Specifying the Format for Printed Results	3-8
Displaying Size Standards and Analysis Parameters	3-10
Displaying Sample and Dye Information	3-11
Setting Sample File Sort Order	3-13
Setting Dye Indicator Preferences	3-14
Defining Folder Locations	3-16
Introduction	3-16
Storing Matrix Files	3-16

4 Analyzing Sample Files

Overview	4-1
In This Chapter	4-1
About Sample Files	4-2
Sample File Generation on 377	4-2
What Files Contain	4-2
Sample Files	4-2
How GeneScan Analyzes Sample Files	4-3
Converting Macintosh Computer Sample Files	4-3
About Converting Files	4-3

Converting Sample Files for Use on a Different Platform	4-4
Introduction	4-4
Converting Macintosh Files to Windows NT-Based Computer Files	4-4
Converting Windows NT-Based Computer Files to Macintosh Files	4-5
Opening Sample Files.	4-7
Introduction	4-7
Procedure	4-7
About the Sample File Window	4-8
What It Displays.	4-8
Five Views	4-8
Sample Results View.	4-9
What It Displays.	4-9
Displaying the View	4-9
Example of Sample Results View	4-10
Description of Columns	4-10
Differences from the Results Display	4-11
Sample Info View	4-11
What It Displays.	4-11
Displaying the View	4-12
Example of Sample Info View	4-13
Description of Information.	4-14
Size Curve View	4-20
What It Displays.	4-20
Displaying the View	4-20
Example of Size Curve View	4-21
Curves Described	4-21
Raw Data View	4-22
What It Displays.	4-22
Displaying the View	4-22
Example of Raw Data View	4-23
What to Evaluate	4-23
EPT Data View	4-24
What It Displays.	4-24
Displaying the View	4-24

EPT Data View Example.	4-25
Colored Lines Described.	4-25
Analyzing a Sample File	4-26
Introduction.	4-26
Procedure	4-26
Installing a New Matrix File.	4-27

5 *Working with Analysis Parameters*

Overview	5-1
In This Chapter	5-1
About the Analysis Parameters	5-2
What They Are	5-2
Why They Are Necessary	5-2
When to Specify a Parameter	5-2
Analysis Parameter Files Provided	5-3
When Analysis Parameters Are Used	5-3
Sizecaller Algorithm Flowchart.	5-4
Flowchart	5-4
Setting Analysis Parameters	5-5
Default Settings.	5-5
Displaying Analysis Parameters	5-6
Analysis Parameters	5-7
Analysis Range Options	5-7
Data Processing Options.	5-7
Peak Detection Options.	5-8
Sizecall Range Options.	5-10
Sizecalling Method Options	5-10
Baselining Option	5-12
Auto-Analysis Only Option	5-12
Using Analysis Parameter Files.	5-13
In This Section	5-13
Assigning the Same Analysis Parameters to All Files	5-13
Assigning Different Parameters to Single Samples	5-15
Displaying Default Parameters	5-16

Creating Custom Analysis Parameter Files	5–16
Changing an Existing Analysis Parameters File	5–17
Deleting Custom Analysis Parameters	5–17

6 *Making a Matrix File*

Overview	6–1
In This Chapter.	6–1
About Matrix Files	6–2
Introduction	6–2
Matrix File Definition	6–2
Multicomponent Definition	6–2
Why Is a Matrix File Necessary.	6–3
When to Create a Matrix File.	6–3
Sample Files Using Matrix File	6–4
Installing a Matrix File.	6–4
When to Assign a Matrix File	6–4
Limitations to Matrix Files	6–4
When to Create a New Matrix File	6–5
Considerations Before Making a Matrix File	6–6
Where to Store Matrix Files.	6–6
Process of Creating a New Matrix File.	6–7
Process Diagram	6–7
Steps to Create a New Matrix	6–7
Loading and Running Dye Standards for the ABI PRISM 310	6–9
Introduction	6–9
Creating a GeneScan Sample Sheet	6–9
Creating a GeneScan Injection List	6–10
Starting the Data Collection Software	6–11
If the Run was Cancelled	6–11
Run Time	6–11
Assigning Matrix File	6–11
Loading and Running Dye Standards for the ABI PRISM 377	6–12
Introduction	6–12
Creating a GeneScan Sample Sheet	6–12

Loading Matrix Standards.	6-13
Running the Matrix Standards	6-14
Generating Matrix Sample Files for the ABI PRISM 377 Instrument	6-15
Who Should Use This Step.	6-15
Gel Handling on Windows	6-15
Verify Tracking	6-15
Generating Sample Files.	6-15
Choosing a Scan Range for the Matrix Calculation	6-17
Introduction.	6-17
Viewing the Raw Data	6-17
What to Look For in the Raw Data Display	6-18
Choosing a Scan Range	6-18
Eliminating Primer Peaks	6-19
Generating a New Matrix File	6-20
Procedure	6-20
Matrix File Example.	6-21
Saving and Naming the Matrix File.	6-22
Introduction.	6-22
Naming Considerations.	6-22
Saving the Matrix File	6-22
Where to Store the Matrix File	6-22
When Data Collection and Analysis Are Performed on Different Computers.	6-22
Assigning the Matrix File to Sample Files	6-23
Introduction.	6-23
Procedure	6-23
Evaluating the Matrix File.	6-25
Introduction.	6-25
Procedure	6-25
Causes for Bad Matrix Files	6-26
If an Error Message Appears.	6-26
Two Causes of Bad Matrix Files.	6-26

7 *Working with Size Standards*

Overview	7-1
In This Chapter.	7-1
About Size Standards	7-2
What Are Size Standards	7-2
Advantages of Using Size Standard	7-2
Size Standards Provided.	7-2
When to Define Size Standards	7-2
If Split Peaks Appear	7-2
Defining the Size Standard	7-3
Two Ways to Define the Size Standard	7-3
Using the New Command	7-3
Using the Analysis Control Window	7-7
Using Size Standards	7-9
In This Section	7-9
Verifying Size Calculations	7-9
Changing the Number of Peaks Detected	7-9
Editing the Size Standard Definition	7-9
Editing an Existing Size Standard	7-10
Deleting an Existing Size Standard	7-12
Analyzing Samples Using the Same Size Standard	7-13
Selecting Separate Size Standards for Samples	7-14

8 *Evaluating Analysis Results*

Overview	8-1
In This Chapter.	8-1
Process of Evaluating Analysis Results	8-2
Evaluating Analysis Results	8-2
Ways to Display Analysis Results.	8-3
Two Ways to Display Analysis Results	8-3
About the Results Display Window	8-4
Introduction	8-4
Displaying the Window	8-4

Results Control Window Callouts Described	8-5
Using the Results Control Window	8-7
Selecting Display Format	8-7
Selecting Electropherogram Panels	8-7
Selecting Samples to Display	8-8
Creating Tiled Electropherogram Displays.	8-10
Unselecting Samples.	8-11
Removing Samples	8-12
Displaying the Results	8-12
Printing the Results.	8-12
Changing How the Results Are Displayed and Printed.	8-13
Procedure	8-13
About the Sample Results View.	8-15
About the View	8-15
Updating the Results	8-16
Re-analyzing the Data.	8-16
Saving and Renaming the Results Control Format	8-17
Introduction.	8-17
Important Considerations	8-17
Saving the Display Format	8-17
Working with a Previously Saved Display	8-18
Renaming the Current Results Display.	8-19

9 Evaluating Electropherograms

Overview	9-1
In This Chapter	9-1
About Electropherogram and Tabular Data Displays	9-2
Introduction.	9-2
How the Window Is Divided.	9-2
What Tabular Data Contains.	9-2
How Electropherogram Panels Are Sized.	9-2
For More Information	9-3
Displaying Electropherogram and Tabular Data	9-4
Procedure	9-4

Example of Tabular Data and Electropherogram	9-5
Table Describing Columns	9-5
Why Some Peaks May Be Visible Only in an Electropherogram . . .	9-6
Why Some Peak Areas May Have a Negative Number	9-6
Highlighting Information	9-6
Changing the Highlight Color	9-7
Adjusting Window Size	9-7
Hiding Selected Rows of Data	9-7
Limiting the Rows to Display	9-8
Displaying Electropherogram Data	9-9
Definition	9-9
Base Pairs Versus Data Points	9-9
Procedure to Display Data	9-9
Electropherogram Example	9-10
Electropherogram Callouts Described	9-10
Working with Electropherogram Data	9-12
In This Section	9-12
Displaying X- and Y-Axis Positions	9-13
Moving the Electropherogram	9-13
Changing the Dye Color	9-14
Displaying Peak Positions	9-14
Highlighting Peaks	9-15
Using Legends to Change the Display	9-15
Scrolling the Display	9-16
Zooming In and Out	9-17
Showing Off-Scale Data	9-18
Electropherogram Displaying Off-Scale Data	9-19
Electropherogram Displaying the Flat-Topped Effect	9-19
Showing Data by Fragment Size	9-20
Changing the Horizontal Scale	9-21
Changing the Vertical Scale	9-22
Assigning Standard or Custom Colors	9-23
Defining Custom Colors in Electropherograms	9-24
Introduction	9-24
Why Change Colors in the Electropherogram	9-24

Saving the Display Format	9-24
Defining Custom Colors	9-25
Defining Individual Plot Colors	9-26
Changing the Dye Scale in Electropherograms	9-28
What the Dye Scale Defines	9-28
Increasing the Dye Scale Example	9-28
Changing the Dye Scale of an Electropherogram	9-29
Changing the Dye Scale Preferences	9-30
Process of Verifying Results	9-31
Introduction	9-31
Steps to Verify Size Calculation	9-31
Verifying Size Calculations	9-33
Introduction	9-33
Verifying for the GeneScan-350 Standard	9-33
Evaluating for Multiple Size Standards	9-35
Using the Analysis Log	9-36
What Is the Analysis Log	9-36
Displaying the Analysis Log	9-36
What to Evaluate	9-37
Removing Information from the Analysis Log	9-37
Closing the Analysis Log	9-37
Verifying Peak Detection	9-38
Introduction	9-38
Verifying Peak Detection	9-38

10 Saving, Archiving, and Copying Files

Overview	10-1
In This Section	10-1
Why Save GeneScan Files	10-2
Reasons for Saving Files	10-2
Saving GeneScan Files	10-3
Introduction	10-3
Saving Projects	10-3

Saving Sample Files	10-3
Saving Results Displays	10-3
Archiving Sample Files	10-4
When to Archive Sample Files	10-4
Procedure	10-4
Transferring Data to Other Applications	10-5
Genotyper Software	10-5
Cutting and Pasting Tabular Data	10-5
Creating a Text File	10-5

11 Printing Results

Overview	11-1
In This Chapter	11-1
About Printing	11-1
Ways You Can Print	11-1
If You Get Unexpected Results	11-1
Printing Run Results Automatically	11-2
Introduction	11-2
From the Data Collection Software	11-2
Procedure	11-2
Printing Selected Sample Files	11-3
Introduction	11-3
Setting Printing Options	11-3
Printing from the Results Control Window	11-3
Printing from the File Menu	11-4

A Creating GeneScan Analysis Modules

Overview	A-1
In This Chapter	A-1
About GeneScan Analysis Modules	A-2
Introduction	A-2
Analysis Module Format	A-2
Files Referenced by the Module	A-2

Default Settings.	A-2
Relationship Between a Module and the Files	A-4
Never Modify the Files	A-4
Creating GeneScan Analysis Modules.	A-5
Summary of the Procedure	A-5
Creating a Size Standard File	A-6
Creating an Analysis Parameter File.	A-10

B Sizecalling Methods

Overview	B-1
In This Appendix.	B-1
Least Squares Method	B-2
About This Method.	B-2
Least Squares Sizecalling Examples.	B-2
Advantages	B-3
Cubic Spline Interpolation Method	B-4
About This Method.	B-4
Possible Local Sizing Inaccuracy	B-4
Local Southern Method	B-5
About This Method.	B-5
The Equation.	B-5
How This Method Works	B-6
Global Southern Method	B-7
About This Method.	B-7
The Equations	B-7
How This Method Works	B-8

C GeneScan Size Standards

Overview	C-1
About the Size Standards	C-1
Size Standards Included	C-1
GeneScan 120 Size Standard.	C-2
About This Size Standard	C-2

Special Uses	C-2
How It Is Prepared	C-2
Fragment Lengths	C-2
Denaturing Electropherogram	C-2
Electropherogram of GeneScan 120 LIZ	C-3
GeneScan 350 All Size Standard	C-4
About This Size Standard.	C-4
How It Is Prepared	C-4
GeneScan 350 Molecular Lengths	C-4
Running Under Denaturing Conditions	C-5
Electropherogram of GeneScan 350	C-6
Double-Stranded GeneScan 500 Fragments	C-6
GeneScan 350 377 Size Standard	C-7
About This Size Standard.	C-7
GeneScan 350 377 Molecular Lengths	C-7
Electropherogram of GeneScan 350 377	C-7
GeneScan 350-250 Size Standard.	C-8
About This Size Standard.	C-8
GeneScan 350-250 Molecular Lengths	C-8
Electropherogram of GeneScan 350-250	C-8
GeneScan-400HD Size Standard	C-9
About This Size Standard.	C-9
Special Uses	C-9
How It Is Prepared	C-9
Fragment Lengths	C-9
Denaturing Electropherogram	C-9
Electropherogram of GeneScan-400HD	C-10
GeneScan 500 All Size Standard	C-11
About This Size Standard.	C-11
How It Is Prepared	C-11
GeneScan 500 Molecular Lengths	C-11
Running Under Denaturing Conditions	C-11
Electropherogram of GeneScan 500	C-12
Double-Stranded GeneScan 500 Fragments	C-12

GeneScan 500 377 Size Standard	C-13
About This Size Standard	C-13
GeneScan 500 377 Molecular Lengths	C-13
Electropherogram of GeneScan 500 377	C-13
GeneScan 500-250 Size Standard	C-14
About This Size Standard	C-14
GeneScan 500-250 Molecular Lengths	C-14
Electropherogram of GeneScan 500-250	C-15

D Troubleshooting the GeneScan Software

Introduction	D-1
In This Appendix	D-1
Troubleshooting Projects and Results	D-2
Table Description	D-2
Troubleshooting Gel Data	D-4
Table Description	D-4
Troubleshooting Genotyping Results	D-7
Table Description	D-7
GeneScan Analysis Software Error Messages	D-8
Introduction	D-8
Analysis Log Error Messages	D-8
Error Messages When Defining Size Standards	D-9

E GeneScan Analysis Software Files

Table of Files	E-1
--------------------------	-----

F Technical Support

Contacting Technical Support	F-1
To Contact Technical Support by E-mail	F-1
Hours for Telephone Technical Support	F-2
To Contact Technical Support by Telephone or Fax	F-2

To Reach Technical Support Through the Internet	F-6
To Obtain Documents on Demand	F-6

G License and Warranty

Applera Corporation Software License and Limited Product Warranty . . .	G-1
Warranty Statement	G-1
Copyright	G-1
License	G-1
Restrictions.	G-2
Limited Warranty	G-2
Limitation of Liability	G-2
Term	G-2
Miscellaneous.	G-3

Glossary

Index

GeneScan Analysis Software Overview

1

Overview

Introduction This chapter provides a general introduction to the ABI PRISM® GeneScan® Analysis Software, information about the organization of this manual, and instructions on how to get help from Applied Biosystems.

In This Chapter Topics in this chapter include the following:

Topic	See Page
About the GeneScan Analysis Software	1-2
GeneScan Software on the 310 Instrument	1-3
GeneScan Analysis on the 377 Instrument	1-4
GeneScan Analysis on the 3100 Instrument	1-5
GeneScan Analysis on the 3700 Instrument	1-6
Registering the Software	1-8
Hardware and Software Requirements	1-9
Installing the GeneScan Analysis Software	1-12
Starting the GeneScan Analysis Software for the First Time	1-13
Removing the Software	1-15

About the GeneScan Analysis Software


What Does the Software Do	The GeneScan Analysis Software performs DNA fragment analysis, which separates a mixture of DNA fragments according to their lengths, provides a profile of the separation, and estimates the lengths and sizes of the fragments.
----------------------------------	---

Instruments	GeneScan Analysis Software Version 3.7 works on the following instruments:
--------------------	--

- ◆ ABI PRISM[®] 310 Genetic Analyzer
 - ◆ ABI PRISM[®] 377 DNA Sequencer
 - ◆ ABI PRISM[®] 3100 Genetic Analyzer
 - ◆ ABI PRISM[®] 3700 DNA Analyzer
-


GeneScan Software on the 310 Instrument

Flowchart




GeneScan Analysis on the 377 Instrument

Flowchart




GeneScan Analysis on the 3100 Instrument

Flowchart



GeneScan Analysis on the 3700 Instrument

Flowchart



Sequence Collector Database Option With the Sequence Collector database system, data is collected by the Windows NT-based computer and exported to a Sequence Collector database on a networked server.

The data can later be viewed, edited, and analyzed using the GeneScan Analysis Software version 3.7 for Windows NT®, either using the same Windows NT-based computer used to collect the data, or a different Windows NT-based computer with access to the Sequence Collector database. The data can also be viewed and edited (but not analyzed) using the GeneScan Analysis Software on a Macintosh® computer.

Related Manuals GeneScan Analysis Software is part of a suite of Applied Biosystems hardware and software products.

If the information you need is not in this manual, it may be in one of the other manuals listed in the table below.

For more information about...	See...	Part Number
ABI PRISM® 3700 DNA Analyzer	<i>ABI PRISM 3700 DNA Analyzer User's Manual</i>	4325941
ABI PRISM® 3100 Genetic Analyzer	<i>ABI PRISM® 3100 Genetic Analyzer User's Manual</i>	4315834
ABI PRISM® 377 DNA Sequencer	<i>ABI PRISM® 377 DNA Analyzer User's Manual</i>	4325703
ABI PRISM® 310 Genetic Analyzer	<i>ABI PRISM® 310 Genetic Analyzer User's Manual</i>	4317588
specific GeneScan chemistry protocols, designing experiments, and preparing samples	♦ GeneScan® Reference Guide, Chemistry Reference for the ABI PRISM 310 Genetic Analyzer, or ♦ The protocols that accompany Applied Biosystems reagent kits.	4303189
attaching the new matrix to an ABI PRISM 377 gel file	<i>Gel Processor User's Manual</i>	
accessing and managing a Sequence Collector database	<i>Sequence Collector User Guide NT</i>	4319527
installation and administration of Sequence Collector	<i>Sequence Collector Installation and System Administration Guide</i>	4319526

Registering the Software

License and Warranty	Before you begin, read Appendix G, “License and Warranty.” This appendix explains your rights and responsibilities regarding the software.
-----------------------------	--

Registering Your Software	To register your copy of the GeneScan Analysis Software, complete the registration card (included in this software package) and return it to Applied Biosystems.
----------------------------------	--

Note Registering the software enables us to send you notification of software updates and any other future information that may be specific to GeneScan Analysis Software owners.

IMPORTANT Your product registration number is located on the Registration card. Be sure to record this number here before you return the Registration card.

Registration Number:

Hardware and Software Requirements

Introduction The GeneScan Analysis Software can be installed on the Windows NT®-based computer connected to your ABI PRISM® instrument or on any other Windows NT-based computer that meets the minimum requirements stated below. The software can be installed on a computer used for analysis only, and on one used for both data collection and analysis.

Computers Connected to Applied Biosystems Instruments	The following table provides important information about computers connected to Applied Biosystems instruments:	
	If...	Then...
	you received this software with a newly purchased instrument	the GeneScan Analysis Software was installed by your Applied Biosystems Customer Support Engineer as part of the installation and setup of the instrument. The system requirements for that computer are described in the instrument user's manual.
	you are updating the GeneScan Analysis Software, or if you are installing the software on a computer other than the one supplied with your instrument	ensure your computer meets the minimum requirements provided below before loading the software. IMPORTANT Your computer MUST meet the requirements listed to run the new GeneScan Analysis Software.

System Requirements

Below are the system requirements to run the GeneScan Analysis Software v. 3.7 for Windows NT® platform on your instrument or analysis computer.

IMPORTANT Applied Biosystems strongly recommends using the computer supplied with your instrument for running GeneScan Analysis software. The software was optimized and tested on systems similar to that shipped with your instrument. Running the software on systems that do not meet the following requirements may cause data loss or other significant problems.

Note These are the minimum requirements. In general, the more memory, the larger the screen size, and the more processing power you have, the better.

System Component	System Requirements	
	For use with 310 and 377 Instruments	For use with 3100 and 3700 Instruments
Model	Dell 733 GX 100 Medium Desktop	Dell WorkStation 220
Processor Speed/Bus	733 MHz/133 MHz	733 MHz/133 MHz
CD-ROM drive	Any	Any
Operating system	Windows NT v. 4. with Service Pack 5	Windows NT v. 4. with Service Pack 5
RAM	The minimum memory requirement is 128 MB, although 256 MB is recommended.	256 MB RAM

Minimum System Recommendations

System Component	Recommendations
Monitor	A 17-inch monitor or larger with 1024 x 768 resolution
Disk Space	Storage requirements depend primarily on the quantity of data to be generated and stored. It is common to store many sample files on the analysis computer. Each sample file is approximately 150–250 KB.
Printer	A PostScript-compatible color printer

**Hard Drive
Partitions**

Use the following information to determine on which drive to install software. During installation, the hard drive of the Windows NT-based computer was partitioned to create the following drives:


When installing...	Use drive...
programs for general use	C.
◆ Data Collection program database software	D.
◆ GeneScan Analysis software	

Installing the GeneScan Analysis Software

Preparation To prepare for the installation:

Step	Action
1	Check that you have at least 40 MB of free disk space to accommodate the GeneScan Analysis Software.
2	Exit all programs that are running.
3	Turn off any virus protection software that you may have running.



Procedure To install the GeneScan Analysis Software from a CD-ROM:

Step	Action
1	Insert the GeneScan v. 3.7 software CD-ROM into the computer's CD-ROM drive. It will automatically launch the setup when you insert the CD and install the software.
2	<p>If it does not run the setup automatically, you can either:</p> <ul style="list-style-type: none">◆ Double-click Setup (.exe)  or◆ Click Start, point to Run, then browse to the CD-ROM drive and click OK <p>Follow the instructions to install the software.</p> <p>The Setup Complete window opens when the setup has finished installing the software on your computer.</p>


Starting the GeneScan Analysis Software for the First Time

Procedure Follow this procedure when you start the GeneScan Analysis Software the first time or when you start the software after moving the application from the Applied Biosystems folder.

To start the GeneScan Analysis Software for the first time:

Step	Action
1	Before opening the GeneScan Analysis Software for the first time: a. Click Start , point to Settings , and click Printers . b. Right-click the printer you expect to use for GeneScan Analysis Software data. c. Click the Close button in the Printers window.
2	Go to the D:\AppliedBio\GeneScan\Bin folder and double-click the GeneScan icon  to start the software. The GeneScan Analysis Software startup screen opens. 

To start the GeneScan Analysis Software for the first time: *(continued)*

Step	Action						
3	<p>The Product Registration dialog box opens.</p> <p>Enter your registration information into the three fields and click OK.</p> <p>Note The registration code is the number you recorded on page 1-8.</p> <div data-bbox="475 358 946 682"></div> <p>This dialog box opens the first time you start the GeneScan Analysis Software, and any time that you move the software to a different disk or partition.</p>						
4	<p>When the Product Registration dialog box closes, the GeneScan Analysis Software menu bar opens.</p> <table><tr><th>For information about...</th><th>See Page</th></tr><tr><td>Creating a Project</td><td>2-9</td></tr><tr><td>Opening Sample Files</td><td>4-7</td></tr></table>	For information about...	See Page	Creating a Project	2-9	Opening Sample Files	4-7
For information about...	See Page						
Creating a Project	2-9						
Opening Sample Files	4-7						

Removing the Software


This section describes how to remove the GeneScan Analysis Software v3.7 from your Windows NT-based computer. The uninstall process deletes all folders and files installed by the GeneScan Analysis Software Installer.

If Files or Folders Have Been Moved

If you have moved the GeneScan Analysis Software files or folders from their original installed locations they may not be found and deleted by the uninstall operation.

Any files that have been added to the application folders, such as those created when the applications are run, are not deleted by the uninstall operation.

Procedure To remove installed GeneScan Analysis Software:

Step	Action
1	<p>Click Start, point to Settings, click Control Panel, and then double-click Add/Remove Programs.</p> <p>The following dialog box opens:</p> 
2	On the Install/Uninstall tab, select GeneScan and click Add/Remove .

To remove installed GeneScan Analysis Software: *(continued)*

Step	Action
3	At the conclusion of the remove operation, an alert box opens with the message whether or not the remove was successful.
4	<p>If files have been moved or added to the GeneScan Analysis Software or AppliedBio folders, the remove operation will be reported as unsuccessful.</p> <p>If this happens do one of the following:</p> <ul style="list-style-type: none">• Examine and delete the remaining files yourself, or• On the Start menu point to Remove GeneScan Analysis v3.7, then select Programs, and Applied Biosystems.

Creating a Project

2

Overview


In This Chapter Topics in this chapter include the following:

Topics	See Page
Creating Projects and Auto-Analysis on 310 Instruments	2-2
Creating Projects and Auto-Analysis on 377 Instruments	2-3
Setting Up for Automatic Analysis on 3100 and 3700 Instruments	2-4
Working with Project Files	2-8
Finding Missing Sample Files	2-14

Creating Projects and Auto-Analysis on 310 Instruments

Process Using the 310 Instrument


The following diagram illustrates data analysis using the ABI PRISM® 310 Genetic Analyzer.



Creating Projects and Auto-Analysis on 377 Instruments

Process Using the 377 Instrument

The following diagram illustrates the process of automatic analysis and project creation for data from the ABI PRISM® 377 DNA Sequencer.



Setting Up for Automatic Analysis on 3100 and 3700 Instruments

Introduction To set up the GeneScan® Analysis Software for automatic analysis after data-collection, you must have previously defined analysis parameters and size standards.

For more information refer to the following sections:

- ◆ “Chapter 5, “Working with Analysis Parameters.”.
- ◆ “Defining the Size Standard” on page 7-3.

Procedure To set up for automatic analysis after using the 3100 and 3700 Data Collection software:

Step	Action
Complete the following steps in the 3100 and 3700 Data Collection software.	
1	Set the GeneScan Run default preferences to auto-analyze and use the pop-up menu to locate and select the GeneScan Analysis Software.
2	<p>In the GeneScan Sample Sheet, for each sample to be analyzed, do the following:</p> <ul style="list-style-type: none">◆ Enter the sample name. This field must be completed for the samples to be active in the Injection or Run Sheet.◆ Indicate which dye is the standard.◆ Select the check box labeled Pres (Present) for each dye/sample you want auto-analyzed.◆ Select any additional check boxes. <p>IMPORTANT If you plan to use the Genotyper® software, you must complete the Sample Info box correctly. For more information, refer to the <i>Genotyper User's Manual</i>.</p>

To set up for automatic analysis after using the 3100 and 3700 Data Collection software: *(continued)*

Step	Action						
3	<p>In the Injection List or Run Sheet for each applicable sample, select AutoAnalyze, and choose one of the following from one of the appropriate pop-up menus:</p> <ul style="list-style-type: none">◆ Matrix file◆ Analysis parameters◆ Size standard <p>The following table lists considerations when choosing a size standard.</p> <table><tr><th>If you choose a...</th><th>Then GeneScan...</th></tr><tr><td>Size Standard definition file in the Collection Injection List or Run Sheet, but you do not specify which dye is the standard in the Sample Sheet</td><td>performs analysis using the dye specified in the Auto-Analysis defaults.</td></tr><tr><td>dye as the standard in the Sample Sheet, but do not specify a size standard definition file in the Collection Injection List or Run Sheet</td><td>does not perform size calling.</td></tr></table>	If you choose a...	Then GeneScan...	Size Standard definition file in the Collection Injection List or Run Sheet, but you do not specify which dye is the standard in the Sample Sheet	performs analysis using the dye specified in the Auto-Analysis defaults.	dye as the standard in the Sample Sheet, but do not specify a size standard definition file in the Collection Injection List or Run Sheet	does not perform size calling.
If you choose a...	Then GeneScan...						
Size Standard definition file in the Collection Injection List or Run Sheet, but you do not specify which dye is the standard in the Sample Sheet	performs analysis using the dye specified in the Auto-Analysis defaults.						
dye as the standard in the Sample Sheet, but do not specify a size standard definition file in the Collection Injection List or Run Sheet	does not perform size calling.						
4	Select Auto Print to print automatically.						
Complete the following steps:							
1	<p>If you are using the Auto-Analysis defaults, choose Auto-Analysis Defaults from the Settings menu.</p> <p>The Auto-Analysis Defaults dialog box opens.</p> <p>The parameters set in the Auto-Analysis Defaults dialog box apply when you have:</p> <ul style="list-style-type: none">◆ Specified <Analysis Defaults> in the Data Collection software.◆ Selected the check box labeled Always Override Collection Settings.						

To set up for automatic analysis after using the 3100 and 3700 Data Collection software: *(continued)*

Step	Action						
2	<p>Select the check box labeled Always Override Collection Settings to have the new parameters take precedence over the following files specified in the 3100 and 3700 Data Collection programs:</p> <ul style="list-style-type: none"> ◆ Size standard ◆ Analysis parameters ◆ Dye standard 						
3	<p>Select a new Size Standard if:</p> <ul style="list-style-type: none"> ◆ Analysis Defaults is specified for the size standard in the data collection settings, or ◆ To override the specified size standard. 						
4	<p>Select the Dye pop-up menu to specify the dye that represents the internal size standard you are running with the samples.</p>						
5	<p>Select the analysis parameters you want to use from the Parameters pop-up menu.</p> <p>To use the default parameters, select Analysis Parameters.</p>						
6	<p>In the Auto-Print section, there are the following options:</p> <table border="1"> <thead> <tr> <th>Choose...</th><th>To print...</th></tr> </thead> <tbody> <tr> <td>Show Electropherograms</td><td> <p>electropherograms.</p> <p>Select the appropriate radio button to specify whether the electropherograms for the four dyes appear:</p> <ul style="list-style-type: none"> ◆ Together in one panel (overlaid), or ◆ In separate panels (tiled). </td></tr> <tr> <td>Show Tabular data</td><td>tabular data.</td></tr> </tbody> </table>	Choose...	To print...	Show Electropherograms	<p>electropherograms.</p> <p>Select the appropriate radio button to specify whether the electropherograms for the four dyes appear:</p> <ul style="list-style-type: none"> ◆ Together in one panel (overlaid), or ◆ In separate panels (tiled). 	Show Tabular data	tabular data.
Choose...	To print...						
Show Electropherograms	<p>electropherograms.</p> <p>Select the appropriate radio button to specify whether the electropherograms for the four dyes appear:</p> <ul style="list-style-type: none"> ◆ Together in one panel (overlaid), or ◆ In separate panels (tiled). 						
Show Tabular data	tabular data.						
7	<p>Click OK.</p> <p>The GeneScan Analysis Software is now prepared to perform automatic analysis.</p> <p>Note The settings that you specify are initial settings. Make two or three trial runs, fine tuning the parameters with each run, to determine which parameters work the best for a particular protocol.</p>						

What Is a Project

A project is a file containing references to a set of sample files that you want to analyze and display together. The project contains Analysis Control (Ctrl+1) and Results Control (Ctrl+2) windows that allow you to analyze specific dye/samples and display the results of analysis.

Why Create a Project

You can create a new project and add any combination of sample files, allowing you to analyze and display samples from different runs.

Adding a sample file to the project sets up a link between the project and the sample file. The file itself is not imported into the project.

For more information, see “Creating a New Project” on page 2-9.

Where to Store Projects



Keep the sample files in the same location on the hard disk relative to the project file so the GeneScan® Analysis Software can locate them when the project is opened.

If...	Then...
sample files are moved and the software does not find them when you open the project	see “Finding Missing Sample Files” on page 2-14.

Working with Project Files






Opening an Existing Project

To open an existing project:

Step	Action
1	<p>Select Open (Ctrl+O) from the File menu.</p> <p>The Open Existing dialog box opens.</p> 
2	<p>Click the Project icon.</p> <p>An Open dialog box opens.</p> 
3	<p>In the dialog box, find and select the project you want to open.</p>
4	<p>Click Open.</p> <p>The project opens in an Analysis Control window.</p>


Creating a New Project

To create a new project:


Step	Action						
1	<p>Select New (Ctrl+N) from the File menu.</p> <p>The Create New dialog box opens.</p> <div><p>Create New:</p><div><div> Project</div><div> Analysis Parameters</div><div> Size Standard</div><div> Matrix</div><div></div></div></div>						
2	<p>Click the Project icon.</p> <p>An untitled Analysis Control window opens.</p> <p>Note For more information on using the Analysis Control window, refer to page 3-6.</p>						
3	<p>There are two ways to add files to a project:</p> <table><tr><th>To add...</th><th>Select...</th></tr><tr><td>sample files you select to the open project</td><td><p>Add Sample Files (Ctrl+B) from the Project menu.</p><p>The Add Sample dialog box opens (refer to “Using the Add Sample Dialog Box” below).</p></td></tr><tr><td>a sample currently open to the open project</td><td><p>Add “file name” from the Project menu.</p></td></tr></table>	To add...	Select...	sample files you select to the open project	<p>Add Sample Files (Ctrl+B) from the Project menu.</p> <p>The Add Sample dialog box opens (refer to “Using the Add Sample Dialog Box” below).</p>	a sample currently open to the open project	<p>Add “file name” from the Project menu.</p>
To add...	Select...						
sample files you select to the open project	<p>Add Sample Files (Ctrl+B) from the Project menu.</p> <p>The Add Sample dialog box opens (refer to “Using the Add Sample Dialog Box” below).</p>						
a sample currently open to the open project	<p>Add “file name” from the Project menu.</p>						

Using the Add Sample Dialog Box

To use the Add Sample dialog box:

Step	Action
1	<p>When the Add Sample Files dialog box opens, find the folder containing the samples that you want to add from the Look in pop-up menu.</p> 

To use the Add Sample dialog box: *(continued)*

Step	Action										
2	<p>You can take the following action:</p> <table> <tr> <th>If you want to...</th><th>Then...</th></tr> <tr> <td>select a single sample file</td><td>double-click the file or select the file and click Add.</td></tr> <tr> <td>select all the sample files</td><td>click Add All.</td></tr> <tr> <td>add a random selection of sample files</td><td>click each file name while pressing the Ctrl key.</td></tr> <tr> <td>add a continuous list of sample files</td><td> <p>a. Click the first sample that you want to add.</p> <p>b. Press the Shift key and click the last sample you want to add.</p> <p>All the files between the first and last file are selected.</p> </td></tr> </table>	If you want to...	Then...	select a single sample file	double-click the file or select the file and click Add .	select all the sample files	click Add All .	add a random selection of sample files	click each file name while pressing the Ctrl key.	add a continuous list of sample files	<p>a. Click the first sample that you want to add.</p> <p>b. Press the Shift key and click the last sample you want to add.</p> <p>All the files between the first and last file are selected.</p>
If you want to...	Then...										
select a single sample file	double-click the file or select the file and click Add .										
select all the sample files	click Add All .										
add a random selection of sample files	click each file name while pressing the Ctrl key.										
add a continuous list of sample files	<p>a. Click the first sample that you want to add.</p> <p>b. Press the Shift key and click the last sample you want to add.</p> <p>All the files between the first and last file are selected.</p>										
3	<p>Click Finish when you have added all the sample files.</p> <table> <tr> <th>If the...</th><th>Then...</th></tr> <tr> <td>sample files appear in the Analysis Control window</td><td>refer to “Analyzing Sample Files: Using the Analysis Control Window” on page 3-6.</td></tr> <tr> <td>locked files alert appears (see below)</td><td>go to “Unlocking Sample Files” below.</td></tr> </table> <div data-bbox="548 982 1220 1190">  </div>	If the...	Then...	sample files appear in the Analysis Control window	refer to “Analyzing Sample Files: Using the Analysis Control Window” on page 3-6.	locked files alert appears (see below)	go to “Unlocking Sample Files” below.				
If the...	Then...										
sample files appear in the Analysis Control window	refer to “Analyzing Sample Files: Using the Analysis Control Window” on page 3-6.										
locked files alert appears (see below)	go to “Unlocking Sample Files” below.										


Unlocking Sample Files

Adding Locked Sample Files to a Project

If...	Then...
the sample files added to a project are locked	the GeneScan Analysis Software does not allow changes to them. You cannot analyze locked files. When you add locked files, an alert appears.

Procedure


To unlock files:

Step	Action				
1	Select Save Project (Ctrl+S) from the File menu and close the project.				
2	Click the Start button, and then point to Programs .				
3	Click Windows NT Explorer .				
4	Select the applicable sample files and select Properties from the File menu. The Properties dialog box opens.				
5	In the Attributes section, unselect the Read-only check box and click OK . 				
6	Close the window and reopen the project. <table><tr><th>If...</th><th>Then...</th></tr><tr><td>you did not close the project before unlocking the files</td><td>once the files are unlocked, close and open the project so the GeneScan Analysis Software recognizes that the sample files are unlocked.</td></tr></table>	If...	Then...	you did not close the project before unlocking the files	once the files are unlocked, close and open the project so the GeneScan Analysis Software recognizes that the sample files are unlocked.
If...	Then...				
you did not close the project before unlocking the files	once the files are unlocked, close and open the project so the GeneScan Analysis Software recognizes that the sample files are unlocked.				

Removing Samples from a Project

Note A removed sample file is not deleted from the hard disk. The reference is removed from the project.

To remove sample files from a project:

Step	Action
1	Select the file or files in the Analysis Control window (Ctrl+1) or the Results Control window (Ctrl+2) that you want to remove. Note Ctrl+click to select multiple files.
2	Select Remove Sample Files from the Project menu, or press the Delete key. A warning dialog box opens.  A screenshot of a Windows-style warning dialog box titled "GeneScan 3.7". It features a yellow warning triangle icon on the left. The text inside the dialog reads: "1 Sample File(s) are selected. Are you sure you want to remove them?". At the bottom right, there are two buttons: "Cancel" and "Remove". The "Remove" button is highlighted with a dashed border.
3	Click Remove .

Finding Missing Sample Files

When Are Files Considered Lost The project and related sample files are usually located in the Run folder created by the Data Collection software.

If...	Then...
the sample files or the project are moved so they are no longer in the same relative position	<p>the GeneScan Analysis Software might not be able to locate the sample files when the project is opened.</p> <p>Usually, this occurs only if the sample files are moved to another disk drive, another server on a network, or another disk partition on the hard drive.</p>

When an Alert Appears The following table describes what happens when the software does not locate the sample files associated with a project:

If...	Then...
the GeneScan Analysis Software does not locate the sample files associated with a project	an alert box opens and the sample file names appear dimmed when the project opens.

Searching for Missing Sample Files You can re-establish the links between the sample files and the project by choosing Find Missing Sample Files from the Project menu and choosing one of the following options from the submenu:

To find missing files:

Choose...	If you...	Description
Fast Search	suspect that the GeneScan Analysis Software could not find the missing sample files because they are located on an unmounted external storage device, or diskette.	Once the volume is mounted, or the diskette containing the files is inserted, Fast Search finds them immediately.

To find missing files: *(continued)*

Choose...	If you...	Description
Search a Folder	know what folder contains the missing sample files.	Specify a folder in which to search for missing sample files. The GeneScan Analysis Software then immediately locates the sample files and re-establishes links to the project.
Exhaustive Search	do not know where any of the specified missing sample files are located.	The GeneScan Analysis Software searches all mounted disk drives, and available servers. When the files are found, the software re-establishes links to the project.

Re-Establishing Links

To re-establish the links with sample files:

Step	Action
1	Click the dimmed file name to select it. When the missing file is selected the Find “file name” command becomes active. The name of the missing file appears inside the quotation marks.
2	Select Find “file name” . A file dialog box opens.
3	In the file dialog box, locate the proper folder and file.
4	Click Open or double-click the file name.

Procedure To set up for automatic analysis after data collection, complete the following steps in the data collection software:

Step	Action						
1	Set the GeneScan Run default preferences to auto-analyze and use the pop-up menu to locate and select the GeneScan Analysis Software.						
2	<p>In the GeneScan Sample Sheet, for each sample to be analyzed, do the following:</p> <ul style="list-style-type: none"> ◆ Enter the sample name. This field must be completed for the samples to be active in the Injection or Run Sheet. ◆ Indicate which dye is the standard. ◆ Select the check box labeled Pres (Present) for each dye/sample you want auto-analyzed. ◆ Select any additional check boxes. <p>IMPORTANT If you plan to use the Genotyper® software, then you must complete the Sample Info box correctly. For more information, refer to the <i>Genotyper User's Manual</i>.</p>						
3	<p>In the Injection List or Run Sheet for each applicable sample, select AutoAnalyze, and choose from the appropriate pop-up menus one of the following:</p> <ul style="list-style-type: none"> ◆ Matrix file or dye set name ◆ Analysis parameters ◆ Size standard <p>The following table lists considerations when choosing a size standard.</p> <table border="1"> <thead> <tr> <th>If you choose a...</th><th>Then GeneScan...</th></tr> </thead> <tbody> <tr> <td>Size Standard definition file in the Collection Injection List or Run Sheet, but you do not specify which dye is the standard in the Sample Sheet</td><td>performs analysis using the dye specified in the Auto-Analysis defaults.</td></tr> <tr> <td>dye as the standard in the Sample Sheet, but do not specify a size standard definition file in the Collection Injection List or Run Sheet</td><td>does not perform size calling.</td></tr> </tbody> </table>	If you choose a...	Then GeneScan...	Size Standard definition file in the Collection Injection List or Run Sheet, but you do not specify which dye is the standard in the Sample Sheet	performs analysis using the dye specified in the Auto-Analysis defaults.	dye as the standard in the Sample Sheet, but do not specify a size standard definition file in the Collection Injection List or Run Sheet	does not perform size calling.
If you choose a...	Then GeneScan...						
Size Standard definition file in the Collection Injection List or Run Sheet, but you do not specify which dye is the standard in the Sample Sheet	performs analysis using the dye specified in the Auto-Analysis defaults.						
dye as the standard in the Sample Sheet, but do not specify a size standard definition file in the Collection Injection List or Run Sheet	does not perform size calling.						
4	Select Auto Print to print automatically.						

Analyzing Project Files

3

Overview

In This Chapter Topics in this chapter include the following:

Topics	See Page
Analyzing Project Files: About the Analysis Control Window	3-2
Analyzing Sample Files: Using the Analysis Control Window	3-6
Defining Folder Locations	3-16


Analyzing Project Files: About the Analysis Control Window

Introduction When you open a project, the Analysis Control window opens. The Analysis Control window is the main window of a project. You can use this window to specify the following for each sample in the project:

- ◆ Dye that represents the size standard you ran with the sample
- ◆ Size standard
- ◆ Analysis parameters
- ◆ Specific dyes to be analyzed
- ◆ Format the document and print the results automatically

For more information, refer to “Analyzing Sample Files: Using the Analysis Control Window” on page 3-6.

Analysis Control Window Example The following is an example of the Analysis Control window:



**Analysis Control
Window Callouts
Described**

The following table describes the callouts shown for the Analysis Control window in the above figure:

Analysis Control window callouts:

Callout	Description						
1	Use these columns to choose the dye colors to analyze and to specify which is the size standard.						
2	Diamonds mark the standards. For more information, refer to “What Are Size Standards” on page 7-2.						
3	Ctrl+click a dye/sample field to specify that dye/sample as the standard.						
4	From the arrow pop-up menu , you can: <ul style="list-style-type: none">◆ Choose Collection Setting.◆ Choose a user-defined size standard.◆ Define a new size standard for that sample. For more information, refer to: <table><tr><th>Topic</th><th>See Page</th></tr><tr><td>Defining the Size Standard</td><td>7-3</td></tr><tr><td>Using Size Standards</td><td>7-9</td></tr></table>	Topic	See Page	Defining the Size Standard	7-3	Using Size Standards	7-9
Topic	See Page						
Defining the Size Standard	7-3						
Using Size Standards	7-9						
5	From the arrow pop-up menu , you can: <ul style="list-style-type: none">◆ Choose Collection Setting.◆ Choose an analysis parameters file. For more information, refer to: <table><tr><th>Topic</th><th>See Page</th></tr><tr><td>About the Analysis Parameters</td><td>5-2</td></tr><tr><td>Using Analysis Parameter Files</td><td>5-13</td></tr></table>	Topic	See Page	About the Analysis Parameters	5-2	Using Analysis Parameter Files	5-13
Topic	See Page						
About the Analysis Parameters	5-2						
Using Analysis Parameter Files	5-13						

Analysis Control window callouts: *(continued)*

Callout	Description						
6	Double-click the size standard text field or the analysis parameters text field to edit the size standard or the analysis parameters.						
	For more information, refer to:						
	<table><tr><th>Topic</th><th>See Page</th></tr><tr><td>Editing an Existing Size Standard</td><td>7-10</td></tr><tr><td>Changing an Existing Analysis Parameters File</td><td>5-17</td></tr></table>	Topic	See Page	Editing an Existing Size Standard	7-10	Changing an Existing Analysis Parameters File	5-17
	Topic	See Page					
Editing an Existing Size Standard	7-10						
Changing an Existing Analysis Parameters File	5-17						
7	<p>A notation appears in this Information Display field when you move the cursor over a sample file name or over a dye color field.</p> <p>Note For more information on how to customize this field, refer to “Displaying Sample and Dye Information” on page 3-11.</p>						


Customizing the Display

The following table explains how you can customize the display by changing the settings:

Note These preferences also apply to the Results Control window. Refer to Chapter 8, “Evaluating Analysis Results.”

To change...	Choose...	For more information
information displayed in the information display field	Project Options from the Settings menu and Sample Info Display from the submenu.	“Displaying Sample and Dye Information” on page 3-11.
the sorting of sample files	Project Options from the Settings menu and Sample File Sorting from the submenu.	“Setting Sample File Sort Order” on page 3-13.
dye indicator code	Preferences from the Settings menu and Dye Indicators from the submenu.	“Setting Dye Indicator Preferences” on page 3-14.
dye color		

Using the Analysis Control Window

To...	Then...	Result
select all the samples	<p>click the upper-left cell.</p> <p>Click here</p> 	All the columns in the Analysis Control window are selected.
select all of one dye color	click the column heading for that dye color.	The column is highlighted for the color selected.
select all dyes for one sample	click the row number.	All the colors in the row are selected.
change all the standards in the size standards column to the same setting	<ol style="list-style-type: none"> Click the arrow in the column heading. Select a file from the pop-up menu. 	The same size standard is displayed for all the samples.
change a size standard in one of the rows	<ol style="list-style-type: none"> Click the arrow in the row. Select a file from the pop-up menu. 	The size standard for the selected row changes.
change all the parameters in the parameters column to the same setting	<ol style="list-style-type: none"> Click the arrow in the row heading. Select a file from the pop-up menu. 	The same parameter is displayed for all the samples.
change a parameter in one of the rows	<ol style="list-style-type: none"> Click the arrow in the row. Select a file from the pop-up menu. 	The parameter for the selected row changes.
apply a choice to selected fields in the size standards or parameters column	<ol style="list-style-type: none"> Click the row in the column containing the information you want to apply and drag down. Select Fill Down (Ctrl+D) from the Edit menu. 	The value in the selected rows changes to the value in the first row selected.

Analyzing Sample Files: Using the Analysis Control Window

Introduction This section describes using the Analysis Control window to perform the following tasks:

Topic	See Page
Accessing Sample Files	3-6
Analyzing Sample Files	3-7
Specifying the Format for Printed Results	3-8
Displaying Size Standards and Analysis Parameters	3-10
Displaying Sample and Dye Information	3-11
Setting Sample File Sort Order	3-13
Setting Dye Indicator Preferences	3-14

Accessing Sample Files There are two ways to access sample files contained in a project from the Analysis Control window.

Note Sample files that are dimmed can not be found by the project. To find missing files, refer to “Finding Missing Sample Files” on page 2-14.

You can...	Then...	
double-click a sample file name	If the sample file is...	Then that...
	open	Sample File window becomes active.
	not open	sample file opens to its Sample Results view.
select a sample file and select one of the five display modes from the Sample menu	the Sample File window appears in the display mode selected.	

Analyzing Sample Files

The Analysis Control window (Ctrl+1) allows you to analyze multiple samples easily. You choose dyes to analyze, dye standard, size standard, and analysis parameters for each sample file, and then analyze the sample files using these settings.

To analyze sample files:

Step	Action								
1	Click the dye color fields for each sample you want to analyze as follows:								
	<table><tr><th>To...</th><th>Click the...</th></tr><tr><td>select a dye for all samples</td><td>colored column header for that dye.</td></tr><tr><td>select all dyes for a single sample file</td><td>index number at the left end of the row in which the sample file appears.</td></tr><tr><td>all dyes for all samples</td><td>area above the row index numbers.</td></tr></table>	To...	Click the...	select a dye for all samples	colored column header for that dye.	select all dyes for a single sample file	index number at the left end of the row in which the sample file appears.	all dyes for all samples	area above the row index numbers.
	To...	Click the...							
	select a dye for all samples	colored column header for that dye.							
	select all dyes for a single sample file	index number at the left end of the row in which the sample file appears.							
all dyes for all samples	area above the row index numbers.								
2	Identify the sample containing the standard as follows:								
	<table><tr><th>To...</th><th>Ctrl+click the...</th></tr><tr><td>identify each sample that contains a size standard</td><td>colored field that represents the size standard.</td></tr><tr><td>select the same dye as the size standard for all samples</td><td>colored column heading for that dye.</td></tr></table>	To...	Ctrl+click the...	identify each sample that contains a size standard	colored field that represents the size standard.	select the same dye as the size standard for all samples	colored column heading for that dye.		
	To...	Ctrl+click the...							
	identify each sample that contains a size standard	colored field that represents the size standard.							
select the same dye as the size standard for all samples	colored column heading for that dye.								
A diamond appears in the field to identify the dye color as the size standard.									
3	Select a defined size standard setting from the pop-up menu in the Size Standard column as follows:								
	<table><tr><th>To...</th><th>Refer to...</th></tr><tr><td>define a new size standard</td><td>"Using Size Standards" on page 7-9.</td></tr><tr><td>edit a size standard</td><td>"Editing an Existing Size Standard" on page 7-10.</td></tr></table>	To...	Refer to...	define a new size standard	"Using Size Standards" on page 7-9.	edit a size standard	"Editing an Existing Size Standard" on page 7-10.		
	To...	Refer to...							
	define a new size standard	"Using Size Standards" on page 7-9.							
edit a size standard	"Editing an Existing Size Standard" on page 7-10.								
4	To install a new matrix, select a set or all of the samples in the Sample File column and choose Install New Matrix from the Sample menu.								

To analyze sample files: *(continued)*


Step	Action				
5	Select a parameter setting from the pop-up menu in the Parameters column as follows:				
	<table><tr><th>To...</th><th>See...</th></tr><tr><td>use the default analysis parameters to specify different parameters</td><td>"About the Analysis Parameters" on page 5-2.</td></tr></table>	To...	See...	use the default analysis parameters to specify different parameters	"About the Analysis Parameters" on page 5-2.
	To...	See...			
use the default analysis parameters to specify different parameters	"About the Analysis Parameters" on page 5-2.				
6	Select the Print Results check box to print the results automatically. For information on print set-up, refer to "Specifying the Format for Printed Results" on page 3-8.				
7	Click Analyze .				
8	To verify the results, refer to "Process of Verifying Results" on page 9-31.				

**Specifying the
Format for Printed
Results**







To specify the format for printed results:

Step	Action
1	In the Analysis Control window (Ctrl+1), select the check box labeled Print Results . When this check box is selected, the Print Setup button becomes active.

To specify the format for printed results: *(continued)*

Step	Action
2	<p>Click the Print Setup button.</p> <p>The Auto Print Setup dialog box opens (see below).</p> <p>All the dyes selected for analysis are also selected for printing.</p> <p>Sample Dye color fields</p> <div></div> <p>Sample Information Display field</p>
3	<p>Moving the cursor over a Sample File name or over a dye color field, a notation appears in the Sample Information Display field.</p>
4	<p>Click the sample dye color fields to specify any sample you do not want to print as shown above.</p>

To specify the format for printed results: *(continued)*

Step	Action						
5	Select the format by clicking either or both of the buttons at the right of the window.						
	<table><tr><th>Click this button...</th><th>To print...</th></tr><tr><td></td><td>electropherograms for the samples and dyes selected for analysis. Select the appropriate radio button to specify whether the electropherograms for all the dyes appear:<ul style="list-style-type: none">◆ Together in one panel (overlaid), or◆ In separate panels (tiled).</td></tr><tr><td></td><td>tabular data.</td></tr></table>	Click this button...	To print...		electropherograms for the samples and dyes selected for analysis. Select the appropriate radio button to specify whether the electropherograms for all the dyes appear: <ul style="list-style-type: none">◆ Together in one panel (overlaid), or◆ In separate panels (tiled).		tabular data.
	Click this button...	To print...					
	electropherograms for the samples and dyes selected for analysis. Select the appropriate radio button to specify whether the electropherograms for all the dyes appear: <ul style="list-style-type: none">◆ Together in one panel (overlaid), or◆ In separate panels (tiled).						
	tabular data.						
6	Click OK .						

Displaying Size Standards and Analysis Parameters

Use the Analysis Control window to open, review, or change size standard and analysis parameters.

To display these files:

You can...	Then...	For more information see...
double-click the field containing the size standard or the analysis parameters file	the Size Standard or Analysis Parameter window opens.	<ul style="list-style-type: none">◆ “Using Analysis Parameter Files” on page 5-13.◆ “Defining the Size Standard” on page 7-3.

**Displaying Sample
and Dye
Information**


How to Display Sample and Dye Information

Move the cursor over a dye color field or over a sample file name field to display information about the samples and dyes.

How to Specify the Information Displayed

The following procedure describes how to specify the information displayed when moving the cursor over a dye color field or a sample file name field:

To specify the information displayed:

Step	Action										
1	<p>Select Project Options from the Settings menu and Sample Info Display from the submenu.</p> <p>The Sample Info Display dialog box opens.</p> 										
2	<p>Select the check boxes in the Dye/Sample Info & Legend section to control what appears when you move the cursor over the dye color fields in the Control windows.</p> <p>The following table describes the check boxes:</p> <table><tr><th>If you select this check box...</th><th>This appears...</th></tr><tr><td>File Name</td><td>Sample file name.</td></tr><tr><td>Sample Name</td><td>Name of the sample file from the sample file.</td></tr><tr><td>Sample Info</td><td>Sample information from the sample file.</td></tr><tr><td>Comment</td><td>Comment from the sample file.</td></tr></table>	If you select this check box...	This appears...	File Name	Sample file name.	Sample Name	Name of the sample file from the sample file.	Sample Info	Sample information from the sample file.	Comment	Comment from the sample file.
If you select this check box...	This appears...										
File Name	Sample file name.										
Sample Name	Name of the sample file from the sample file.										
Sample Info	Sample information from the sample file.										
Comment	Comment from the sample file.										


To specify the information displayed: *(continued)*

Step	Action												
3	In the Sample File Info heading, select the check boxes for the information that you want to display when you move the cursor over the Sample File name field.												
	The following table describes the buttons:												
	<table><tr><th>Select this button...</th><th>To display...</th></tr><tr><td>Show User Name</td><td>user name from the sample file.</td></tr><tr><td>Show Instrument Name</td><td>instrument name from the sample file.</td></tr><tr><td>Show Path Name</td><td>path and name of the folder where the file is located.</td></tr><tr><td>Show Run Date</td><td>run date and start time from the sample file.</td></tr><tr><td>Show Creation Date</td><td>date and time the sample file was created.</td></tr></table>	Select this button...	To display...	Show User Name	user name from the sample file.	Show Instrument Name	instrument name from the sample file.	Show Path Name	path and name of the folder where the file is located.	Show Run Date	run date and start time from the sample file.	Show Creation Date	date and time the sample file was created.
	Select this button...	To display...											
	Show User Name	user name from the sample file.											
	Show Instrument Name	instrument name from the sample file.											
	Show Path Name	path and name of the folder where the file is located.											
Show Run Date	run date and start time from the sample file.												
Show Creation Date	date and time the sample file was created.												
4	Select the check box labeled Save as Defaults to have the options you choose saved as the default settings.												
5	Click OK .												

**Setting Sample
File Sort Order**

Use this option to specify how sample files are sorted using three criteria. If no sorting option is specified, then the program sorts by sample number.


To set sample file sort order:

Step	Action
1	<p>Select Project Options from the Settings menu and Sample File Sorting from the submenu.</p> <p>The Sample File Sorting dialog box opens.</p> 
2	<p>From the pop-up menus, select from the following items:</p> <ul style="list-style-type: none">◆ File Name◆ Directory◆ Sample Number◆ User Name◆ Instrument Name◆ Run Date◆ Creation Date◆ As Added (sorts the files in the order that they were added) <p>The precedence indicates the sorting level.</p>
3	<p>Select a button for each item to indicate whether to sort in ascending or descending order.</p>
4	<p>Select the check box labeled Save as Defaults to have the options you choose saved as the default settings.</p>
5	<p>Click OK.</p>

Setting Dye Indicator Preferences

The following procedure describes how to change the defaults that determine what dye colors appear on the screen and on printed results. Setting default dye and plot colors sets the colors used for both the Control windows and the Results displays:

To set default dye and plot colors:

Step	Action						
1	<p>Select Preferences from the Settings menu and Dye Indicators from the submenu.</p> <p>The Preferences window opens.</p> <p>If the Preferences window is already displayed, select Dye Indicators from the pop-up menu.</p> <div><div>Vertical scroll bar</div></div> <p>Note Use the scroll bar to see the orange dye color.</p>						
2	<p>The following table describes the Dye Color and Plot Color columns:</p> <table><tr><th>Item</th><th>Description</th></tr><tr><td>Dye Color column</td><td>Shows the colors that represent the dyes in the Control window list. The dye color is also identified in the left color legend in the Results display.</td></tr><tr><td>Plot Color column</td><td>Shows the colors used for plotting the data in the electropherograms.</td></tr></table>	Item	Description	Dye Color column	Shows the colors that represent the dyes in the Control window list. The dye color is also identified in the left color legend in the Results display.	Plot Color column	Shows the colors used for plotting the data in the electropherograms.
Item	Description						
Dye Color column	Shows the colors that represent the dyes in the Control window list. The dye color is also identified in the left color legend in the Results display.						
Plot Color column	Shows the colors used for plotting the data in the electropherograms.						
3	<p>Use the vertical scroll bar to change the dye color and plot color for a fifth dye.</p>						

To set default dye and plot colors: *(continued)*

Step	Action						
4	To change a code, type a different character in the appropriate entry field in the Code column.						
5	Take the following action to change a color: <table><tr><th>To...</th><th>Then...</th></tr><tr><td>change a color</td><td>select a new color from the pop-up menu.</td></tr><tr><td>define a new color</td><td>Select Other from the pop-up menu. A color picker opens.</td></tr></table>	To...	Then...	change a color	select a new color from the pop-up menu.	define a new color	Select Other from the pop-up menu. A color picker opens.
To...	Then...						
change a color	select a new color from the pop-up menu.						
define a new color	Select Other from the pop-up menu. A color picker opens.						
6	Click OK when finished changing dye indicator preferences.						

Defining Folder Locations

Introduction The GeneScan Analysis Software looks in the designated folders for the:

- ◆ Size Standard file
- ◆ Analysis Parameter file
- ◆ Matrix file

When saving one of these files for the first time, the default folder locations for saving the files are those same designated folders.

Storing Matrix Files Store matrix files, that are intended for use by Data Collection software to assign to collection runs, in the AppliedBio folder. The AppliedBio folder is located on the computer on which the Data Collection software is installed.

If Data Collection and Analysis are Performed on Different Computers

Make a copy of a matrix and store it as follows. This is useful when data collection and analysis are performed on different computers.

Store a copy in the...	For use by the...
AppliedBio folder	Data Collection software.
GS Matrix Folder	GeneScan Analysis Software.

Note The ABI PRISM® instrument Data Collection software uses the files installed by the GeneScan Analysis Software in the AppliedBio folder. When you run the analysis software, the program also creates several files (such as a Preference file) and an Analysis Log.

Analyzing Sample Files

4

Overview

In This Chapter Topics in this chapter include the following:

Topics	See Page
About Sample Files	4-2
Converting Macintosh Computer Sample Files	4-3
Converting Sample Files for Use on a Different Platform	4-4
Opening Sample Files	4-7
About the Sample File Window	4-8
Sample Results View	4-9
Sample Info View	4-11
Size Curve View	4-20
Raw Data View	4-22
EPT Data View	4-24
Analyzing a Sample File	4-26

About Sample Files

Sample File Generation on 377 The Gel Processor bundled with 377 Data Collection software creates Sample files after extracting lanes. The information from each lane in the gel file is tracked and extracted, and the resulting Sample files are placed in their respective sample folder. If you change tracking, lane assignment, or Sample Sheet information, you have to regenerate the Sample files.

The software consults Sample Sheet information to determine whether a lane is used (contains sample). The lane tracker uses this information to assign lane numbers to the tracker lines. In addition, the Gel Processor Software only extracts those lanes identified as Used.

What Files Contain Sample files contain electrophoresis data collected on ABI PRISM® instruments. Unanalyzed data contain raw data.

Sample Files Refer to the following table on how Sample files are generated:

Instrument	Software That Generates Sample Files	Notes
ABI PRISM 310	310 Data Collection	Sample file is created after each injection. Data Collection software invokes GeneScan software for auto-analysis
ABI PRISM® 377 DNA Sequencer	Gel Processor software	Gel file is created by 377 Data Collection software. Gel Processor tracks gels and extracts data into Sample files. Gel Processor invokes GeneScan for auto-analysis.
ABI PRISM® 3100 Genetic Analyzer ABI PRISM® 3700 DNA Analyzer	Data Extractor	Sample files are extracted by Data Extractor after each run is completed. If set for auto-analysis, the software is invoked to analyze the data.

GeneScan software is used to manually analyze, edit and view sample files generated on any ABI PRISM instrument.

How GeneScan Analyzes Sample Files

The GeneScan® Analysis Software performs the following steps in analyzing sample files:

Step	Action
1	Processes the raw data signals to generate analyzed data signal and then uses the analyzed signals to detect the signal peaks associated with DNA fragments.
2	Performs size calling by identifying the peaks of the in-lane size standard found in each sample.
3	<p>Determines the fragment size of each experimental peak within the sample based on the size calling curve generated using the size standard peaks, the selected size calling method, and by comparing it to the pre-defined size standard file.</p> <p>The algorithmic steps of the process from raw data to analyzed data are as follows:</p> <ul style="list-style-type: none">◆ Baselineing◆ Smoothing, if any◆ Peak detection

Converting Macintosh Computer Sample Files

About Converting Files

When you insert the GeneScan Analysis Software version 3.7 for Windows NT® platform CD-ROM into a Macintosh® computer's CD-ROM drive, a folder appears that contains two applications, Mac to Win and Win to Mac.

Use these applications to change sample files created on the Macintosh computer to files that can be read by a Windows NT-based computer, and to change files that were created on a Windows NT-based computer so that they can be read by a Macintosh computer.



Converting Sample Files for Use on a Different Platform

Introduction Use the conversion utilities to change sample files created on the Macintosh® computer to files that can be read by a Windows NT™-based computer, and to change files that were created on a Windows NT-based computer so that they can be read by a Macintosh computer.


IMPORTANT This utility will run only on a Macintosh computer.

Converting Macintosh Files to Windows NT-Based Computer Files

To convert Macintosh sample files for use on the Windows NT-based computer:

Step	Action
1	Insert the <i>ABI PRISM® GeneScan Analysis Software v. 3.7 for Windows NT® platform</i> CD-ROM into the Macintosh computer where your Macintosh sample files are stored.
2	Double-click the CD-ROM icon.
3	Double-click the Sample File Mac to Win icon to start the application. <div> Sample File Mac to Win</div> <p>The following dialog box is displayed:</p> <div></div>

To convert Macintosh sample files for use on the Windows NT-based computer: (continued)



Step	Action
4	<p>Click Run.</p> <p>The following directory dialog box is displayed:</p> 
5	<p>Navigate to the folder that contains the sample files that you want to convert and click Choose.</p> <p>The program will perform the task and automatically quit. The converted sample files will have the extension .ab1.</p>

**Converting
Windows
NT-Based
Computer Files to
Macintosh Files**

To convert Windows NT-based computer sample files for use on the Macintosh computer:

Step	Action
1	Insert the <i>GeneScan Analysis Software v. 3.7 for Windows NT platform</i> CD-ROM into the Macintosh computer where your Windows NT-based sample files are stored.
2	Double-click the CD-ROM icon.

To convert Windows NT-based computer sample files for use on the Macintosh computer: *(continued)*

Step	Action				
3	<p>Click the Sample File Win to Mac icon to start the application.</p> <div> Sample File Win to Mac</div> <p>The following directory dialog box is displayed:</p> <div></div>				
4	<p>Navigate to the folder that contains the sample files that you want to convert and click the Choose button.</p> <table><tr><th>If...</th><th>Then...</th></tr><tr><td>there are no problems, the program will perform the task and automatically quit</td><td>When you open the folder, you can double-click the sample files to open them using GeneScan Analysis Software on a Macintosh computer.</td></tr></table>	If...	Then...	there are no problems, the program will perform the task and automatically quit	When you open the folder, you can double-click the sample files to open them using GeneScan Analysis Software on a Macintosh computer.
If...	Then...				
there are no problems, the program will perform the task and automatically quit	When you open the folder, you can double-click the sample files to open them using GeneScan Analysis Software on a Macintosh computer.				


Opening Sample Files

Introduction Sample files can be opened as separate files outside of projects, and display related information about each sample file.


To open sample files:

If you are interested in...	Then...
one or two sample files	it is often more convenient to open sample files individually and analyze or view the data without opening an entire project.
multiple sample files	use a project. For information on opening Sample Files from within a Project, see "Accessing Sample Files" on page 3-6

Procedure To open a sample file as a separate file:

Step	Action
1	<p>Select Open (Ctrl+O) from the File menu. The Open Existing dialog box appears.</p>  <p>Note You can also double-click the sample file name in the folder containing the files. If the GeneScan Analysis Software is not running, the software starts and opens the sample file.</p>

To open a sample file as a separate file: *(continued)*

Step	Action
2	Click the Sample icon. An Open dialog box appears. 
3	In the dialog box, navigate to the folder and select the sample file that you want to open.
4	Click Open . The Sample File window appears. For more information on the Sample File window, refer to page 4-8.

About the Sample File Window

What It Displays You can use the display modes in the Sample File window to review the analyzed and raw data, and all pertinent data collection, sizing and sample description information from a single window. The Sample Results view appears as the default.


Five Views The five views of the Sample File window are:

View	See Page
Sample Results View	4-9
Sample Info View	4-11
Size Curve View	4-20
Raw Data View	4-22
EPT Data View	4-24

Sample Results View


What It Displays The Sample Results View displays the sample file's analyzed data in both electropherogram and tabular data form.

Displaying the View The following table describes ways to display the Sample Results View:

To display the view...	Do this...	
from the Sample File window	You can either	Result
	click the button for the Sample Results view at the bottom left of the Sample File window. 	The Sample Results View appears. Refer to “Example of Sample Results View” on page 4-10.
	select Sample Results (Ctrl+E) from the Sample menu.	
from a project window	a. Select a sample or Ctrl+click to select multiple samples in the Analysis Control (Ctrl+1) or the Results Control (Ctrl+2) window. b. Select Sample Results (Ctrl+E) from the Sample menu.	

**Example of
Sample Results
View**

The following is an example of the Sample Results view:



**Description of
Columns**

The following table describes the columns in the above figure:

This column	Identifies
Dye/Sample Peak	Dye color and Peak number.
Minutes	The time, in minutes, from the start of the run to the time the fragment was detected.
Size	<div>The number of base pairs in the fragment.</div> <div>This value is calculated automatically only if you:</div> <div><div>◆ Run the size standard in the same lane or injection as the sample, and</div><div>◆ Perform size calling.</div></div>
Peak Height	Signal size.
Peak Area	Area of the detected peak.
Data point	Data point of the fragment at its maximum peak height.

Differences from the Results Display

The Sample Results view displays the same electropherogram and tabular data as the Results Display, with the following differences:

- ◆ One sample file is displayed.
 - ◆ Show or hide dye/sample data by clicking the buttons below the electropherogram.
 - ◆ Cannot display legends.
 - ◆ Cannot use custom plot colors.
-
-

Sample Info View


What It Displays

Displays the following sample file information:

Information displayed	See Page
Run Information	4-14
Data Collection Settings	4-14
Sample Information	4-14
Gel Information	4-15
Analysis Records	4-15
Dyes Within Analysis Records	4-16
Data Collection Settings	4-17
Gel Information (Polymer)	4-17
Sample Information	4-18
Analysis Records	4-18
Dyes Within Analysis Records	4-19

Displaying the View

The following table describes ways to display the Sample Info View:

To display the view...	Do this...	
from the Sample File window	You can either	Result
	click the button for the Sample Info view at the bottom left of the Sample File window. 	The Sample Info view appears. Refer to “Example of Sample Info View” on page 4-13.
	select Sample Info (Ctrl+I) from the Sample menu.	
from a project window	a. Select a sample or Ctrl+click to select multiple samples in the Analysis Control (Ctrl+1) or the Results Control (Ctrl+2) window.	
	b. Select Sample Info (Ctrl+I) from the Sample menu. The information is organized in four panels. Click the triangles to expand or collapse the panels to display specific information.	

Example of Sample Info View

The following is an example of the Sample File Information window in Sample Info View:

13C6.fsa

Sample File Information

▼ **Run Information**

User Name: JFred	Run Date: Sat, Dec 16, 2000
Instrument: 377-96 #	Start Time: 2:33:15 AM
Data Coll. Version: ABI PRISM 377-96 Collection 3.0.0b1c5	Run Duration: 2 Hrs 30 Mins 38 Secs
	Total Points: 5708

▼ **Data Collection Settings**

Module File:	E.P. Voltage: 0 Volts
Matrix File: 5Dye Matrix.mtx	E.P. Current: 0 mAmps
Parameters:	E.P. Power: 0 Watts
Size Standard:	Temperature: 0 °C
	Laser Power: 0 mWatts

▼ **Gel Information**

Gel Type:	Gel File: 5Dye 377NT Gel File.gel
Gel Percent: 0.0 %	Gel Thickness: 0.0 mm
Well-To-Read Distance: 36.0 cm	Lane Number: 13
Number of Channels: 194	Channel Averaging: 3
Number of Lanes: 34	Lane Extraction Range: 0 - 5708 pts

▼ **Sample Information**

Sample Name: C6

Dye	Sample Info	Comment
B:	Blue Info	Comment
G:	Green Info	Comment
Y:	Yellow Info	Comment
R:	Red Info	Comment
O:	Orange Info	Comment

▼ **Analysis Records**

B: Analyzed 2:18:25 PM Fri, Jan 05, 2001

Parameters: GS500Analysis.gsp	Standard: GS 500 All.szs	Peak Totals Found In: - Sample: 23 - Dye Std: 35 Std Defined: 16 Std Matched: 16
Analysis Range: 0 - 5707 pts	Dye Std: 0	
Baselined: Yes Multi-Componented: Yes	Size Method: Local Southern Method	
Data Smoothing: None	Size Range: Full Range	
Peak Detection Threshold: 50	Std Peak Det. Threshold: 50	
Peak Detection Min. Half-Width: 2	Baseline Window Size: 261	
Polynomial Degree: 3	Peak Win: 19	Slope Threshold for Pk Start: 0.0 End: 0.0

Description of Information This information is generally applied to Data Collection software. The following tables list the information in the Sample Info View:

Run Information

Information found under the header	Information entered in
<ul style="list-style-type: none"> ◆ User Name ◆ Instrument ◆ Data Collection software version ◆ Run date and start time ◆ Tube ◆ Run Date ◆ Start Time ◆ Run Duration ◆ Total Points 	Data Collection software Run file and run information

Data Collection Settings

Information found under the header	Information is entered in
<ul style="list-style-type: none"> ◆ Module File ◆ Matrix ◆ Analysis Parameters ◆ Size Standard ◆ Run Voltage, Injection Voltage, and Injection Duration ◆ Temperature ◆ Laser power 	Data Collection software Run file and run information

Sample Information

Information found under header	Information is entered in
Sample Name	Data Collection software Plate Setup Record
Sample data and comment for each dye color	

Gel Information

This information applies only to samples from the ABI PRISM 377 instrument.

Information found under header	Information is entered in
<ul style="list-style-type: none">◆ Gel Type◆ Gel File name◆ Gel Percent◆ Gel Thickness◆ Well-To-Read Distance◆ Number of Channels◆ Lane Number◆ Channel Averaging◆ Lane Extraction Range	Data Collection software

Analysis Records

Information found under header	Information is entered in
Date and time each color was analyzed, and more panel-display arrows	Analysis information

Dyes Within Analysis Records

Information found under header	Information is entered
Analysis parameters file and range analyzed	Analysis Settings
Whether baselined or multicomponented	
Data smoothing	
Peak detection threshold and minimum half-width	
Polynomial degree	
Peak window	
Size standard file	Analysis Control window
The dye color used for the standard, Sizing method, and range	
Standard peak detection threshold	Analysis Settings
Baseline window size	
Slope threshold for peak start	
Value for end	
Total number of peaks: found in sample and dye standard, defined in standard matched with standard peaks	Analysis results

Data Collection Settings

Information found under the header	Information is inserted from
<ul style="list-style-type: none"> ◆ Module File ◆ Matrix File ◆ Analysis Parameters ◆ Size Standard ◆ Electrophoresis voltage, current, and power (ABI 373 and ABI PRISM 377 only) ◆ Run Voltage, Injection Voltage, and Injection Duration (ABI PRISM 310 only) ◆ Temperature ◆ Laser power 	<p>ABI 373 or ABI PRISM 377 gel file (data collection Run File and run information).</p> <p>This information is embedded in the gel file.</p> <p>ABI PRISM 310 data collection Run file and run information.</p>

Gel Information (Polymer)

Information found under the header	Information is inserted from
ABI PRISM 310	
◆ Gel type	Optional (user entered)
◆ Length to detector	Data Collection information
◆ Lot # and expiration date	
ABI 373 and ABI PRISM 377	
◆ Gel type	Optional (user entered).
◆ Name of gel file	Gel file (data collection Run file) This information is embedded in the gel file.
◆ Gel percentage and thickness	Optional (user entered)
◆ Well-to-read (separation) distance	
◆ Number of channels	Gel file (data collection Run file)
◆ Number of Lanes	This information is embedded in the gel file.

Information found under the header	Information is inserted from										
<p>♦ Channel Averaging</p> <p>Note Zero (0) indicates use of pre-averaging offscale data.</p> <p>The following is an example of how pre-averaging offscale data appears:</p> <table border="1"> <thead> <tr> <th>Extraction method</th><th>Displayed</th></tr> </thead> <tbody> <tr> <td>3 channel averaging</td><td>3</td></tr> <tr> <td>3 channel averaging with pre-averaging offscale</td><td>30</td></tr> <tr> <td>3 channel averaging-weighted</td><td>- 3</td></tr> <tr> <td>3 channel averaging-weighted, pre-averaging offscale</td><td>- 30</td></tr> </tbody> </table> <p>♦ Range of data points extracted</p>	Extraction method	Displayed	3 channel averaging	3	3 channel averaging with pre-averaging offscale	30	3 channel averaging-weighted	- 3	3 channel averaging-weighted, pre-averaging offscale	- 30	Gel processing parameters
Extraction method	Displayed										
3 channel averaging	3										
3 channel averaging with pre-averaging offscale	30										
3 channel averaging-weighted	- 3										
3 channel averaging-weighted, pre-averaging offscale	- 30										

Sample Information

Information found under header	Information is inserted from
Sample data and comment for each dye color	<p>♦ ABI PRISM 310 Sample Sheet</p> <p>♦ ABI 373 or ABI PRISM 377 gel file (Sample Sheet)</p> <p>This information is embedded in the gel file.</p>

Analysis Records

Information found under header	Information is inserted from
Date and time each color was analyzed, and more panel-display arrows	Analysis information


Dyes Within Analysis Records

Information found under header	Information is inserted from
◆ Analysis parameters file and range analyzed	Analysis Settings
◆ Whether baselined or multicomponented	
◆ Data smoothing	
◆ Peak detection threshold and minimum half-width	
◆ Size standard file	Analysis Control window
◆ The dye color used for the standard, Sizing method, and range	
◆ Standard peak detection threshold	Analysis Settings
◆ Split Peak correction	
◆ Total number of peaks: found in sample and dye standard, defined in standard matched with standard peaks	Analysis results


Size Curve View

What It Displays Displays sizing curves for sample files. The size curve is a measure of how well the internal size standard matches the standard definition, and whether or not it is linear.

Displaying the View The following table describes ways to display the Size Curve View:

To display the view...	Do this...	
from the Sample File window	You can either	Result
	<div>click the button for the Size Curve view at the bottom left of the Sample File window.</div> <div></div> <div>select Size Curve (Ctrl+U) from the Sample menu.</div>	<div>The Size Curve view appears.</div> <div>Refer to “Example of Size Curve View” on page 4-21.</div>
from a project window	<div>a. Select a sample or Ctrl+click to select multiple samples in the Analysis Control (Ctrl+1) or the Results Control (Ctrl+2) window.</div> <div>b. Select Size Curve (Ctrl+U) from the Sample menu.</div> <div>For each selected sample, a Sample File window opens and displays its Size Curve view.</div>	

Example of Size Curve View The Size Curve view displays two curves, as shown in the figure below. For a description of the curves, refer to “Curves Described” below.



Curves Described The following table describes the curves in the above figure:


Note Sizing errors due to anomalous mobilities may be displayed as nonlinear.

This curve	Represents
Red curve	The sizecalling curve, based on the sizecalling method used to analyze the data.
Black curve	<p>The best-fit least squares curve, which the GeneScan Analysis Software calculates for all samples, regardless of the size calling method.</p> <p>This curve is provided to help evaluate the linearity of the sizing curve.</p> <p>When the sizing curve and best-fit curve match, they overlap so you see only the size curve.</p>

Raw Data View


What It Displays It displays the raw data collected for a sample. This information is stored in the sample file.

Displaying the View The following table describes ways to display the Raw Data View:

To display the view...	Do this...	
from the Sample File window	You can either	Result
	<p>click the button for the Raw Data view at the bottom left of the Sample File window.</p>  <p>select Raw Data (Ctrl+R) from the Sample menu.</p>	<p>The Raw Data View appears.</p> <p>Refer to “Example of Raw Data View” on page 4-23.</p>
from a project window	<p>a. Select a sample or Ctrl+click to select multiple samples in the Analysis Control (Ctrl+1) or the Results Control (Ctrl+2) window.</p> <p>b. Select Raw Data (Ctrl+R) from the Sample menu.</p>	

Example of Raw Data View

The following is an example of the Sample File window in Raw Data View:



What to Evaluate

Use the Raw Data view to evaluate:


- ◆ Problems or noise in the baseline that could result in poor size calling.
- ◆ Start and stop points for analysis.

For information on changing the horizontal (refer to page 9-21) or vertical scale of the data (refer to page 9-22).

EPT Data View


What It Displays The EPT Data View (electrophoresis power and temperature) displays this information collected for a sample, and it is stored in the sample file.

Displaying the View The following table describes ways to display the EPT Data View:

To display the view....	Do this...	
from a Sample File window	You can either	Result
	click the button for the EPT Data view at the bottom left of the Sample File window. 	The EPT Data view appears. Refer to “EPT Data View Example” on page 4-25.
	select EPT Data (Ctrl+M) from the Sample menu.	
from a project window	a. Select a sample or Ctrl+click to select multiple samples in the Analysis Control (Ctrl+1) or the Results Control (Ctrl+2) window. b. Select EPT Data (Ctrl+M) from the Sample menu, or click the EPT data button.	

**EPT Data View
Example**

The following is an example of the Sample File window in EPT Data View:



**Colored Lines
Described**


The following table describes the lines in the above figure:

Line	Description
Blue	Electric voltage in volts/10
Black	Electric power in watts
Red	Run temperature in °C
Green	Electric current in mA (milliamps)

Analyzing a Sample File

Introduction The GeneScan Analysis Software analyzes raw data stored in sample files according to parameters and standards that you select. You can use the analyzed data to detect peaks associated with DNA fragments and identify those peaks with an established size standard.

Procedure To analyze a sample file:

Step	Action										
1	<p>Select Analyze “Sample File Name” from the Sample menu (Ctrl+Y).</p> <p>The Analyze Sample file dialog box appears.</p> 										
2	<p>From the Analyze Sample file dialog box, select one of the following options:</p> <table><tr><th>Choose...</th><th>To select...</th></tr><tr><td>Dyes to Analyze check boxes</td><td>any number of dyes to analyze.</td></tr><tr><td>Analyze Parameters pop-up menu</td><td>from the default parameters or any parameter files in the folder location specified in the application preferences.</td></tr><tr><td>Size Standard pop-up menu</td><td>from the default standard, or any standard files in the folder location that you specify in the application preferences.</td></tr><tr><td>Standard Dye pop-up menu</td><td>the inline standard dye.</td></tr></table>	Choose...	To select...	Dyes to Analyze check boxes	any number of dyes to analyze.	Analyze Parameters pop-up menu	from the default parameters or any parameter files in the folder location specified in the application preferences.	Size Standard pop-up menu	from the default standard, or any standard files in the folder location that you specify in the application preferences.	Standard Dye pop-up menu	the inline standard dye.
Choose...	To select...										
Dyes to Analyze check boxes	any number of dyes to analyze.										
Analyze Parameters pop-up menu	from the default parameters or any parameter files in the folder location specified in the application preferences.										
Size Standard pop-up menu	from the default standard, or any standard files in the folder location that you specify in the application preferences.										
Standard Dye pop-up menu	the inline standard dye.										
3	<p>After analysis, evaluate the results.</p> <p>For more information on evaluating the results, refer to Chapter 8, “Evaluating Analysis Results.”</p>										

Installing a New Matrix File

Use the following procedure to install a new matrix file for the Sample file that you want to analyze.

For information on attaching the new matrix to an ABI 373 or ABI PRISM 377 gel file, refer to “Gel Processor User’s Manual.”

To install a new matrix file:

Step	Action
1	Choose Install New Matrix from the Sample menu. A directory dialog box appears. The Folder Preferences settings determine where the GeneScan Analysis Software looks for the matrix file. For more information, refer to “Defining Folder Locations” on page 3-16.
2	Select the new matrix file in the dialog box and click Open . A message appears when the matrix is successfully assigned.
3	Re-Analyze the Sample file. Applying a new matrix file clears previous analysis information, so you must re-analyze the file.

Working with Analysis Parameters

5

Overview

In This Chapter Topics in this chapter include the following:

Topics	See Page
About the Analysis Parameters	5-2
Sizecaller Algorithm Flowchart	5-4
Setting Analysis Parameters	5-5
Using Analysis Parameter Files	5-13

About the Analysis Parameters

What They Are Analysis parameter files are used by 3100 and 3700 Data Collection software to inform the software of the parameters to use to automatically analyze fragment data. Use the analysis parameters provided, or follow the procedure in “Creating GeneScan Analysis Modules” on page A-1 to create new analysis parameter files. The auto-analysis feature is part of the 3100 and 3700 Data Collection software.

These parameters are also used for analysis using the GeneScan® Analysis Software.

Why They Are Necessary Analysis parameter files are required because the auto-analysis feature in 3100 and 3700 software has no user interface, so the parameters for analysis cannot be directly configured in the analysis software. The GeneScan Analysis Software creates the analysis parameters, and the software saves the parameters in a folder that can be read by the Data Collection system software that performs auto-analysis.

The analysis parameter file can also be used for manual analysis of fragment data. It is then used by auto-analysis while it is processing run data.

When to Specify a Parameter Before performing a run, specify the analysis parameter to use for analysis of each sample. Do this when preparing the plate record for the plate that contains the samples on 3100 and 3700 instruments and when preparing the run sheet on 310 and 377 instruments.

Note If there is no analysis parameter specified for a particular sample in the plate record, the sample will not be analyzed by auto-analysis. Use the GeneScan Analysis Software to analyze the data.

**Analysis
Parameter Files
Provided**

The following analysis parameter files are provided:

Analysis parameter file	Suggested parameters for using...
GS120Analysis.gsp	GeneScan 120 Size Standard
GS350Analysis.gsp	GeneScan 350 Size Standard
GS400Cubic Analysis.gsp	GeneScan 400 HD Size Standard
GS400HDAnalysis.gsp	GeneScan 400 HD Size Standard
GS400Ord2Analysis.gsp	GeneScan 400 HD Size Standard
GS500Analysis.gsp	GeneScan 500 Size Standard


**When Analysis
Parameters Are
Used**

Analysis parameter files are opened and read by the Data Extractor program as each sample is extracted.

Note If you make a change to an existing analysis parameter file (without changing its file name), this change will affect any sample extracted after the change is saved.

Sizecaller Algorithm Flowchart

Flowchart The following flowchart shows how the sizecaller algorithm works:



Setting Analysis Parameters

Default Settings The following table describes the default settings for the analysis parameters. Some of these parameters will be different depending on the analysis setting you select. This is an example of one of the settings.

Parameter	Default setting																										
Analysis Range	Full Range																										
Data Processing	Smooth Options—the Light button is selected																										
Peak Detection	Select Analysis Parameters from the Settings menu. The Analysis Parameters dialog box opens. <table><tr><th>Parameter</th><th>Default setting</th></tr><tr><td rowspan="7">Peak Amplitude Thresholds</td><td><table><tr><th>Dye</th><th>Setting</th></tr><tr><td>Blue</td><td>50</td></tr><tr><td>Green</td><td>50</td></tr><tr><td>Yellow</td><td>50</td></tr><tr><td>Red</td><td>50</td></tr><tr><td>Orange</td><td>50</td></tr></table></td></tr><tr><td>Min. Peak Half Width</td><td>2 points</td></tr><tr><td>Degree of Polynomial differentiation</td><td>3</td></tr><tr><td>Peak Window Size</td><td>19 points</td></tr><tr><td>Slope Threshold for Peak Start</td><td>0.0</td></tr><tr><td>Slope Threshold for Peak End</td><td>0.0</td></tr></table>	Parameter	Default setting	Peak Amplitude Thresholds	<table><tr><th>Dye</th><th>Setting</th></tr><tr><td>Blue</td><td>50</td></tr><tr><td>Green</td><td>50</td></tr><tr><td>Yellow</td><td>50</td></tr><tr><td>Red</td><td>50</td></tr><tr><td>Orange</td><td>50</td></tr></table>	Dye	Setting	Blue	50	Green	50	Yellow	50	Red	50	Orange	50	Min. Peak Half Width	2 points	Degree of Polynomial differentiation	3	Peak Window Size	19 points	Slope Threshold for Peak Start	0.0	Slope Threshold for Peak End	0.0
Parameter	Default setting																										
Peak Amplitude Thresholds	<table><tr><th>Dye</th><th>Setting</th></tr><tr><td>Blue</td><td>50</td></tr><tr><td>Green</td><td>50</td></tr><tr><td>Yellow</td><td>50</td></tr><tr><td>Red</td><td>50</td></tr><tr><td>Orange</td><td>50</td></tr></table>	Dye	Setting		Blue	50	Green	50	Yellow	50	Red	50	Orange	50													
	Dye	Setting																									
	Blue	50																									
	Green	50																									
	Yellow	50																									
	Red	50																									
	Orange	50																									
Min. Peak Half Width	2 points																										
Degree of Polynomial differentiation	3																										
Peak Window Size	19 points																										
Slope Threshold for Peak Start	0.0																										
Slope Threshold for Peak End	0.0																										
Full Range	Full Range																										
Sizecalling Method	Local Southern Method																										
Baselining	Baseline Window Size 251 points																										
Auto Analysis	Size Standard <None>																										

**Displaying
Analysis
Parameters**

Select Analysis Parameters from the Settings menu. The Analysis Parameters dialog box opens.

Note To display the orange Peak Amplitude Threshold, use the scroll bar under the values and scroll to the right.

Analysis Parameters

Analysis Range

☒ Full Range

☐ This Range (Data Points)

Start:

Stop:

Data Processing

Smooth Options

☒ None

☐ Light

☐ Heavy

Size Call Range

☒ Full Range

☐ This Range (Base Pairs)

Min:

Max:

Size Calling Method

☐ 2nd Order Least Squares

☐ 3rd Order Least Squares

☐ Cubic Spline Interpolation

☒ Local Southern Method

☐ Global Southern Method

Peak Detection

Peak Amplitude Thresholds

B: Y:

G: R:

Min. Peak Half Width: Pts

Polynomial Degree:

Peak Window Size: Pts

Slope Threshold for Peak Start:

Slope Threshold for Peak End:

Baselining

BaseLine Window Size: Pts

Auto Analysis Only

Size Standard:

Cancel OK

Analysis Parameters

There are seven analysis parameters:

Analysis parameter	See Page
Analysis Range Options	5-7
Data Processing Options	5-7
Peak Detection Options	5-8
Sizecall Range Options	5-10
Sizecalling Method Options	5-10
Baselining Option	5-12
Auto-Analysis Only Option	5-12

Analysis Range Options

The following are the Analysis Range options:

Item	Description
Full Range button	Use to analyze all the data collected on the genetic analysis instrument for each sample.
This Range (Data Points) button	Enter Start and Stop data point numbers in the entry fields in order to specify only a limited range to be analyzed for each sample. This affects what is displayed in the results display. Normally, set the analysis range to start after the primer peak.

Data Processing Options

The Smooth Options only affect the appearance of the analyzed electropherograms.

Note Since the tabulated peaks are calculated from unsmoothed data, they might not be consistent with the smoothed display. For example, the reported peak heights might not be the same as those visible from the smoothed electropherogram. It is recommended that you select the None button.

Peak Detection Options

The Peak Detection options locate peaks at the positive-to-negative zero crossings of the first derivative of the baselined electropherogram. The peak detector computes the first derivative at a data point i by fitting a polynomial to a window centered on i .

Peak Detection Parameter Options Described

The following table describes the options:

Item	Description	For example
Peak Amplitude Thresholds	<p>The sizecaller reports to the user only those peaks whose heights are at least the Peak Amplitude Threshold for that dye.</p> <p>Set the dye amplitude threshold at a level that allows the software to detect peaks, but eliminate noise.</p> <p>For each dye, the GeneScan Analysis Software detects peaks above the threshold entered in the entry field.</p>	<p>If you leave the default value of 50, peaks with amplitude above 50 are analyzed and appear in the tabular data.</p> <p>Lower amplitude peaks still appear in the electropherogram, but are not analyzed and do not appear in the tabular data.</p>
Minimum Peak Half Width	<p>Defines what constitutes a peak.</p> <p>Use to specify the smallest full width at half maximum for peak detection.</p> <p>The range is 2–99.</p> <p>A typical number might be 3 for microsatellites, or 10 for SSCPs.</p>	<p>If this number is large, the software ignores noise spikes.</p> <p>If the peaks in the data are narrow, set the value to a low number.</p> <p>Experiment with this value to determine the best number for the data.</p>

The following table describes the options: *(continued)*

Item	Description	For example
Polynomial Degree	Sets the degree of the polynomial.	These parameters control the sensitivity of this process. Sensitivity increases with the polynomial degree and decreases with the window size.
	Min. setting	
	2	
Peak Window Size	Max. setting	Use polynomials of degree 2 or 3 for well-isolated peaks, such as those from a size standard, and a degree 4 for finer control.
	5	
	Sets the width of the window.	
Min. setting	Max. setting	For degree 4, the Peak Window Size should be 1 to 2 times the full width at half maximum of the peaks that you wish to detect. These parameters cannot be set for each color independently.
1 above the Degree of Polynomial differentiation setting.	Number of scans between peaks.	
Slope Threshold for Peak Start Slope Threshold for Peak End	Determines where a peak starts and stops.	For example, a peak ends when the first derivative again exceeds the Slope Threshold for Peak End. Slope Threshold for peak start must be non-negative and Slope Threshold for peak end must be nonpositive. Values other than 0 will move the extent of the peak toward its center.

The following table describes the options: *(continued)*

Item	Description	For example
Baseline Window Size	<p>Controls the slope of the baseliner.</p> <p>The sizecaller determines a baseline value for each scan i.</p> <p>Basically, the sizecaller sets the baseline to the lowest electropherogram value that it sees in a window of size β centered on scan i.</p>	<p>To achieve symmetry, you should set β to an odd value, though even values are acceptable.</p> <p>Small β will cause the baseline to creep up into peaks.</p> <p>Large β will create a baseline that does not touch all peaks; that is, the peaks will not be baseline resolved after the sizecaller subtracts the baseline from the electropherogram.</p>

Sizecall Range Options About the Sizecall Range Options

Use the Sizecall Range parameter options to specify the range of size fragment (in base pairs) to be included in the peak tabular data.

Sizecall Range Parameter Options Described

Item	Description
Full Range button	Select this choice to report all peaks within the range of the matched size standard.
This Range (Base Pairs) button	Select this choice to limit the reported peaks by fragment size.

Sizecalling Method Options About Sizecalling Method Options

Click a button to select the desired sizecalling method. The GeneScan Analysis Software uses these methods to determine the molecular length of an unknown fragment.

Sizecalling Method Parameter Options Described

Item	Description
2nd Order Least Squares and 3rd Order Least Squares	<p>Both Least Squares Methods use regression analysis to build a best-fit sizecalling curve.</p> <p>For information on the Sizecalling Methods, refer to page B-1.</p>
Cubic Spline Interpolation	<p>Forces the sizing curve through all the known points of the selected GeneScan size standard.</p> <p>For information on the Cubic Spline Interpolation Method, refer to page B-4.</p>
Local Southern Method	<p>Determines the sizes of fragments by using the reciprocal relationship between fragment length and mobility.</p> <p>For information on the Local Southern Method, refer to page B-5.</p>
Global Southern Method	<p>Similar to the Least Squares Method in that it compensates for standard fragments that may run anomalously.</p> <p>For information on the Global Southern Method, refer to page B-7.</p>

Baselining Option About the Baselining Option

The Baselining option controls the scope of the baseliner. Use this option to set the size *Beta* of the Baseline Window. The sizecaller computes a baseline for the electropherogram of each dye independently.

How the Baselining Option Works

A baseline comprises a value at each data point *i*. Basically, the baseline value at each data point *i*, is the lowest electropherogram value in a window whose width *Beta* is set using the Baselining option, and centered at each data point *i*.

More accurately, the baseline computed in this manner is intermediate. The real baseline value at each data point *i*, is the highest intermediate value, again in a window whose width *Beta* is set using the Baselining options and centered at each data point *i*. The sizecaller baselines an electropherogram by subtracting the baseline from the raw electropherogram.

If the Baseline Window Is Too Small or Too Large

The following table describes what happens if the baseline window is either too small or too large:

Using...	Causes...
a small baseline window size	the baseline to creep into the peaks, resulting in shorter peaks in the analyzed data.
a large baseline window size	the baseline to ride too low, resulting in elevated and possibly not baseline-resolved peaks.

Auto-Analysis Only Option	When performing auto-analysis, select the Size Standard to use from this pop-up menu.
--------------------------------------	---


Using Analysis Parameter Files

In This Section This section contains the following topics:

Topic	See Page
Assigning the Same Analysis Parameters to All Files	5-13
Assigning Different Parameters to Single Samples	5-15
Displaying Default Parameters	5-16
Creating Custom Analysis Parameter Files	5-16
Changing an Existing Analysis Parameters File	5-17
Deleting Custom Analysis Parameters	5-17

Assigning the Same Analysis Parameters to All Files

To analyze samples using the same analysis parameters:


Step	Action
1	If the Analysis Control window is not displayed, then select Analysis Control (Ctrl+1) from the Windows menu.
2	<p>Click the arrow in the Parameters column heading and choose parameters from the pop-up menu.</p> <p>Your menu choice applies to all fields in the column.</p> <p>Note Selecting this pop-up menu automatically selects ALL samples. You cannot assign parameters to only a subset of samples.</p> <div></div>

To analyze samples using the same analysis parameters: *(continued)*

Step	Action
3	The pop-up menu contains the following options:

Assigning Different Parameters to Single Samples

To apply separate analysis parameters to selected samples:

Step	Action						
1	If the Analysis Control window is not displayed, then select Analysis Control (Ctrl+1) from the Windows menu.						
2	<p>Click the arrow in the Parameters column for the sample parameter settings that you want to change.</p> <p>A pop-up menu opens.</p> 						
3	<p>The pop-up menu contains the following options:</p> <table border="1"> <thead> <tr> <th>Choose...</th><th>To...</th></tr> </thead> <tbody> <tr> <td>Analysis Parameters</td><td>apply the parameters that are stored as preferences in the software.</td></tr> <tr> <td>Define New</td><td>display the Analysis Control dialog box. For information on completing the fields, see "About the Analysis Parameters" on page 5-2.</td></tr> </tbody> </table>	Choose...	To...	Analysis Parameters	apply the parameters that are stored as preferences in the software.	Define New	display the Analysis Control dialog box. For information on completing the fields, see "About the Analysis Parameters" on page 5-2.
Choose...	To...						
Analysis Parameters	apply the parameters that are stored as preferences in the software.						
Define New	display the Analysis Control dialog box. For information on completing the fields, see "About the Analysis Parameters" on page 5-2.						
4	<p>Repeat step 2 and step 3 for each sample.</p> <p>Note You can also use the Cut, Copy, and Paste commands from the Edit menu.</p>						


Displaying Default Parameters

There are two ways to display the default parameters:

You can...	Then...
select Analysis Parameters from the Settings menu.	the Analysis Parameters dialog box opens with the default parameters. For information on defining the Analysis Parameters, refer to page 5-2.
double-click Analysis Parameters from any entry in the Parameters column in the Analysis Control window.	


Creating Custom Analysis Parameter Files

To create a custom analysis parameter file:

Step	Action
1	<p>Select New from the File menu.</p> <p>Note You can also select the Define New option from the pop-up menu in the Parameters column of the Analysis Control window.</p> <p>The Create New dialog box opens.</p> 
2	<p>Click the Analysis Parameters icon.</p> <p>The Analysis Parameters dialog box opens.</p>
3	<p>Change the parameters as necessary.</p> <p>For more information on the analysis parameters, refer to page 5-2.</p>
4	<p>Choose Save (Ctrl+S) from the File menu.</p> <p>A dialog box opens.</p>
5	<p>Enter a descriptive name and click Save.</p> <p>The file now opens in the pop-up menu for analysis parameters in the GeneScan Analysis Control window.</p> <p>You can also select the file in the Data Collection software for automatic analysis.</p>

Changing an Existing Analysis Parameters File

To change an existing analysis parameter file:

Step	Action				
1	<p>Select Open (Ctrl+O) from the File menu.</p> <p>Note You can also double-click the file name in the Parameters column of the Analysis Control window.</p> <p>The Open Existing dialog box opens.</p> 				
2	<p>Click the Analysis Parameters icon.</p> <p>An Open directory dialog box opens.</p> <table border="1"> <thead> <tr> <th>If...</th><th>Then...</th></tr> </thead> <tbody> <tr> <td>the Params folder does not appear</td><td> <p>navigate to the folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params</p> </td></tr> </tbody> </table>	If...	Then...	the Params folder does not appear	<p>navigate to the folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params</p>
If...	Then...				
the Params folder does not appear	<p>navigate to the folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params</p>				
3	Select a file that you want to change and click Open .				
4	Make the changes and close the window by clicking OK .				

Deleting Custom Analysis Parameters

To delete a custom analysis parameters file:

Step	Action
1	Click the Start button, and then point to Programs .
2	<p>Click Windows NT® Explorer and find the file in the folder you specified as the location for the analysis parameters files.</p> <p>Note The custom analysis parameters file is in the folder that is normally called the GS Parameters Folder. This folder is inside the ABI PRISM GeneScan folder.</p>
3	Drag the custom analysis parameters file to the Recycle Bin.

Making a Matrix File

Overview

This chapter describes the processes of making a matrix file, loading and running samples and evaluating the file. You have the option of selecting four or five dyes depending on the application when creating a new matrix for data collection.

Note This chapter applies only to the 310 and 377 instruments.

In This Chapter Topics in this chapter include the following:

Topics	See page
About Matrix Files	6-2
Process of Creating a New Matrix File	6-7
Loading and Running Dye Standards for the ABI PRISM 310	6-9
Loading and Running Dye Standards for the ABI PRISM 377	6-12
Generating Matrix Sample Files for the ABI PRISM 377 Instrument	6-15
Choosing a Scan Range for the Matrix Calculation	6-17
Generating a New Matrix File	6-20
Saving and Naming the Matrix File	6-22
Assigning the Matrix File to Sample Files	6-23
Evaluating the Matrix File	6-25
Causes for Bad Matrix Files	6-26

About Matrix Files

Introduction There are three dye-labeling chemistries currently available to prepare nucleic acid samples to use the GeneScan® Analysis Software on ABI PRISM® instruments:

- ◆ Fluorescent NHS-Ester
- ◆ Fluorescent dNTP
- ◆ Fluorescent Phosphoramidite

Each chemistry has a set of dye labels that fluoresce at different wavelengths when excited by a laser.

During data collection on the...	The wavelengths are separated...
310 or 377, 377XL, or 96-lane upgrade instrument	by a spectrograph into a known spectral pattern across a detection system with the sequencer.

Matrix File Definition Matrix files are mathematical matrices that correct for spectral overlap of fluorescent emission spectra data collected from ABI PRISM® instruments.


A matrix file allows you to account for spectral overlap when analyzing Sample files.

Multicomponent Definition This process of eliminating the bleed-through caused by spectral overlaps is called multicomponenting.

Applying a matrix file to raw data allows you to generate multicomponented data.

Why Is a Matrix File Necessary

A matrix file is necessary because the four or five dyes used to label the fragments fluoresce at different wavelengths and may have spectral overlaps, as shown below:



When to Create a Matrix File


Create a matrix file for each dye set used from that particular instrument before analyzing fragment data.

You may have to create new matrix files for different gel compositions or unusual run conditions.


Sample Files Using Matrix File

The figures below show examples of data analyzed with and without a matrix file.

You can see that peak data from a Sample file analyzed without a matrix file displays the expected peak, along with extra peaks in other dye colors, or bleed-through from other dye colors.



Sample file analyzed without a matrix file



Sample file analyzed with a matrix file

Installing a Matrix File

You can install a matrix file into a gel file or into a Sample file.

Normally the matrix from a matrix file is installed within a gel file or a Sample file automatically upon generation during or after a run. Additionally a matrix can be manually installed into a 377/310 Sample file from within the GeneScan Analysis Software.

When to Assign a Matrix File

Before you can successfully analyze Sample files using the GeneScan Analysis Software, you must make a new matrix file or assign an existing one to a set of Sample files.

Limitations to Matrix Files

You can only assign a matrix file to Sample files generated on the same instrument, under the same electrophoresis, gel matrix and buffer conditions, and using the same dye set.

Note If you are using a fifth dye, then you need to create a new matrix file for that dye.

**When to Create a
New Matrix File**

Create a new matrix file in the following conditions:

- ◆ For each dye set:
 - NHS-Esters
 - Phosphoramidite set
 - Fluorescent dNTPs
 - ◆ Whenever you change the dye set you use to label sample fragments, for example, if you are using the fifth dye.
 - ◆ When you use gel materials or buffers with pH values that differ greatly from the pH value of the gel material or buffer on which the existing matrix files were generated.
 - ◆ When you use dyes other than those provided by Applied Biosystems.
 - ◆ When you run the same gel on a different instrument.
 - ◆ When you see multiple unexpected peaks of different colors under an expected peak.
 - ◆ When you recalibrate your CCD camera (310 and 377 instruments) and the change is greater than 3 pixels from the original pixel position.
 - ◆ When you replace the CCD camera (ABI PRISM 310 and ABI PRISM 377).
-

Considerations Before Making a Matrix File

The following table lists some of the considerations before making a matrix file:

Consideration	Comment
How much dye matrix standard to load?	With the ABI PRISM 377, loading more than 3 μ L, produces too much signal. Any amount that results in a signal over 4,000 FUs is too strong.
Which lanes to load with the dye matrix standards?	For gel electrophoresis, load the matrix standards with an empty lane between each sample to avoid contamination of the individual dyes by residual material leaking adjacent samples.
What exact gel data will be used for matrix creation?	After generating a gel image, for ABI PRISM 377 instrument, check that the tracking of the gel file is adequate.

Where to Store Matrix Files


Store matrix files intended for use by Data Collection software in:

D:\AppliedBio\Shared\Analysis\Sizecaller\Matrix\

If Data Collection and Analysis are installed on different computers, the location is the same. Remember to copy the matrix file from the analysis computer to the Data Collection computer.

Process of Creating a New Matrix File

Process Diagram The following diagram shows the procedure for making a new matrix file:



For sample preparation and loading information, refer to the appropriate instrument user manual

Steps to Create a New Matrix

The following table lists the steps to create a new matrix file:

Step	Process	See Page
1	Loading and Running Dye Standards for the ABI PRISM 310	6-9
2	Loading and Running Dye Standards for the ABI PRISM 377	6-12
3	Choosing a Scan Range for the Matrix Calculation	6-17
4	Generating a New Matrix File	6-20

The following table lists the steps to create a new matrix file: *(continued)*

Step	Process	See Page
5	Saving and Naming the Matrix File	6-22
6	Assigning the Matrix File to Sample Files	6-23
7	Evaluating the Matrix File	6-25
8	Causes for Bad Matrix Files	6-26

Loading and Running Dye Standards for the ABI PRISM 310


Introduction This section describes how to do the following tasks before running matrix standards on the ABI PRISM 310:

Topic	See Page
Creating a GeneScan Sample Sheet	6-9
Creating a GeneScan Injection List	6-10
Starting the Data Collection Software	6-11

Note When loading the matrix standards on an instrument, note which colors you load in which autosampler positions.

Creating a GeneScan Sample Sheet

To create a GeneScan Sample Sheet:

Step	Action
1	Open the ABI PRISM 310 Data Collection software.
2	Choose New from the File menu.
3	Click the GeneScan Sample Sheet 48 or 96 tube icon. The GeneScan Sample Sheet opens. 
4	Enter the appropriate colors (Blue, Green, Yellow, Red, or Orange) in the positions for each sample name.
5	Enter any additional information about the sample in the Comments column.


To create a GeneScan Sample Sheet: *(continued)*

Step	Action
6	Use the Save As command and save the Sample Sheet to the Sample Sheet folder.

Creating a GeneScan Injection List

The Sample Sheet is imported into the Injection List, which defines the sample names and the initial injection order.

To create a GeneScan Injection List:

Step	Action
1	Choose New from the File menu.
2	Click the GeneScan Injection List icon. The GeneScan Injection List opens.
	
3	Select a Sample Sheet from the Sample Sheet pop-up menu.
4	Select the appropriate module in the Module pop-up menu for lines 1 through 4, for example, A1, A3, A5, and A7.
5	Unselect Auto-Anlz (analysis).
6	Unselect Auto-Print .

Starting the Data Collection Software

To start the Data Collection software:

Step	Action	
1	Click the Run button.	
	The following windows open:	
	Window	Description
	Raw Data window	Shows the real-time chromatogram of the run.
	Log Window	Shows the real-time written record of run events.
2	Choose Status from the Window menu.	
	The current run is in italics in the Injection List.	
	You can monitor activities such as electrophoresis current, laser power, running time, and gel temperature.	

If the Run was Cancelled

If the run was cancelled and you are using Data Collection Software Version 1.0.4 or later (Macintosh or Windows version), the sample file is saved if you skip to the next sample or cancel a run.

Run Time

Run time is approximately 30 minutes for the GS STR POP-4 module, so the total run time will be about 120 minutes.

Assigning Matrix File

For information on “Assigning the Matrix File to Sample Files,” refer to page 6-23.

Loading and Running Dye Standards for the ABI PRISM 377


Introduction This section describes how to do the following tasks using the 377 and the 377 with XL Upgrade instruments:

Topic	See page
Creating a GeneScan Sample Sheet	6-12
Loading Matrix Standards	6-13
Running the Matrix Standards	6-14

Note When loading the matrix standards on an instrument, note which colors you load in which lanes for gel-based systems.

Creating a GeneScan Sample Sheet The GeneScan Sample Sheet assigns sample and dye information to the appropriate lane.
To create a GeneScan Sample Sheet:

Step	Action
1	Open the ABI PRISM 377 Data Collection software. Define or verify the data collection preferences.
2	Choose New from the File menu.
3	Click the GeneScan Sample Sheet icon. The GeneScan Sample Sheet opens.



To create a GeneScan Sample Sheet: *(continued)*

Step	Action
4	Enter the individual colors in the appropriate lanes where the matrix standards are loaded. Note It is important to fill out the Sample Sheet completely.
5	Enter any additional information about the sample in the Comments column.
6	Use the Save As command and save the Sample Sheet to the Sample Sheet folder.

Loading Matrix Standards


To load matrix standards:

Step	Action									
1	For denaturing gels load: <ul style="list-style-type: none">◆ 0.5–2 µL of matrix standard per lane.◆ 8 lanes with different colors, leaving an empty lane between each lane of matrix standard.									
2	<p>Complete the information in the data collection Run sheet, making sure to choose the appropriate PreRun and Run modules.</p> <p>Take the following action:</p> <table><tr><th>For this matrix...</th><th>Choose modules that use...</th><th>Module file</th></tr><tr><td>Dye Primer matrix</td><td>Virtual Filter A</td><td>GS 36A...</td></tr><tr><td>Fluorescent Amidite</td><td>Virtual Filter C</td><td>GS36C...</td></tr></table>	For this matrix...	Choose modules that use...	Module file	Dye Primer matrix	Virtual Filter A	GS 36A...	Fluorescent Amidite	Virtual Filter C	GS36C...
For this matrix...	Choose modules that use...	Module file								
Dye Primer matrix	Virtual Filter A	GS 36A...								
Fluorescent Amidite	Virtual Filter C	GS36C...								
3	Electrophorese samples according to conditions specified in your instrument manual.									

Running the Matrix Standards

Run the matrix standards under the precise conditions you want to generate a matrix file.

To run the matrix standards:

Step	Action									
1	<p>Complete the information in the data collection Run sheet.</p> <p>Refer to the figure below to select the appropriate PreRun and Run Modules.</p> 									
2	<p>Take the following action:</p> <table><tr><th>For this matrix...</th><th>Choose modules that use...</th><th>Module file</th></tr><tr><td>Dye Primer or dNTP matrix</td><td>Virtual Filter A</td><td>GS 36A...</td></tr><tr><td>GS Amidite matrix</td><td>Virtual Filter C</td><td>GS36C...</td></tr></table>	For this matrix...	Choose modules that use...	Module file	Dye Primer or dNTP matrix	Virtual Filter A	GS 36A...	GS Amidite matrix	Virtual Filter C	GS36C...
For this matrix...	Choose modules that use...	Module file								
Dye Primer or dNTP matrix	Virtual Filter A	GS 36A...								
GS Amidite matrix	Virtual Filter C	GS36C...								
3	Start the electrophoresis run according to the conditions specified in your instruction manual.									
4	Go to “Generating Matrix Sample Files for the ABI PRISM 377 Instrument” on page 6-15.									

Generating Matrix Sample Files for the ABI PRISM 377 Instrument

Who Should Use This Step Follow these guidelines to use this step:

Who should use this step	Who should not use this step
<p>This step is only necessary if you are using an:</p> <ul style="list-style-type: none">◆ ABI PRISM 377 instrument◆ XL-upgraded instrument◆ 96-lane upgrade	<p>The ABI PRISM 310 Genetic Analyzer automatically processes collection data and generates Sample files when the run completes.</p> <p>If you ran your Dye Matrix Standards on the ABI PRISM 310, go to “Choosing a Scan Range for the Matrix Calculation” on page 6-17.</p>

Gel Handling on Windows The 377 Data Collection for Windows generates gel files. There is a new program on Windows called Gel Processor. This program handles all gel file-related tasks, including tracking and sample file generation. For more information on the Gel Processor, refer to the Gel Processor user's manual.

Verify Tracking Before generating Sample files, verify that the lanes were set and tracked correctly. Once you have successfully completed your run of Dye Matrix Standards, open the gel file in the Gel Processor and visually verify the positions of the tracker lines.

Generating Sample Files To generate Sample files from within Gel Processor:

Step	Action
1	<p>Click the appropriate Lane Indicators at the top of the Gel window to make sure the tracker line is optimally aligned over each band in all lanes.</p> <p>If any of the lanes are not properly tracked, use the tracker line editing tools to align the tracker lines in each lane. Refer to the <i>Gel Processor User's Manual</i>.</p>
2	<p>If you change any of the lane assignments in any way, save the changes.</p>

To generate Sample files from within Gel Processor: *(continued)*

Step	Action
3	<p>From the Gel menu, choose Track & Extract Lanes.</p> <p>Note If the gel is already tracked, choose Extract Lanes from the Gel menu.</p> <p>The project Analysis Control window opens containing each of the Sample files.</p>

Choosing a Scan Range for the Matrix Calculation

Introduction Depending on how well your Matrix Standards run, it may be necessary for you to choose a specific range of data points to be considered for your matrix calculation.

In order to choose appropriate values for the Scan range, you must first view the Sample file raw data from each of the matrix standard files, so you can decide where to choose the start and stop points for the scan range.

Viewing the Raw Data When you have multiple Sample files, raw data can be accessed more easily through a project's Analysis Control window. Raw data can provide useful information about the Sample files you have created.

Note You can view Sample files without opening a project. However, this procedure is easier if you use a project to organize the Sample files.

To view raw data:

Step	Action
1	<p>Use the following steps to create a project for the Dye Matrix Standards:</p> <ol style="list-style-type: none">Choose New from the File menu. The Create New dialog box opens.Click the Project icon. An untitled Analysis Control window opens.Choose Add Sample Files from the Project menu.Find and open your matrix run folder.Select the five Sample files representing the blue, green, yellow, red, and orange dye-labeled runs, and click Add.Click Done after the Sample files are transferred. <p>For more information, refer to "Creating a New Project" on page 2-9.</p>
2	<p>From the Analysis Control window, select the matrix standard Sample files by clicking on the first Sample file, holding down the mouse button, and releasing on the last Sample file.</p>

To view raw data: *(continued)*

Step	Action
3	<p>From the Project menu, choose Raw Data (Ctrl+R).</p> <p>Electropherograms displaying raw data from the matrix standard Sample files appear.</p> <p>Note For the ABI PRISM 377 instrument, you can also view raw data from the gel display by selecting one of the lanes containing Dye Matrix Standard and looking at the Slice View to the left of the gel image.</p> <p>For more information about viewing the gel image, refer to the <i>Gel Processor User's Manual</i>.</p>

What to Look For in the Raw Data Display

In the raw data display of the Sample files verify the following:

- ◆ Data peaks are present in all four of the matrix standards.
- ◆ There are no anomalies.
- ◆ The baseline is stable.
- ◆ Peaks should be on-scale—no more than 4000 relative fluorescent units—and the peaks of the dye of interest should have a value of at least 200.

If peak data does not show these characteristics, refer to “Causes for Bad Matrix Files” on page 6-26, for possible interpretations of your peak data.

Choosing a Scan Range

To choose a scan range:

Step	Action
1	<p>Move the cursor well away from the primer peak, in a region at the beginning of the run and in a flat part of the baseline, and record the scan numbers.</p> <p>Note When choosing the start point, do not include primer peaks in the scan range (refer to “Eliminating Primer Peaks”). Also, the region for both the start and stop points should be flat points at the baseline.</p>
2	<p>Record the data point values for both the start and stop points of the scan range.</p> <p>You will need to enter these values in the next step when generating the new matrix file (refer to page 6-20).</p>

To choose a scan range: *(continued)*

Step	Action
3	<p>Close the raw data boxes and the project by clicking in the upper left-hand corner of the window.</p> <p>Note Holding down the Alt key while clicking in the upper left-hand corner of the window will close the windows simultaneously.</p>


Eliminating Primer Peaks

Both the primer peaks and the data peaks are displayed when viewing the raw data of your matrix standards. Any time you run dye-labeled samples on a gel (377 instrument), or capillary (310 instrument), you always have excess dye-labeled primer in the reaction. The primer peak displays as the first peak, usually off-scale because it is in molar excess.

Eliminate the primer peak when making a matrix, by choosing the start point after the primer peak in a flat area with a stable baseline.

Generating a New Matrix File

Procedure To generate a new matrix file:


Step	Action
1	Choose New from the File menu. The Create New dialog box opens.
2	Click the Matrix icon. The Make New Matrix dialog box opens. 
3	Choose the number of dyes from the Number of Dyes pop-up menu. If 5 dyes are selected, a button is added to the bottom of the list.
4	The B, G, Y, R, and O buttons represent dye colors. a. Click a button to display a pop-up menu. b. Use the pop-up menu to access a Sample file to link to each of the dye-labeled primers. c. Choose the Sample file that represents the dye color for that button.
5	Enter the start point that you determined when choosing a scan range in the Start at field. Refer to “Choosing a Scan Range” on page 6-18.

To generate a new matrix file: *(continued)*

Step	Action
6	Enter the total number of data points to include to calculate the matrix in the Points field. Note You must have at least five peaks to make a matrix. In most cases, leave the default value, unless you must exclude a portion of your data because of artifacts or bleed-through.
7	Click OK . This generates a new matrix file.

Matrix File Example

The following is an example of the Matrix Values window that opens showing the values used to calculate the overlap correction.



	B	G	Y	R	O
B	1.0000	1.4527	0.5962	0.6376	0.1484
G	0.8133	1.0000	1.5131	1.5981	2.0000
Y	0.6811	0.7925	1.0000	1.4584	2.0000
R	0.4637	0.5135	0.8426	1.0000	2.0000
O	0.0596	0.0746	0.1621	0.1728	1.0000

For each dye, the value where the dye fluorescence is read by the appropriate filter is 1.000. The adjacent colors show the amount of overlap for which the system must compensate. The adjacent values, in most cases, should be less than 1.000, but equal to or greater than 0.0000.

Saving and Naming the Matrix File

Introduction The matrix file is instrument-specific. You cannot apply a matrix file you made on the ABI PRISM 377 to data you collected on an ABI PRISM 310, nor can you apply a matrix file made on an ABI PRISM 377 to a Sample or gel file made on another ABI PRISM 377. In other words, you cannot apply matrix files created on one instrument to other instruments of the same model.

Naming Considerations When naming a matrix file, consider including the following information in the name:

Item	For example
Instrument type	ABI PRISM 377, or ABI PRISM 310.
Filter set used	D, C, F, G5 and E5
Gel conditions	native or denaturing.

Saving the Matrix File To save the matrix file:

Step	Action
1	Select Save from the File menu. A dialog box opens.
2	Enter a descriptive name for the new matrix file and click Save .

Where to Store the Matrix File Store matrix files that are intended for use by data collection to assign to collection runs in:

D:\AppliedBio\Shared\Analysis\Sizecaller\Matrix\


When Data Collection and Analysis Are Performed on Different Computers Be sure to copy matrix files generated on the analysis computer (the computer running GeneScan analysis) to the Data Collection computer, if different. This will ensure that the correct matrix is stored in the GeneScan sample file. The proper matrix is required for accurate analysis of 310 and 377 sample files.

Assigning the Matrix File to Sample Files


Introduction After generating the new matrix file, assign it to all the Sample files that you want to analyze.

IMPORTANT After assigning your matrix file to Sample files, refer to “Evaluating the Matrix File” on page 6-25.

Procedure To assign a matrix file to Sample files and review them:

Step	Action
Assigning a matrix file to Sample files:	
1	From the project that contains your matrix standard Sample files, open the Analysis Control window. 
2	Select the matrix standard Sample files.
3	Select Assign New Matrix from the Project menu.
4	Select the matrix file you just created.
5	Select numbers 1 , 2 , 3 , and 4 on the left side of the window to highlight the colors for each row.
6	Select Set Analysis Parameters from the Settings menu.
7	Enter the appropriate range for the Analysis Range and click OK to return to the Analysis Control window.
8	Click Analyze .

To assign a matrix file to Sample files and review them: *(continued)*

Step	Action
Reviewing the results:	
1	<p>Choose Results Control from the Windows menu.</p> <p>The Results Control window opens as shown below.</p>  <p style="text-align: right;">Dye/Samples</p> <p style="text-align: right;"># of Panels pull-down menu</p>
2	Select 4 from the # of Panels pull-down menu (see above).
3	Click 1 under Dye/Samples (see above).
4	Click 1 on the Sample Files side of the Results Control window.
5	For the rest of the matrices, click Analyze and repeat steps 1 and 2.

Evaluating the Matrix File

Introduction After creating a new matrix file and assigning it to select Sample files, the next step is evaluating the quality of the matrix file. The quality of the matrix file has a direct impact on the quality of the results data.

Procedure To evaluate the matrix file:

Step	Action	
1	Analyze the Sample files used to make the matrix.	
2	Display Results data for all the Dye Matrix Standard Sample files on one screen, showing only electropherogram data.	
3	For each displayed Sample file:	
	You should see...	If not...
	that the only visible peaks represent the color of the Dye Matrix Standard run in that lane, or for that injection (ABI PRISM 310).	you probably have a bad matrix file. For instructions on how to identify and correct problems with bad matrix files see, "Causes for Bad Matrix Files" on page 6-26.
	all other lines should be relatively flat. This indicates that the matrix properly compensated for the spectral overlap. For example, for the blue matrix standard Sample file, you should only see blue.	
	sharp, well-defined, singularly colored peak data.	

Causes for Bad Matrix Files


If an Error Message Appears

There are two possible causes for the error messages shown in the following table:


For this cause...	Take this action...
designated the wrong files.	reassign the matrix files. Refer to “Assigning the Matrix File to Sample Files” on page 6-23.
signal is too weak to make a matrix.	rerun the matrix standards.

Two Causes of Bad Matrix Files

The following table lists two common causes of bad matrix files:

Problem	Cause	What to do
Artifact peaks of different colors under the true peaks. Refer to the figure below.	Loading too much dye when running matrix standards, resulting in dye bleed-through.	Complete another run and recreate the matrix.
		
Noisy baseline.	If the matrix subtracts too much of a particular color from the sample, then the baseline may become too elevated, resulting in false peaks. Refer to the figure below.	Complete another matrix run and make sure you do not have any off-scale data.

The following table lists two common causes of bad matrix files: *(continued)*



Working with Size Standards

7

Overview

In This Chapter Topics in this chapter include the following:

Topics	See Page
About Size Standards	7-2
Defining the Size Standard	7-3
Using Size Standards	7-9

About Size Standards

What Are Size Standards	Size standards are specific DNA fragments of known sizes. After defining the peaks of a size standard, the GeneScan® Analysis Software matches this definition to the internal size standard included with the run. The software assigns the defined size values to the appropriate peaks of the internal size standard, and uses this information with the selected sizecalling method to size all unknown fragments.
--------------------------------	--

Advantages of Using Size Standard	Running an internal size standard results in determining accurate and precise molecular length. This is because the internal size standard and the unknown fragments undergo exactly the same electrophoretic forces. The GeneScan Analysis Software can then compensate for band-shift artifacts caused by variations in the gel and in the sample from lane to lane or injection to injection.
--	--

Size Standards Provided	Applied Biosystems provides several fluorescently labeled size standards, which are described in Appendix C.
--------------------------------	--

You can also label and use other fragments if they better suit the fragment sizes with which you are working.

When to Define Size Standards	Normally, a size standard is defined using the GeneScan Analysis Software after running the size standard with samples on the instrument. The software detects peaks for a selected dye color in a selected sample file and allows you to define the peak sizes. You can save the defined size standard in a file and use it to automatically analyze other samples run with the same size standard and under the same conditions.
--------------------------------------	--

For a list of the size standards that are supplied with the GeneScan Analysis Software, refer to Appendix C, “GeneScan Size Standards.”

If Split Peaks Appear	Split peaks might appear in size standards in which both DNA strands are labeled. For some peaks, the two strands migrate at different rates when running under denaturing conditions, and they appear as two peaks approximately half the height of normal, non-split, peaks. One peak of the two runs is true to size. Assign a size to that peak for the size standard definition, and assign zero to the other peak.
------------------------------	--

Defining the Size Standard



Two Ways to Define the Size Standard

There are two ways to define a new size standard:


Topic	See Page
Using the New Command	7-3
Using the Analysis Control Window	7-7

Using the New Command


To use the New command to define a new size standard:

Step	Action
1	<p>Select New from the File menu.</p> <p>The Create New dialog box opens.</p> <div>The 'Create New' dialog box is a light gray window. It contains the text 'Create New:' at the top. Below this text are four icons: 'Project' (a folder with a document), 'Analysis Parameters' (a document with a red star), 'Size Standard' (a document with a yellow bar chart), and 'Matrix' (a document with a red grid). To the right of these icons is a 'Cancel' button.</div>
2	<p>Click the Size Standard icon.</p> <p>The following dialog box opens:</p> <div>The 'Select Sample File' dialog box is a blue window. It has a title bar with a question mark and a close button. The main text says 'Select a Sample File to use as a Template for the New Size Standard Definition'. Below this is a 'Look in:' dropdown menu showing 'AppliedBio'. To the right of the dropdown are three icons: a folder, a document, and a grid. Below the dropdown is a list box containing 'GelProcessor'. At the bottom, there is a 'File name:' text box, a 'Files of type:' dropdown menu showing 'Sample Files (*.FSA)', and two buttons: 'Open' and 'Cancel'.</div>

To use the New command to define a new size standard: *(continued)*

Step	Action						
3	<p>Select the sample file that contains the dye standard you want to use as a template for the new size standard and click Open.</p> <p>The Select Dye and Analysis Parameters dialog box opens.</p> 						
4	<p>From the Dye pop-up menu, select the code that represents the dye label of the size standard in the selected sample File.</p>						
5	<p>From the Analysis Parameters pop-up menu, select the analysis parameters to use.</p> <p>The pop-up menu contains the following options:</p> <table border="1"> <thead> <tr> <th>Item</th><th>Description</th></tr> </thead> <tbody> <tr> <td>Analysis Parameters</td><td>Applies the parameters that are stored as preferences in the software.</td></tr> <tr> <td>custom parameters that are listed at the bottom of the menu</td><td> <p>These are files that you defined and they are located in the Params folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params.</p> </td></tr> </tbody> </table>	Item	Description	Analysis Parameters	Applies the parameters that are stored as preferences in the software.	custom parameters that are listed at the bottom of the menu	<p>These are files that you defined and they are located in the Params folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params.</p>
Item	Description						
Analysis Parameters	Applies the parameters that are stored as preferences in the software.						
custom parameters that are listed at the bottom of the menu	<p>These are files that you defined and they are located in the Params folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params.</p>						

To use the New command to define a new size standard: *(continued)*

Step	Action																												
6	<p>Click OK.</p> <p>A window opens (see below) showing the electropherogram and a table of peaks for the dye color and sample selected.</p> <p>You should be able to recognize the peak pattern of the size standard in the electropherogram.</p> <p>Note You can only change the peak size value in the right column of the table. You cannot change or rearrange the peak numbers.</p> <table><tr><th>If...</th><th>Then...</th></tr><tr><td>too many peaks appear in the electropherogram or the baseline is too high</td><td>you may need to adjust the analysis parameters. See “Using Analysis Parameter Files” on page 5-13.</td></tr></table> <p>The software assigns a number to each peak found in the electropherogram in order, from left to right.</p>  <p>The screenshot shows a software window titled "GS 350 2400 36 STD.szs". It contains an electropherogram plot with a red line on a black background. The x-axis is labeled from 1400 to 4400 in increments of 200. The y-axis is labeled from 0 to 152 in increments of 38. Below the plot is a table with three columns: "Peak", "Data Points", and "Size". The table contains seven rows of data. To the right of the table is a section titled "Template File Info:" which includes fields for "File:", "Run Date:", "Run Time:", "Dye:", and "Parameters:".</p> <table><tr><th>Peak</th><th>Data Points</th><th>Size</th></tr><tr><td>1</td><td>1470</td><td>75</td></tr><tr><td>2</td><td>1691</td><td>100</td></tr><tr><td>3</td><td>2067</td><td>139</td></tr><tr><td>4</td><td>2167</td><td>150</td></tr><tr><td>5</td><td>2273</td><td>160</td></tr><tr><td>6</td><td>2702</td><td>200</td></tr><tr><td>7</td><td>3253</td><td>250</td></tr></table> <p>Template File Info: File: 01*1 Run Date: Wed, Jul 31, 1996 Run Time: 3:45:39 PM Dye: R Parameters: 2400 36 STD AP</p>	If...	Then...	too many peaks appear in the electropherogram or the baseline is too high	you may need to adjust the analysis parameters. See “Using Analysis Parameter Files” on page 5-13.	Peak	Data Points	Size	1	1470	75	2	1691	100	3	2067	139	4	2167	150	5	2273	160	6	2702	200	7	3253	250
If...	Then...																												
too many peaks appear in the electropherogram or the baseline is too high	you may need to adjust the analysis parameters. See “Using Analysis Parameter Files” on page 5-13.																												
Peak	Data Points	Size																											
1	1470	75																											
2	1691	100																											
3	2067	139																											
4	2167	150																											
5	2273	160																											
6	2702	200																											
7	3253	250																											

To use the New command to define a new size standard: *(continued)*

Step	Action								
7	<p>Specify the peaks of the size standards and their sizes as follows:</p> <table><tr><th>Step</th><th>Action</th></tr><tr><td>a.</td><td><p>Click the peak you want to define either in the electropherogram or in the table.</p><p>Use the Zoom In (Ctrl+ Plus sign) and Zoom Out (Ctrl+ Minus sign) commands from the View menu to zoom the electropherogram for easier viewing.</p><ul style="list-style-type: none">– If you click a peak in the electropherogram, the corresponding row in the table is highlighted.– If you click a peak number in the table, the corresponding peak in the electropherogram is highlighted.</td></tr><tr><td>b.</td><td><p>Type the value for the selected peak in the corresponding Size field in the table.</p><p>Refer to Appendix C, “GeneScan Size Standards,” for values and peaks patterns.</p><p>Note Leave a zero in the Size field to ignore a peak for the size standard definition.</p></td></tr><tr><td>c.</td><td><p>Press Enter to automatically move to the next size standard peak.</p></td></tr></table>	Step	Action	a.	<p>Click the peak you want to define either in the electropherogram or in the table.</p> <p>Use the Zoom In (Ctrl+ Plus sign) and Zoom Out (Ctrl+ Minus sign) commands from the View menu to zoom the electropherogram for easier viewing.</p> <ul style="list-style-type: none">– If you click a peak in the electropherogram, the corresponding row in the table is highlighted.– If you click a peak number in the table, the corresponding peak in the electropherogram is highlighted.	b.	<p>Type the value for the selected peak in the corresponding Size field in the table.</p> <p>Refer to Appendix C, “GeneScan Size Standards,” for values and peaks patterns.</p> <p>Note Leave a zero in the Size field to ignore a peak for the size standard definition.</p>	c.	<p>Press Enter to automatically move to the next size standard peak.</p>
Step	Action								
a.	<p>Click the peak you want to define either in the electropherogram or in the table.</p> <p>Use the Zoom In (Ctrl+ Plus sign) and Zoom Out (Ctrl+ Minus sign) commands from the View menu to zoom the electropherogram for easier viewing.</p> <ul style="list-style-type: none">– If you click a peak in the electropherogram, the corresponding row in the table is highlighted.– If you click a peak number in the table, the corresponding peak in the electropherogram is highlighted.								
b.	<p>Type the value for the selected peak in the corresponding Size field in the table.</p> <p>Refer to Appendix C, “GeneScan Size Standards,” for values and peaks patterns.</p> <p>Note Leave a zero in the Size field to ignore a peak for the size standard definition.</p>								
c.	<p>Press Enter to automatically move to the next size standard peak.</p>								
8	<p>When you finish defining the peaks, save the size standard by selecting Save As from the File menu.</p> <p>Note You can also click the Close button.</p>								
9	<p>Enter a descriptive name for the size standard and click Save.</p> <p>Note Run conditions are not stored in the Size Standard file. Use a name that clearly defines the size standard for future use.</p> <p>This file is automatically saved in the SizeStandards folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards</p>								

Using the Analysis Control Window

To define a new size standard:

Step	Action								
1	Open an existing project or create a new project.								
	<table><tr><th>For information on...</th><th>See Page</th></tr><tr><td>Opening an Existing Project</td><td>2-8</td></tr><tr><td>Creating a New Project</td><td>2-9</td></tr></table>	For information on...	See Page	Opening an Existing Project	2-8	Creating a New Project	2-9		
	For information on...	See Page							
	Opening an Existing Project	2-8							
Creating a New Project	2-9								
The Analysis Control window should open. If it does not, select Analysis Control (Ctrl+1) from the Windows menu.									
2	In the Analysis Control window, find the row that contains the sample for which you want to define the size standard.								
3	In that row, Ctrl+click the dye color cell that represents your size standard.								
	A diamond symbol (◆) appears in the cell, identifying it as the size standard.								
4	Click the arrow in the Parameters field of the same row, and select an option from the pop-up menu that opens.								
	The pop-up menu contains the following options:								
	<table><tr><th>Item</th><th>Description</th></tr><tr><td>Analysis Parameters</td><td>Applies the parameters that are stored as preferences in the software.</td></tr><tr><td>filename</td><td>Applies the settings specified in the data collection run file.</td></tr><tr><td>Custom parameters that are listed at the bottom of the menu</td><td>These are files that you defined and they are located in the Params folder. The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params</td></tr></table>	Item	Description	Analysis Parameters	Applies the parameters that are stored as preferences in the software.	filename	Applies the settings specified in the data collection run file.	Custom parameters that are listed at the bottom of the menu	These are files that you defined and they are located in the Params folder. The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params
	Item	Description							
Analysis Parameters	Applies the parameters that are stored as preferences in the software.								
filename	Applies the settings specified in the data collection run file.								
Custom parameters that are listed at the bottom of the menu	These are files that you defined and they are located in the Params folder. The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params								
5	In the same row, click the arrow in the Size Standard field, and select Define New from the pop-up menu.								
	A window opens showing the electropherogram and a table of peaks for the dye color and sample you selected.								
	You should be able to recognize the peak pattern of the size standard in the electropherogram.								

To define a new size standard: *(continued)*

Step	Action
6	Follow step 7 to step 9 on page 7-6. The name of the size standard appears in the Size Standard pop-up menu in the Analysis Control window.

Using Size Standards

In This Section This section contains the following topics:

Topic	See Page
Changing the Number of Peaks Detected	7-9
Editing the Size Standard Definition	7-9
Using the Open Command to Edit an Existing Size Standard	7-10
Using the Analysis Control Window to Edit an Existing Size Standard	7-11
Deleting an Existing Size Standard	7-12
Analyzing Samples Using the Same Size Standard	7-13
Selecting Separate Size Standards for Samples	7-14

Verifying Size Calculations For information, see “Verifying Size Calculations” on page 9-33.

Changing the Number of Peaks Detected Use the Analysis Parameters dialog box to change the number of peaks detected in the Define New Standard window.

For information, refer to “Peak Detection Options” on page 5-8.

Editing the Size Standard Definition The following procedure describes how to edit the size standard definition (Collection Setting standard) that is embedded in the sample file. Double-clicking the size standard definition will not open the file.

To edit the size standard definition:

Step	Action
1	Choose Open from the File menu. The Open Existing dialog box opens. The Collection Setting does not change.
2	Click the Size Standard icon. An Open dialog box opens.

To edit the size standard definition: *(continued)*


Step	Action
3	Navigate to the SizeStandards folder, select the file that you want to edit, and click Open . The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards
4	Make any changes to the file and click the Close button.
5	Choose the edited file from the Size Standard pop-up menu to apply the changes.

Editing an Existing Size Standard


The two ways to edit a previously defined size standard are by:

- ◆ Using the Open Command to Edit an Existing Size Standard
- ◆ Using the Analysis Control Window to Edit an Existing Size Standard

Using the Open Command to Edit an Existing Size Standard

Step	Action
1	Choose Open (Ctrl+O) from the File menu. The Open Existing dialog box opens. 
2	Click the Size Standard icon. An Open dialog box opens.
3	Navigate to the SizeStandards folder, select the size standard file that you want to modify, and click Open . The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards
4	Edit the peak size values by following steps 7 to 9 on page 7-6.
5	Select Save (Ctrl+S) from the File menu to update an existing file, or Save As to specify a new name.

Using the Analysis Control Window to Edit an Existing Size Standard

Step	Action
1	If the Analysis Control window is not open, select Analysis Control (Ctrl+1) from the Windows menu.
2	Click the arrow in the Size Standard column for the sample that you want to change. 
3	Select a size standard from the pop-up menu.
4	Edit the peak size values by following steps 7 to 9 on page 7-6
5	Select Save (Ctrl+S), or Save As from the File menu to save changes and specify a new file name.

Deleting an Existing Size Standard


The following procedure describes how to delete a user-defined size standard from the Params Folder so that it no longer appears in the Size Standard pop-up menus. The size standard is permanently removed, and you must redefine it to use it again.

To delete an existing size standard:

Step	Action
1	Click the Start button, and then point to Programs .
2	Click Windows NT® Explorer to find and open the SizeStandards folder. The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards
3	Select the size standard file that you want to delete.
4	Drag the file to the Recycle Bin . Note You can also drag the size standard to another folder for storage.

Analyzing Samples
Using the Same
Size Standard

To select the same size standard to analyze all samples:


Step	Action								
1	If the Analysis Control window is not displayed, select Analysis Control (Ctrl+1) from the Windows menu.								
2	<p>Click the arrow in the Size Standard column heading and choose a size standard file from the pop-up menu.</p> <p>Your menu choice applies to all fields in the column.</p> <p>Note Alternatively, you can choose a value from the pop-up menu, click the header to select the entire column, and select Fill Down (Ctrl+D) from the Edit menu.</p> <div></div>								
3	<p>The pop-up menu contains the following options:</p> <table><thead><tr><th>Item</th><th>Description</th></tr></thead><tbody><tr><td>None</td><td>Apply no standard definition.</td></tr><tr><td>Collection Setting</td><td><p>Apply the size standard specified in the Data Collection software, which is embedded in the sample file.</p><p>For information on editing this file, see “Editing the Size Standard Definition” on page 7-9.</p></td></tr><tr><td>Custom standards that are listed at the bottom of the menu</td><td><p>These are files that you defined and they are located in the SizeStandards folder.</p><p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards</p></td></tr></tbody></table>	Item	Description	None	Apply no standard definition.	Collection Setting	<p>Apply the size standard specified in the Data Collection software, which is embedded in the sample file.</p> <p>For information on editing this file, see “Editing the Size Standard Definition” on page 7-9.</p>	Custom standards that are listed at the bottom of the menu	<p>These are files that you defined and they are located in the SizeStandards folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards</p>
Item	Description								
None	Apply no standard definition.								
Collection Setting	<p>Apply the size standard specified in the Data Collection software, which is embedded in the sample file.</p> <p>For information on editing this file, see “Editing the Size Standard Definition” on page 7-9.</p>								
Custom standards that are listed at the bottom of the menu	<p>These are files that you defined and they are located in the SizeStandards folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards</p>								

Selecting Separate Size Standards for Samples

To apply separate size standards to selected samples, click the arrow in the Size Standard column for the sample that you want to change (refer to the figure below) to open the pop-up menu.

For information on using the pop-up menu, see step 3 above.

Pop-up menu



Evaluating Analysis Results

8

Overview

In This Chapter Topics in this chapter include the following:

Topic	See Page
Process of Evaluating Analysis Results	8-2
Ways to Display Analysis Results	8-3
About the Results Display Window	8-4
Using the Results Control Window	8-7
Changing How the Results Are Displayed and Printed	8-13
About the Sample Results View	8-15
Updating the Results	8-16
Saving and Renaming the Results Control Format	8-17

Process of Evaluating Analysis Results

Evaluating Analysis Results

The following table describes the steps to evaluate the analysis results:

Step	Action	For information, see...
1	Displaying analysis results	"Ways to Display Analysis Results" on page 8-3
2	Using electropherogram and tabular data displays	"About Electropherogram and Tabular Data Displays" on page 9-2
3	Viewing electropherograms	"Displaying Electropherogram Data" on page 9-9
4	Verifying analysis results	"Process of Verifying Results" on page 9-31
5	Displaying other sample file data	◆ "Sample Info View" on page 4-11 ◆ "Size Curve View" on page 4-20 ◆ "Raw Data View" on page 4-22 ◆ "EPT Data View" on page 4-24
6	Using the analysis log	"Using the Analysis Log" on page 9-36
7	Remembering and renaming the results display	"Saving and Renaming the Results Control Format" on page 8-17

Ways to Display Analysis Results

Two Ways to Display Analysis Results


The following table describes two ways to display analysis results:

You can use...	Description	See Page
About the Results Display Window	This window is created from the Results Control window of a project. Use to group and view multiple sample files as electropherogram and tabular data.	8-4
About the Sample Results View	This view displays one sample file at a time, like the Results Display window, but is more convenient to view analysis results from a single sample file. It also allows quick access to supporting information views.	8-15

About the Results Display Window







Introduction The Results Display window is created from the Results Control window of a project. It allows you to group and view multiple sample files as electropherogram or tabular data. You can use the window to show up to eight panels, with multiple dye/samples per panel.

Displaying the Window When you select Results Control (Ctrl+2) from the Windows menu, a window opens like the example shown below. The callouts are described in the table that follows this example.






**Results Control
Window Callouts
Described**

The following table describes the callouts in the previous figure of the Results Control window:

Callout	Description												
1	Click one of the dye color fields (B, G, Y, R, O) to select dye/sample information.												
2	Click to display electropherograms for the selected samples.												
3	Click to show which sample is selected for this panel.												
4	Choose the number of electropherogram panels available for display from the pop-up menu.												
5	Click to display tabular data for the selected samples.												
6	Identifies the sample by row number and dye code. Sample information is displayed as specified in the project options.												
7	<table><tr><td colspan="2">Plot color indicator</td></tr><tr><td><table><tr><th>If you...</th><th>Then...</th></tr><tr><td>double-click the plot color indicator</td><td><p>the Choose a Plot Color dialog box opens.</p><div></div><p>For more information, refer to “Defining Individual Plot Colors” on page 9-26.</p></td></tr><tr><td>Ctrl+double-click the dye color indicator</td><td><p>the plot color indicator returns to the default color.</p><p>For more information, refer to “Setting Dye Indicator Preferences” on page 3-14.</p></td></tr></table></td><td></td></tr><tr><td colspan="2">Note If you change the plot color, a vertical line appears beside the indicator, showing that it has been modified.</td></tr></table>	Plot color indicator		<table><tr><th>If you...</th><th>Then...</th></tr><tr><td>double-click the plot color indicator</td><td><p>the Choose a Plot Color dialog box opens.</p><div></div><p>For more information, refer to “Defining Individual Plot Colors” on page 9-26.</p></td></tr><tr><td>Ctrl+double-click the dye color indicator</td><td><p>the plot color indicator returns to the default color.</p><p>For more information, refer to “Setting Dye Indicator Preferences” on page 3-14.</p></td></tr></table>	If you...	Then...	double-click the plot color indicator	<p>the Choose a Plot Color dialog box opens.</p> <div></div> <p>For more information, refer to “Defining Individual Plot Colors” on page 9-26.</p>	Ctrl+double-click the dye color indicator	<p>the plot color indicator returns to the default color.</p> <p>For more information, refer to “Setting Dye Indicator Preferences” on page 3-14.</p>		Note If you change the plot color, a vertical line appears beside the indicator, showing that it has been modified.	
Plot color indicator													
<table><tr><th>If you...</th><th>Then...</th></tr><tr><td>double-click the plot color indicator</td><td><p>the Choose a Plot Color dialog box opens.</p><div></div><p>For more information, refer to “Defining Individual Plot Colors” on page 9-26.</p></td></tr><tr><td>Ctrl+double-click the dye color indicator</td><td><p>the plot color indicator returns to the default color.</p><p>For more information, refer to “Setting Dye Indicator Preferences” on page 3-14.</p></td></tr></table>	If you...	Then...	double-click the plot color indicator	<p>the Choose a Plot Color dialog box opens.</p> <div></div> <p>For more information, refer to “Defining Individual Plot Colors” on page 9-26.</p>	Ctrl+double-click the dye color indicator	<p>the plot color indicator returns to the default color.</p> <p>For more information, refer to “Setting Dye Indicator Preferences” on page 3-14.</p>							
If you...	Then...												
double-click the plot color indicator	<p>the Choose a Plot Color dialog box opens.</p> <div></div> <p>For more information, refer to “Defining Individual Plot Colors” on page 9-26.</p>												
Ctrl+double-click the dye color indicator	<p>the plot color indicator returns to the default color.</p> <p>For more information, refer to “Setting Dye Indicator Preferences” on page 3-14.</p>												
Note If you change the plot color, a vertical line appears beside the indicator, showing that it has been modified.													



The following table describes the callouts in the previous figure of the Results Control window: *(continued)*

Callout	Description						
8	Dye color indicator						
	<table><tr><th>If you...</th><th>Then...</th></tr><tr><td>double-click the dye color indicator</td><td>the Choose a Dye Scale dialog box opens.  For more information, see “Changing the Dye Scale in Electropherograms” on page 9-28.</td></tr><tr><td>Ctrl+double-click the dye color indicator</td><td>the dye color indicator returns to the default scale. For more information, see “Changing the Dye Scale Preferences” on page 9-30.</td></tr></table>	If you...	Then...	double-click the dye color indicator	the Choose a Dye Scale dialog box opens.  For more information, see “Changing the Dye Scale in Electropherograms” on page 9-28.	Ctrl+double-click the dye color indicator	the dye color indicator returns to the default scale. For more information, see “Changing the Dye Scale Preferences” on page 9-30.
	If you...	Then...					
double-click the dye color indicator	the Choose a Dye Scale dialog box opens.  For more information, see “Changing the Dye Scale in Electropherograms” on page 9-28.						
Ctrl+double-click the dye color indicator	the dye color indicator returns to the default scale. For more information, see “Changing the Dye Scale Preferences” on page 9-30.						
Note If you change the scale, a vertical line appears beside the indicator, showing that it has been modified.							
9	Quick Tile buttons. See “Creating Tiled Electropherogram Displays” on page 8-10.						
10	Display button. See “Displaying the Results” on page 8-12.						
11	Clear Panel button. See “Removing Samples” on page 8-12.						
12	Print button. See “Printing the Results” on page 8-12.						
13	Clear All button. See “Removing Samples” on page 8-12.						

Using the Results Control Window


Selecting Display Format

Choose the results display format by clicking the electropherogram or tabular icons, or both. You must select one of the icons to display or print data.

Click this icon...	To display...	And the...
	electropherogram data	panel information below the icon is enabled. Refer to the figure below.
	tabular data	Tabular data appears in a single table, so electropherogram panel configuration is not relevant.

Selecting Electropherogram Panels

In electropherograms, the data appears in panels. You can overlay up to 16 samples within up to eight panels. Select the number of available panels (up to eight) from the pop-up menu labeled # of Panels.




Note Set the number of panels when using the Quick Tile option. Using this option causes the GeneScan® Analysis Software to automatically change panels as you select samples to display (see “Creating Tiled Electropherogram Displays” on page 8-10).

**Selecting Samples
to Display**

Use the dye color field on the left side of the Results Control window to select the dye to display.

You can select any or all dye colors, including the standard, for each sample file.

To select samples to display:

Step	Action
1	<div>Click the Electropherogram icon.</div> <div></div> <div>The panel below the button is enabled.</div>
2	<div>Click a panel number in the list to the left of the Dye/Samples list.</div> <div>If fewer than eight panels are available and you want to use more, choose a larger number (up to eight) from the # of Panels pop-up menu.</div>

To select samples to display: *(continued)*

Step	Action								
3	<p>Click the dye color fields to select or unselect the corresponding samples.</p> <p>Take the following action:</p> <table border="1"><thead><tr><th>To...</th><th>Click the...</th></tr></thead><tbody><tr><td>select the entire column</td><td>header of the column.</td></tr><tr><td>unselect the entire column</td><td>Clear All button, or click the header if it is already selected.</td></tr><tr><td>select all colors for a sample file</td><td>index number to the left of the color columns.</td></tr></tbody></table>	To...	Click the...	select the entire column	header of the column.	unselect the entire column	Clear All button, or click the header if it is already selected.	select all colors for a sample file	index number to the left of the color columns.
To...	Click the...								
select the entire column	header of the column.								
unselect the entire column	Clear All button, or click the header if it is already selected.								
select all colors for a sample file	index number to the left of the color columns.								

The samples corresponding to the selected dyes appear in the **Dye/Samples** list to the right of the **Sample File** list, as shown in the figure below.

Dye color of sample

Plot color of electropherogram

Row number and dye code followed by dye/sample information

To select samples to display: *(continued)*

Step	Action
4	When you change to a new panel, the dye colors of the samples you selected in other panels appear dark gray to indicate that they have been selected. You can select them again in the current panel.

Creating Tiled Electropherogram Displays

Procedure

To create tiled electropherogram displays:

Step	Action						
1	Choose the number of panels you want to display from # of Panels pop-up menu.						
2	Click the On button under Quick Tile .						
3	Select samples by clicking color fields. For information on: <table><tr><th>Topic</th><th>See Page</th></tr><tr><td>Selecting Samples to Display</td><td>8-8</td></tr><tr><td>Setting the tiled electropherogram preferences</td><td>8-13</td></tr></table> Each time a sample is selected, the program automatically changes to the next panel, so each selection is placed after the one containing your previous selection. After you select a sample for the last panel, the panel displays the first panel again.	Topic	See Page	Selecting Samples to Display	8-8	Setting the tiled electropherogram preferences	8-13
Topic	See Page						
Selecting Samples to Display	8-8						
Setting the tiled electropherogram preferences	8-13						

Example

The following table describes two examples of how to use the Quick Tile feature:


If you...	Then...
have four samples and choose four panels for display	click the column heading for the blue dye to select all four blue dye-labeled samples. The blue dye for each sample file appears in a separate panel.

The following table describes two examples of how to use the Quick Tile feature: *(continued)*

If you...	Then...
click the row index number for the first sample file to select all dye colors for one sample file	each dye for that sample file appears in a separate panel.

Unselecting Samples

To unselect the samples that you have selected for display:

Step	Action
1	<p>If you specified the Electropherogram display, click the panel number to the left of the Dye/Samples list to display the panel containing the samples you want to unselect.</p> <p>Dye color fields Panel numbers</p> 
2	<p>On the left side of the Results Control window, unselect the dye color fields corresponding to the samples you want to remove.</p> <p>The dye/sample identifiers are removed from the Dye/Samples list as you unselect samples.</p>

Removing Samples To remove samples:

Step	Action						
1	Select the panel number on the buttons to the left of the Dye/Samples list.						
2	You can take the following action: <table><tr><th>To remove...</th><th>Click...</th></tr><tr><td>all the samples you selected to display in a panel</td><td>Clear Panel.</td></tr><tr><td>the samples you have selected to display in all panels</td><td>Clear All.</td></tr></table>	To remove...	Click...	all the samples you selected to display in a panel	Clear Panel.	the samples you have selected to display in all panels	Clear All.
To remove...	Click...						
all the samples you selected to display in a panel	Clear Panel.						
the samples you have selected to display in all panels	Clear All.						

- Displaying the Results** To display the results on the screen, you can either:
- ◆ Click the Display button.
 - ◆ Press the Return or Enter key.


Printing the Results To print the results:

Step	Action
1	You can either: <ul style="list-style-type: none">◆ Click the Print button, or◆ Select Print (Ctrl+ P) from the File menu.
2	Click OK in the dialog box that opens.

Changing How the Results Are Displayed and Printed

Procedure You can set certain display preferences that remain in effect each time you display or print results data.

To change how results are displayed and printed:

Step	Action
1	<p>Select Preferences from the Settings menu and Results Display from the submenu.</p> <p>The following dialog box opens:</p> 
	<p>There are three Results Display preference categories:</p> <ul style="list-style-type: none">◆ Default Display Attributes◆ Stacked Electropherogram Panels◆ Printing Preferences

To change how results are displayed and printed: *(continued)*

Step	Action														
2	Set the Default Display Attributes to control the display attributes of new results displays, as follows:														
	<table><tr><th>You can select the...</th><th>For more information...</th></tr><tr><td>Align By Size check box</td><td>"Showing Data by Fragment Size" on page 9-20.</td></tr><tr><td>Show Peak Positions check box</td><td>"Displaying Peak Positions" on page 9-14.</td></tr><tr><td>Show Legends check box</td><td>"Using Legends to Change the Display" on page 9-15.</td></tr><tr><td>Show Offscale Region check box</td><td>"Showing Off-Scale Data" on page 9-18.</td></tr><tr><td>Standard or Custom Plot Colors buttons</td><td>"Defining Custom Colors" on page 9-25.</td></tr><tr><td>Opaque or Transparent Peak Highlighting buttons</td><td>"Highlighting Peaks" on page 9-15.</td></tr></table>	You can select the...	For more information...	Align By Size check box	"Showing Data by Fragment Size" on page 9-20.	Show Peak Positions check box	"Displaying Peak Positions" on page 9-14.	Show Legends check box	"Using Legends to Change the Display" on page 9-15.	Show Offscale Region check box	"Showing Off-Scale Data" on page 9-18.	Standard or Custom Plot Colors buttons	"Defining Custom Colors" on page 9-25.	Opaque or Transparent Peak Highlighting buttons	"Highlighting Peaks" on page 9-15.
	You can select the...	For more information...													
	Align By Size check box	"Showing Data by Fragment Size" on page 9-20.													
	Show Peak Positions check box	"Displaying Peak Positions" on page 9-14.													
	Show Legends check box	"Using Legends to Change the Display" on page 9-15.													
	Show Offscale Region check box	"Showing Off-Scale Data" on page 9-18.													
	Standard or Custom Plot Colors buttons	"Defining Custom Colors" on page 9-25.													
Opaque or Transparent Peak Highlighting buttons	"Highlighting Peaks" on page 9-15.														
3	Set the Stacked Electropherogram Panels , as follows:														
	<table><tr><th>Choose...</th><th>To set...</th></tr><tr><td>Use Common Vertical Scale check box</td><td>all panels in a display so they have the same vertical scale. The common scale is based on the electropherogram with the largest vertical scale.</td></tr><tr><td>Panel Height Resize Limits</td><td>minimum and maximum values for electropherogram panel height in the results display. Use to limit how much the electropherogram panels stretch or shrink to fit the size of the window.</td></tr></table>	Choose...	To set...	Use Common Vertical Scale check box	all panels in a display so they have the same vertical scale. The common scale is based on the electropherogram with the largest vertical scale.	Panel Height Resize Limits	minimum and maximum values for electropherogram panel height in the results display. Use to limit how much the electropherogram panels stretch or shrink to fit the size of the window.								
	Choose...	To set...													
Use Common Vertical Scale check box	all panels in a display so they have the same vertical scale. The common scale is based on the electropherogram with the largest vertical scale.														
Panel Height Resize Limits	minimum and maximum values for electropherogram panel height in the results display. Use to limit how much the electropherogram panels stretch or shrink to fit the size of the window.														

To change how results are displayed and printed: *(continued)*

Step	Action								
4	Set the Printing Preferences , as follows:								
	<table><tr><th>If...</th><th>Then...</th></tr><tr><td>the height of the panels that appears on screen is acceptable</td><td>select the As Shown on Screen button.</td></tr><tr><td>you want to print the electropherogram at a specified height</td><td>click the Fixed at button and enter a value in the field.</td></tr><tr><td>you want to force a page break after the electropherograms have been printed</td><td>select the Page Break before Tabular Data check box.</td></tr></table>	If...	Then...	the height of the panels that appears on screen is acceptable	select the As Shown on Screen button.	you want to print the electropherogram at a specified height	click the Fixed at button and enter a value in the field.	you want to force a page break after the electropherograms have been printed	select the Page Break before Tabular Data check box.
	If...	Then...							
	the height of the panels that appears on screen is acceptable	select the As Shown on Screen button.							
you want to print the electropherogram at a specified height	click the Fixed at button and enter a value in the field.								
you want to force a page break after the electropherograms have been printed	select the Page Break before Tabular Data check box.								
5	Click OK .								

About the Sample Results View

About the View The Sample Results view is displayed within a Sample File window. You can access the window through a sample file or through a project's Analysis Control or Results Control window.

If...	Then...
you are opening a sample file as a stand-alone file	the Sample Results view is the default display within the Sample File window.
For more information on...	See Page
the Sample Results view	4-9

Updating the Results

Re-analyzing the Data	<div><div></div><div>The Results Control and the Sample Results windows are dynamic.</div><div>If you re-analyze your data with either window active, then the software updates this window.</div></div>
------------------------------	--


Saving and Renaming the Results Control Format

Introduction You can use the Results Control window to view multiple sample files in electropherogram and tabular format. The GeneScan Analysis Software allows you to save formats for future use. You can then redisplay or print these formats without having to redefine them again.

Important Considerations The following are important considerations for saving a Results Control format:

- ◆ You must save the project for the display to be available when you open the project again.
- ◆ Remembering a display preserves the combination of windows/panels/data and customized color settings. It does not preserve any zooming you have performed.


Saving the Display Format To save a display format for future viewing:

Step	Action
1	With the Results Control window set for the display, either: <ul style="list-style-type: none">◆ Click the Display button, or◆ Press the Enter key.
2	With the display on the screen, select Remember Display from the Project menu. The Remember Display dialog box opens. 
3	Enter a name for the display and click OK .

**Working with a
Previously Saved
Display**


Use the Previous Displays dialog box to display, print, remove or rename a saved display.

To work with a previously saved displays:

Step	Action										
1	<p>Select Previous Displays from the Project menu.</p> <p>The following is an example of the Previous Displays dialog box:</p> 										
2	<p>Select a display or multiple displays and take the following action:</p> <table><tr><th>To...</th><th>Click an item in the list and click...</th></tr><tr><td>display the saved formats</td><td>Display.</td></tr><tr><td>print the saved formats</td><td>Print. The standard print dialog box opens.</td></tr><tr><td>remove the saved formats</td><td>Remove. An alert opens.</td></tr><tr><td>rename a currently saved format</td><td>Rename. The Rename dialog box opens.</td></tr></table> <p>Note You can only rename one display at a time.</p> <p>See “Renaming the Current Results Display” on page 8-19.</p>	To...	Click an item in the list and click...	display the saved formats	Display.	print the saved formats	Print. The standard print dialog box opens.	remove the saved formats	Remove. An alert opens.	rename a currently saved format	Rename. The Rename dialog box opens.
To...	Click an item in the list and click...										
display the saved formats	Display.										
print the saved formats	Print. The standard print dialog box opens.										
remove the saved formats	Remove. An alert opens.										
rename a currently saved format	Rename. The Rename dialog box opens.										

Renaming the Current Results Display

To rename the Results display that is currently on the screen and to save the display under a different name:

Step	Action
1	Ensure that the display is the active window.
2	<p>Select Rename Display from the Project menu.</p> <p>The following is an example of the Rename Display dialog box:</p> 
3	<p>Enter a new name for the display and click OK.</p> <p>The new name opens in the Previous Displays dialog box.</p>

Evaluating Electropherograms

9

Overview

In This Chapter Topics in this chapter include the following:

Topic	See Page
About Electropherogram and Tabular Data Displays	9-2
Displaying Electropherogram and Tabular Data	9-4
Displaying Electropherogram Data	9-9
Working with Electropherogram Data	9-12
Defining Custom Colors in Electropherograms	9-24
Changing the Dye Scale in Electropherograms	9-28
Process of Verifying Results	9-31
Verifying Size Calculations	9-33
Using the Analysis Log	9-36
Verifying Peak Detection	9-38

About Electropherogram and Tabular Data Displays

Introduction After analyzing the data, you can display the results for each sample in electropherogram and tabular data. You can also customize the electropherogram and tabular data display.

Note Altering the appearance of the electropherograms and the tabular data displays does not change the analyzed data contained in the sample file on which they are based.

How the Window Is Divided When electropherogram and tabular data are displayed together, the window is divided into upper and lower windows.

Window	Contains
Upper window	electropherogram data
Lower window	tabular data
For more information on...	See
customizing the window's appearance by adjusting the relative size of each window	"Adjusting Window Size" on page 9-7.

What Tabular Data Contains If you analyze samples and perform sizecalling, the tabular data contains the estimated sizes (in base pairs) of all detected fragments. Use this information for detailed data analysis and further calculations. The peaks matched to the defined size standard are identified by dots next to the Dye/Sample Peak field.

Sample peaks that are larger (in base pairs) than the largest defined peak in the selected standard are not sized. The corresponding size fields are blank.

Note Tabular data displays only peaks that are detected based on the Dye Amplitude Thresholds and Minimum Peak Half Width setting of the analysis parameters.

How Electropherogram Panels Are Sized In the Results Display window, the GeneScan® Analysis Software sizes all electropherogram panels to fit within the electropherogram portion of the window by using the largest size that fits them into the visible area. You can scroll to see the portion of the display that is not visible.




**For More
Information**

For more information, see the following topics:

Topic	See Page
Displaying Electropherogram and Tabular Data	9-4
Displaying Electropherogram Data	9-9
Working with Electropherogram Data	9-12

Displaying Electropherogram and Tabular Data

Procedure To display electropherogram and tabular data:

Step	Action								
1	Select Results Control (Ctrl+2) from the Windows menu to open the Results Control window.								
2	Click the Electropherogram button  and the Tabular button  .								
3	If applicable, select the number of electropherogram panels from the # of Panels list.								
4	<div>Click the dye color fields to select or unselect the corresponding samples.</div> <div>Dye color fields</div> <div></div> <div>Take the following action:</div> <table><tr><th>To...</th><th>Click the...</th></tr><tr><td>select the entire column</td><td>header.</td></tr><tr><td>unselect the entire column</td><td>Clear All button.</td></tr><tr><td>select all colors for a sample file</td><td>index number to the left of the color columns.</td></tr></table> <div>The samples corresponding to the selected dyes appear in the Dye/Samples list to the right of the Sample File list.</div>	To...	Click the...	select the entire column	header.	unselect the entire column	Clear All button.	select all colors for a sample file	index number to the left of the color columns.
To...	Click the...								
select the entire column	header.								
unselect the entire column	Clear All button.								
select all colors for a sample file	index number to the left of the color columns.								
5	<div>Click Display.</div> <div>The electropherogram and tabular data are displayed in the Results Display window.</div> <div>See also “Working with Electropherogram Data” on page 9-12.</div>								

Example of Tabular Data and Electropherogram

The following is an example of tabular data with a corresponding electropherogram:




Table Describing Columns

The following table describes the columns in the “Example of Tabular Data and Electropherogram” above:

Column heading	Identifies...
Dye/Sample Peak	<ul style="list-style-type: none"> ◆ Sample index number ◆ Dye color ◆ Peak number
Minutes	The time, in minutes, from the start of the run to the time the fragment was detected
Size	<p>The differences in fragment mobility</p> <p>This value is calculated automatically only if you:</p> <ul style="list-style-type: none"> ◆ Run the size standard in the same lane or injection as the sample, and ◆ Perform sizecalling
Peak Height	Signal size (RFU)

The following table describes the columns in the “Example of Tabular Data and Electropherogram” above: *(continued)*

Column heading	Identifies...
Peak Area	Area of the detected peak
Data point	Data point of the fragment at its maximum peak height

**Why Some Peaks
May Be Visible
Only in an
Electropherogram**

Peaks may be visible in the electropherogram and not listed in the tabular data because:

Reason	For more information, see...
The software detects the peaks based on the Peak Amplitude Thresholds and Min Peak Half Width.	“Peak Detection Options” on page 5-8
Electropherograms display the peaks that fall within the range specified by the Sizecall Range parameters that are defined in the Analysis Parameters dialog box.	“Sizecall Range Options” on page 5-10

**Why Some Peak
Areas May Have a
Negative Number**

If negative values appear in peak areas in the electropherogram, it is because a portion of the peak is below the baseline. The GeneScan Analysis Software display does not show the part of the electropherogram that is below the baseline.

**Highlighting
Information**

To highlight information for one peak in the electropherogram and tabular data:

Click...	Then...
the peak in the electropherogram	the peak fills with color and the corresponding row in the tabular data window is highlighted.
the Dye/Sample Peak number in the tabular data window	highlights the corresponding peak in the electropherogram.

Changing the Highlight Color


To change the highlight transparency:

Step	Action
1	Select Peak Highlighting from the View menu.
2	Then select either Opaque or Transparent from the submenu.

Note For more information on highlighting peaks, see page 9-15.

Adjusting Window Size

To adjust the relative size of the electropherogram and tabular windows:

Step	Action
1	Move the cursor to the window divider (the double line between the two windows).
2	When the cursor changes to a bidirectional arrow (), click the window divider line and drag it up or down.

Hiding Selected Rows of Data

To hide selected rows of data:

Step	Action	
1	Take the following action:	
	If you want to...	Then...
	select a row	either: <ul style="list-style-type: none">◆ Click the first field in the row, or◆ Click the corresponding peak in the electropherogram.
	select several rows that are not next to each other	Ctrl+click the rows.
2	Select Hide Selected Rows (Ctrl+H) from the View menu.	

**Limiting the Rows
to Display**

To limit the display to the selected rows of data:


Step	Action
1	Select the rows you want to display.
2	Select Show ONLY Selected Rows (Ctrl+G) from the View menu. Note Select Show All Rows (Ctrl+G) from the View menu to display all of the tabular data after limiting the display.

Displaying Electropherogram Data

Definition Each electropherogram provides a profile of the selected dye samples it represents. The y-axis represents the relative fluorescence of the detected fragments as they occurred over time. The x-axis represents time and can be displayed by data points or base pairs.


Base Pairs Versus Data Points The tick marks on the x-axis can represent size in base pairs instead of data points. This option is only available for runs that include an internal size standard with the sample (see “Showing Data by Fragment Size” on page 9-20).

Procedure to Display Data To display electropherogram data:

Step	Action								
1	Select Results Control (Ctrl+2) from the Windows menu to open the Results Control window.								
2	Click the Electropherogram button.								
3	If applicable, select the number of electropherogram panels from the # of Panels pop-up menu.								
4	<div>Click the dye color fields to select or unselect the corresponding samples.</div> <div>Dye color fields</div> <div></div> <div>Take the following action:</div> <table><tr><th>To...</th><th>Click the...</th></tr><tr><td>select the entire column</td><td>header.</td></tr><tr><td>unselect the entire column</td><td>Clear All button.</td></tr><tr><td>select all colors for a sample file</td><td>index number to the left of the color columns.</td></tr></table> <div>The samples corresponding to the selected dyes appear in the Dye/Samples list to the right of the Sample File list.</div>	To...	Click the...	select the entire column	header.	unselect the entire column	Clear All button.	select all colors for a sample file	index number to the left of the color columns.
To...	Click the...								
select the entire column	header.								
unselect the entire column	Clear All button.								
select all colors for a sample file	index number to the left of the color columns.								
5	<div>Click Display.</div> <div>The electropherogram is displayed in the Results Display window as shown in the example below.</div>								

Electropherogram Example

The following is an example of an electropherogram:



Electropherogram Callouts Described

The following table describes the callouts in the figure above:

Call out	Description	See...
1	Cross hairs	“Displaying X- and Y-Axis Positions” on page 9-13
2	Magnifying glass Use to zoom in a specific area or hold down while pressing the Alt key to zoom out to a smaller scale. Or, click and drag a marque around an area to zoom in to that area.	“Zooming In and Out” on page 9-17
3	Legend Text from the sample file that appear beneath electropherogram panels in the Results Display window.	“Using Legends to Change the Display” on page 9-15

The following table describes the callouts in the figure above: *(continued)*

Call out	Description	See...
4	Dye color indicator	“Dye color indicator” on page 8-6
5	Plot color indicator	“Plot color indicator” on page 8-5
6	Scroll bar Use the scroll bar to scroll horizontally.	“Scrolling the Display” on page 9-16

Working with Electropherogram Data

In This Section This section describes how to perform the following tasks:

Task	See Page
Displaying X- and Y-Axis Positions	9-13
Moving the Electropherogram	9-13
Changing the Dye Color	9-14
Displaying Peak Positions	9-14
Highlighting Peaks	9-15
Using Legends to Change the Display	9-15
Scrolling the Display	9-16
Zooming In and Out	9-17
Showing Off-Scale Data	9-18
Electropherogram Displaying Off-Scale Data	9-19
Electropherogram Displaying the Flat-Topped Effect	9-19
Showing Data by Fragment Size	9-20
Changing the Horizontal Scale	9-21
Changing the Vertical Scale	9-22
Assigning Standard or Custom Colors	9-23

Displaying X- and Y-Axis Positions


The following table describes how to display the x- and y-axis positions:

To...	Then...
display the x- and y-axis positions.	click the cross hairs and select an area in the Electropherogram. The x- and y-axis values appear in the box in the lower left corner of the electropherogram. If tabular data is also displayed, the row in the table is highlighted.

Moving the Electropherogram

Move the associated electropherogram to the front by clicking one of the following in the legend:

- ◆ Dye scale indicator
- ◆ Plot color indicator, or
- ◆ Text




Changing the Dye Color

The following table describes how to change the dye color and how to return to the default dye color:

To...	Then...
change the dye color	double-click the dye color indicator.
return to the default dye color	Ctrl+double-click the dye color indicator.
If...	Then...
you change the dye color or scale	a vertical line appears beside the indicator, showing that it has been modified.

Displaying Peak Positions

Use Show Peak Positions from the View menu to examine how the GeneScan Analysis Software defines peaks by displaying markers that identify the beginning, center, and end of each peak.



Highlighting Peaks Use the Peak Highlighting command to highlight a selected peak with the dye/sample's plot color.

To highlight selected peaks:



Step	Action						
1	Select Peak Highlighting from the Views menus and either Opaque or Transparent from the submenu.						
	<table><tr><th>Use this option...</th><th>To...</th></tr><tr><td>Opaque</td><td>fill the peak with a solid color that can obscure peaks behind the selected peak.</td></tr><tr><td>Transparent</td><td>use a slightly diffused plot color that allows you to view overlapping peaks.</td></tr></table>	Use this option...	To...	Opaque	fill the peak with a solid color that can obscure peaks behind the selected peak.	Transparent	use a slightly diffused plot color that allows you to view overlapping peaks.
	Use this option...	To...					
	Opaque	fill the peak with a solid color that can obscure peaks behind the selected peak.					
Transparent	use a slightly diffused plot color that allows you to view overlapping peaks.						
2	Click a detected peak in an electropherogram. The peak is highlighted.						

Using Legends to Change the Display The following table shows how to use legends to change how electropherograms are displayed:

If you want to...	Then...
show or hide legends	select Show Legends from the View menu.
open sample file windows	double-click the corresponding legend text.
reorganize overlaid electropherograms	a. Display the electropherograms with legends. b. Click either the dye scale indicator, plot color indicator, or the text for the sample you want to move to the front.

Scrolling the Display

The following table describes ways to scroll the display:

Use the	Description	
scroll bar	If you want to...	Then...
	shift the electropherogram to the right or the left	click in the gray region of the scroll bar to the right or left of the scroll box.
	scroll across the electropherogram	click an arrow at the end of the scroll bar.
	control the amount of scroll	drag the scroll bar to the right or the left.
scroller symbols	<p>a. Hold the mouse cursor over either the vertical or horizontal scale of the electropherogram.</p> <p>Either a vertical scroller symbol () or a horizontal scroller symbol () appears.</p> <p>b. Hold down the mouse button and move the mouse in the direction of the information you want to view.</p>	

Zooming In and Out

Out

About Zooming In and Out

By default, the GeneScan Analysis Software scales each electropherogram horizontally to show all peaks detected during the run. While this provides a good overview of the run, some peaks may be quite compressed.

Improving Visibility

To improve visibility, you can change the horizontal scale of the electropherograms by zooming.

Zooming affects:

- ◆ Only the horizontal scale, and zooms the middle portion of the window
- ◆ All displayed electropherogram panels

How to Change the View Scale

If you want to...	Then...
see views with greater detail	<ul style="list-style-type: none">◆ Select Zoom In (Ctrl++) from the View menu, or◆ Click the magnifying glass in the upper-left corner of the window, and drag around a specific area to zoom in on the electropherogram.
see a smaller scale view of the data after zooming in	<ul style="list-style-type: none">◆ Select Zoom Out (Ctrl+ -) from the View menu, or◆ Click the magnifying glass cursor, hold down the Alt key, and click the electropherogram. <p>The data appears in successively smaller scale views.</p>
quickly scale the data so that the entire length fits within the window, again	Select Zoom Out (Full Range) from the View menu.


Showing Off-Scale
Data

This section contains the following information:

Topic	See Page
Procedure	9-18
About Flat-Topped Peaks	9-18
Electropherogram Displaying Off-Scale Data	9-19
Electropherogram Displaying the Flat-Topped Effect	9-19

Procedure

To show off-scale data:


Step	Action				
1	Select Preferences from the Settings menu and Results Display from the submenu. The Results Display Preferences dialog box opens.				
2	Select the Show Offscale Regions check box to highlight with a red bar regions in the electropherogram that contain off-scale data (see “Electropherogram Displaying Off-Scale Data” on page 9-19).  Note Select the Zoom In (Ctrl++) command from the View menu to more clearly show the areas of off-scale data. <table><tr><th>If...</th><th>Then...</th></tr><tr><td>the sample was sized</td><td>the Analysis Log lists the numbers of off-scale regions in the analysis range for each sample file.</td></tr></table> Note You can toggle this command for individual electropherograms by selecting Hide/Show Offscale Regions (Ctrl+`) from the View menu.	If...	Then...	the sample was sized	the Analysis Log lists the numbers of off-scale regions in the analysis range for each sample file.
If...	Then...				
the sample was sized	the Analysis Log lists the numbers of off-scale regions in the analysis range for each sample file.				

About Flat-Topped Peaks

An additional feature is that peaks that contain off-scale data points are drawn in the electropherograms as “flat topped;” that is, the top section of the peak is flat rather than pointed (see “Electropherogram Displaying the Flat-Topped Effect” on page 9-19).

Select Analysis Parameters from the Settings Menu (see “Data Processing Options” on page 5-7).

The following is an example of an electropherogram displaying off-scale data:



Showing Data by Fragment Size

The following table explains how to show data by fragment size:

If...	Then...
you analyze your samples with an internal size standard	<p>use the Align by Size (Ctrl+T) command from the View menu to align the horizontal scale of the electropherograms by fragment size instead of by data point.</p> <p>Note You can display data by size only if you analyzed and performed sizecalling of your samples using a size standard.</p>
For example, if	Then...
you run two identical samples in different runs	<p>the Align by Size command adjusts for run-to-run variations by aligning peaks by size value. This eliminates any apparent differences that were caused by run discrepancies.</p> <p>Note You can display overlaid samples in the same dye in different colors. See “Defining Custom Colors in Electropherograms” on page 9-24.</p>


How to Switch Between Size and Data Point Display

If you want to...	Then...
show data by size	<p>select Align by Size (Ctrl+ T) from the View menu.</p> <p>When the data is aligned by size, the menu command changes to Align by Data Point.</p> <p>Select the command again, to show the data aligned by data point.</p>
set the default peak alignment	<p>select Preferences from the Settings menu, and select Results Display from the submenu.</p> <p>You can use the Results Display Preferences dialog box to set certain preferences that remain in effect each time you display or print results data.</p> <p>For more information, see “Changing How the Results Are Displayed and Printed” on page 8-13.</p>

Changing the Horizontal Scale

Changing the Horizontal Scale for All Electropherograms

To change the scale of the horizontal axis for all electropherograms:

1	Display the electropherogram panels you want to change.
2	<p>Select Horizontal Scale from the View menu.</p> <p>The Horizontal Scale Parameters dialog box opens.</p> <p>You can also move the cursor over the horizontal axis of a displayed electropherogram and double-click.</p> 
3	Enter the increments represented by the tick marks for the horizontal axis in the Tick Spacing box.
4	Enter a range in the entry fields labeled Display from and Display to .
5	Click OK .

Changing the Horizontal Scale for Individual Electropherograms


To change the horizontal scale for individual electropherograms:

Step	Action
1	Display the electropherograms.
2	<p>Move the cursor over the horizontal axis of the panel that you want to change and double-click.</p> <p>The Horizontal Scale Parameters dialog box opens.</p>
3	Enter the increments represented by the tick marks for the horizontal axis in the Tick Spacing box.
4	Enter a range in the entry fields labeled Display from and Display to .
5	<p>Click OK.</p> <p>The horizontal scale changes.</p>

**Changing the
Vertical Scale**


Changing the Vertical Scale for All Electropherograms

To change the vertical scale for all electropherograms:

Step	Action
1	Display the electropherogram panel that you want to change.
2	Select Vertical Scale from the View menu. The Vertical Scale dialog box opens. 
3	Enter the increments represented by the tick marks for the vertical axis in the Tick Spacing box.
4	Enter a range in the entry fields labeled Display from and Display to .
5	Click OK .

Changing the Vertical Scale for Individual Electropherograms

To change the vertical scale for individual electropherograms:

Step	Action
1	Display the electropherograms.
2	<p>Move the cursor over the vertical axis of the panel that you want to change, and double-click.</p> <p>The following dialog box opens:</p> 
3	Enter tick mark increments and a range.
4	<p>Ensure that the Apply to all Electropherogram panels check box is not selected.</p> <p>Select the check box only to apply the changes to all displayed electropherogram panels.</p>
5	<p>Click OK.</p> <p>The vertical scale changes only for the electropherogram panel that you selected.</p>

Assigning Standard or Custom Colors

The following table describes how to assign standard or custom colors and where to look for more information:

To...	Then...
assign standard or custom colors	<p>select the Plot Colors command from the View menu and either Standard or Custom from the submenu.</p> <p>For information on defining custom colors, see "Defining Custom Colors in Electropherograms" on page 9-24.</p>

Defining Custom Colors in Electropherograms

- Introduction

The GeneScan Analysis Software assigns a plot color to each dye/sample added to an electropherogram. Normally, it is the color associated with the individual dye/sample by the Dye Indicators Preferences.
- Note

To change the default dye colors in the Analysis Control window and the Results displays, see “Setting Dye Indicator Preferences” on page 3-14.
- Note

Custom plot colors are not available in the Sample Results view.

- Why Change Colors in the Electropherogram

Change the colors in the electropherogram to:

- ◆ Differentiate between different samples labeled with the same color dye
 - ◆ Improve contrast between different dye colors
 - ◆ Show data in a special color for a presentation
 - ◆ Optimize plot colors for a particular printer



Saving the Display Format

The following table describes the options to save the display format after customizing the display colors:

If you...	And then...
specify saving the display format of a Results Display window after customizing the display	open it at a later time, the custom colors still appear. For more information, see “Saving and Renaming the Results Control Format” on page 8-17.
do not save the display format after manually customizing the colors	display the same results again. The electropherograms are redrawn using default colors.

Defining Custom Colors

To define custom colors for all electropherograms:



Step	Action
1	<p>Select Project Options from the Settings menus and Choose Custom Plot Colors from the submenu.</p> <p>The Custom Plot Colors dialog box opens.</p> 
2	<p>Select new colors from the pop-up menus beside the 16 plot numbers.</p> <p>Note The plot numbers indicate the order of the samples in the electropherogram legend.</p>
3	<p>Select Other from the pop-up menu to specify a color that does not appear in the pop-up menu.</p> <p>The following is an example of a color picker that opens:</p> 

To define custom colors for all electropherograms: *(continued)*

Step	Action
4	Position the cross hairs pointer on the color you want, and click. The color appears in the Custom Color box.
5	To save the custom color, click Add to Custom Colors .
6	Click OK . The Color Picker window closes.
7	In the Custom Plot Colors dialog box, select the check box labeled Save As Defaults to save the customized colors.
8	Click OK .

Defining Individual Plot Colors

To define the individual plot colors:

Step	Action
1	Display the electropherograms with legends.
2	<p>Double-click the plot color indicator next to the sample you want to change.</p>  <p>Plot color indicator</p> <p>The Choose a Plot Color dialog box opens:</p>  <p>Pop-up menu</p>
3	<p>Select a new color from the pop-up menu.</p> <p>If you select Other, then a color picker opens.</p>

To define the individual plot colors: *(continued)*

Step	Action
4	<p>Click OK.</p> <p>The color of the electropherogram for the individual sample changes, and a vertical line appears beside the plot color indicator to signify that it has been modified.</p> <p>Note You can change the plot color in the same way from the Results Control window. When you do so, the dye/sample is plotted with the set color each time you open the applicable Results Display window.</p> <p>Note Press Ctrl and double-click the plot color indicator to reset it to the original color.</p>

Changing the Dye Scale in Electropherograms

What the Dye Scale Defines The dye scale defines how dyes in an electropherogram appear relative to each other. You can compensate for peaks with different intensities by redefining the dye scale.



Note Changing the dye scale affects only the display, not the underlying data.

Increasing the Dye Scale Example The following table describes one way to increase the dye scale:

If...	Then...	Action...
you loaded a smaller amount of green sample in relation to the red sample	the peaks for the green sample might appear half as tall as those of the red sample.	To make it easier to view both samples on the same scale, increase the dye scale value of the green sample to make the peaks appear similar.


Changing the Dye Scale of an Electropherogram

To change the dye scale of an individual electropherogram:

Step	Action
1	Display the electropherograms with legends.
2	<p>Double-click the dye color indicator next to the sample you want to change.</p>  <p>Plot color indicator</p> <p>The Choose a Dye Scale dialog box opens.</p> 
3	<p>Enter a new scale in the dialog box and click OK.</p> <p>The dye scale for the individual sample changes, and a vertical line appears beside the dye color indicator to signify that it is modified.</p> <p>Note You can change the dye color in the same way from the Results Control window. The dye is scaled each time you open the Results Display window.</p>

Changing the Dye
Scale Preferences

To change the dye scale preferences:

Step	Action
1	<p>Select Preferences from the Settings menu and select Results Dye Scales from the submenu.</p> <p>The Preferences dialog box opens with the Results Dye Scales pop-up menu displayed.</p> 
2	<p>Enter a positive number between 0.1 and 100 for each sample relative to any other sample, and click OK.</p> <p>Note Dye scale values do not automatically revert to default values. Change them back to the defaults before examining results of another run.</p>

Process of Verifying Results

Introduction You can use the electropherogram and tabular displays to verify the results of analysis by checking the GeneScan Analysis Software calculated sizes and peaks.

Note The sizecalling of the standard and of sample fragments varies according to the sizecalling method you defined in the Analysis Parameters and the accuracy of the defined standard.

Steps to Verify Size Calculation

To verify size calculations:

Step	Action	See
1	Compare how well multiple size standard electropherograms line up within a Results Display window when aligned by size.	"Verifying Size Calculations" on page 9-33.
2	View the sizing curve calculated by the GeneScan Analysis Software.	"Example of Size Curve View" on page 4-21.
3	Determine how well the defined size standard matches the size standard run with your sample. Use the Peak Total information in the Sample Info view of the Sample File window.	"Sample Info View" on page 4-11.
4	View the Analysis Log , which provides messages for each analyzed sample file. If there is a problem or a questionable condition during sizecalling, a warning message is displayed in the Analysis Log .	"Using the Analysis Log" on page 9-36.
5	Use the Raw Data view to display information about the raw data for a sample. Analyzed sample files contain raw and analyzed data.	"Example of Raw Data View" on page 4-23.

To verify size calculations: *(continued)*

Step	Action	See
6	<p>Use EPT Data to troubleshoot problems caused by poor run conditions, such as:</p> <ul style="list-style-type: none">◆ EP voltage◆ EP current◆ Laser power◆ Run temperature versus time <p>EPT data can be displayed for each sample file.</p>	<p>“EPT Data View Example” on page 4-25.</p>

Verifying Size Calculations

Introduction This section describes the following topics:

- ◆ Verifying for the GeneScan-350 Standard
- ◆ Evaluating for Multiple Size Standards

Verifying for the GeneScan-350 Standard


About the Standard

In the GeneScan-350 standard, the size of the first peak should be approximately 50 bp, the second 75 bp, and so on, assuming that processing started after the 35-bp fragment passed the scan region. If peaks appear to be correctly measured for your run, measurement of the sample fragments that ran with the standard should also be correct.


Note For a complete list of fragment sizes, refer to Appendix C, “GeneScan Size Standards.”

Procedure

To verify the size calculation for the GeneScan-350 standard:


Step	Action
1	In the Results Control window, select the Electropherogram button and the Tabular button.
2	Click the Clear All button to clear the panels.
3	Click the dye color for the size standard you want to view. Dye color fields 
4	Click Display . Note You can also open the Sample File window for the sample file of interest to verify sizecalling.

To verify the size calculation for the GeneScan-350 standard: *(continued)*

Step	Action																																																																		
5	<p>Click each peak in the electropherogram and check the tabular data to ensure it is the correct size (see figure below).</p> <p>Note The peaks matched to the defined size standard are identified by dots next to the Dye/Sample Peak field.</p> <p>Indicates peak matched to defined size standard</p>  <table><tr><th>Dye/Sample Peak</th><th>Minutes</th><th>Size</th><th>Peak Height</th><th>Peak Area</th><th>Data Point</th></tr><tr><td>1R_1</td><td>22.06</td><td>50.00</td><td>288</td><td>4298</td><td>2685</td></tr><tr><td>1R_2</td><td>23.01</td><td>56.75</td><td>18</td><td>222</td><td>2801</td></tr><tr><td>1R_3</td><td>23.15</td><td>57.68</td><td>16</td><td>177</td><td>2817</td></tr><tr><td>1R_4</td><td>23.33</td><td>58.96</td><td>54</td><td>609</td><td>2839</td></tr><tr><td>1R_5</td><td>23.48</td><td>60.00</td><td>263</td><td>3103</td><td>2857</td></tr><tr><td>1R_6</td><td>23.61</td><td>60.91</td><td>33</td><td>430</td><td>2873</td></tr><tr><td>1R_7</td><td>23.84</td><td>62.55</td><td>14</td><td>85</td><td>2902</td></tr><tr><td>1R_8</td><td>23.89</td><td>62.89</td><td>14</td><td>80</td><td>2908</td></tr><tr><td>1R_9</td><td>23.98</td><td>63.51</td><td>10</td><td>69</td><td>2919</td></tr><tr><td>1R_10</td><td>24.17</td><td>64.82</td><td>12</td><td>30</td><td>2942</td></tr></table>	Dye/Sample Peak	Minutes	Size	Peak Height	Peak Area	Data Point	1R_1	22.06	50.00	288	4298	2685	1R_2	23.01	56.75	18	222	2801	1R_3	23.15	57.68	16	177	2817	1R_4	23.33	58.96	54	609	2839	1R_5	23.48	60.00	263	3103	2857	1R_6	23.61	60.91	33	430	2873	1R_7	23.84	62.55	14	85	2902	1R_8	23.89	62.89	14	80	2908	1R_9	23.98	63.51	10	69	2919	1R_10	24.17	64.82	12	30	2942
Dye/Sample Peak	Minutes	Size	Peak Height	Peak Area	Data Point																																																														
1R_1	22.06	50.00	288	4298	2685																																																														
1R_2	23.01	56.75	18	222	2801																																																														
1R_3	23.15	57.68	16	177	2817																																																														
1R_4	23.33	58.96	54	609	2839																																																														
1R_5	23.48	60.00	263	3103	2857																																																														
1R_6	23.61	60.91	33	430	2873																																																														
1R_7	23.84	62.55	14	85	2902																																																														
1R_8	23.89	62.89	14	80	2908																																																														
1R_9	23.98	63.51	10	69	2919																																																														
1R_10	24.17	64.82	12	30	2942																																																														
6	<p>If you find that the sizes were not calculated correctly, you can:</p> <ul style="list-style-type: none">◆ Redefine the size standard, or◆ Change the analysis parameters and re-analyze the affected samples (see “Sizecaller Algorithm Flowchart” on page 5-4)																																																																		

Evaluating for Multiple Size Standards

To evaluate size calculations for multiple size standards:

Step	Action
1	In the Results Control window, select the Electropherogram button.
2	Click the Clear All button to clear the panels.
3	Click the On button to turn on the Quick Tile option.
4	<p>Click the dye colors for the size standards you want to view.</p> <p>Dye color fields</p>  <p>When the Quick Tile option is on, the GeneScan Analysis Software inserts each in a separate panel.</p> <p>Note Click the header of the appropriate dye/sample column to display all standards in the project that are the same color dye.</p>
5	<p>Click Display.</p> <p>The standards appear in tiled electropherogram displays.</p> <p>For information on setting the tiled electropherogram displays, see “Creating Tiled Electropherogram Displays” on page 8-10.</p>
6	<p>Select Align by Size (Ctrl+T) from the View menu if the electropherograms are not already aligned by size.</p> <p>The size standards should line up when aligned by size.</p> <p>Note You can set preferences so that all new displays show data aligned by size by selecting Preferences from the Settings menu and Results Display from the submenu. Click the Align By Size check box.</p> <p>For more information, see “Saving and Renaming the Results Control Format” on page 8-17.</p>

Using the Analysis Log


What Is the Analysis Log

The Analysis Log maintains a running record of analysis performed by the GeneScan Analysis Software. If a problem occurs during analysis of a sample file, the Analysis Log automatically opens in the foreground as an alert.

Displaying the Analysis Log

Select Analysis Log (Ctrl+O) from the Windows menu.

The following is an example of the Analysis Log:



What to Evaluate Evaluate the following:

What to evaluate	What not to evaluate
Potential problems that the GeneScan Analysis Software might have had during sizecalling. Analysis Log will alert you if more than two defined size standard peaks were not matched.	The GeneScan Analysis Software will not alert you to any consecutive peaks at the end of the definition. This is to avoid logging warnings when your sample was not run long enough to include all the defined size standard peaks. This prevents you from having to create a new standard for shorter runs. The Analysis Log will, however, alert you if less than 50% of the defined size standard peaks were not matched, regardless of the peak locations in the definition.

Removing Information from the Analysis Log

To remove information from the Analysis Log:

Step	Action
1	Select the information you want to remove. Note Choose Select All (Ctrl+A) from the Edit menu to select all the information.
2	Select Clear from the Edit menu.


Closing the Analysis Log

- You can either:
- ◆ Click the **Close** button in the upper-left corner, or
 - ◆ Select **Close** (Ctrl+W) from the **File** menu

Verifying Peak Detection

Introduction Use the Show Peak Positions command from the View menu, while the electropherogram and associated tabular data are displayed, to verify results by examining how the GeneScan Analysis Software defined the total area that comprises each peak and the center of the peak.

Verifying Peak Detection To verify peak detection:

Step	Action
1	<p>Select Show Peak Positions from the View menu.</p> <p>Markers appear that identify the beginning, center, and end of each peak.</p> <p>Note For a better view, select Zoom In from the View menu, or use the Zoom tool (see “Zooming In and Out” on page 9-17).</p>  The screenshot shows a window titled 'untitled Display-4' containing an electropherogram. The y-axis represents fluorescence intensity with values from 0 to 531. The x-axis represents time or scan number with a value of 5600. The trace shows several peaks. Two prominent peaks are labeled with 'I' markers at their centers. A smaller peak is labeled with a 'T' marker. The baseline is relatively flat with some noise. At the bottom of the window, there is a status bar showing '1B : Test Sample(cap#1).f5a / someBlueColorInfo' and coordinates 'X: 5393 Y: 285'.
2	<p>Examine the display to ensure that each peak's center, beginning, and end points are correct.</p>

To verify peak detection: *(continued)*

Step	Action
3	<p>Select Hide Peak Positions from the View menu to suppress the display of the peak markers.</p> <p>Note You can also use the Sample Info view to display information about the peaks detected and matched.</p> <p>For more information, see “Description of Information” on page 4-14.</p>

Saving, Archiving, and Copying Files 10

Overview

In This Section Topics in this chapter include the following:

Topics	See Page
Why Save GeneScan Files	10-2
Saving GeneScan Files	10-3
Archiving Sample Files	10-4
Transferring Data to Other Applications	10-5

Why Save GeneScan Files

Reasons for Saving Files The following table explains why you save projects, sample files, and Results Displays. For information on archiving files, see “Archiving Sample Files” on page 10-4.

Save	Because	See
GeneScan® Analysis Software projects	It protects the links to sample files and their preferences. Projects contain links to sample files and preferences regarding display and analysis.	“Saving Projects” on page 10-3.
Sample files	It protects the links to projects and their preferences. Sample files also contain raw data and critical information about the run, settings, and analysis control.	“Saving Sample Files” on page 10-3.
Results Displays	It saves the Results Display settings in projects when you have a display format that suits your needs.	“Saving Results Displays” on page 10-3.

Saving GeneScan Files

Introduction This section describes how to save projects, sample files, and results displays.

Saving Projects **Note** You do not need to save a sample file after analysis. The analyzed data is written directly to the sample file during analysis.

The following table describes the options to save a project:

If you choose...	Then...						
Save Project (Ctrl+S)	you can take the following action:						
	<table><tr><th>If you...</th><th>Then...</th></tr><tr><td>previously saved the project</td><td>it is automatically saved using the same name.</td></tr><tr><td>had not saved the project</td><td>the Save this document as dialog box opens. Select a location for the file, enter a name, and click Save.</td></tr></table>	If you...	Then...	previously saved the project	it is automatically saved using the same name.	had not saved the project	the Save this document as dialog box opens. Select a location for the file, enter a name, and click Save .
	If you...	Then...					
previously saved the project	it is automatically saved using the same name.						
had not saved the project	the Save this document as dialog box opens. Select a location for the file, enter a name, and click Save .						
Save Project As	the Save this document as dialog box opens. Select a location for the file, enter a name, and click Save .						

Saving Sample Files To save sample files, select **Save** (Ctrl+S) from the **File** menu. If you select **Close** from the **File** menu or click the **Close** button when you have not saved the changes, a dialog box opens with a message asking if you want to save them.

Saving Results Displays You can combine electropherograms and tabular data in many ways for display, and the GeneScan Analysis Software allows you to save display combinations and formats for future viewing.

For more information on saving a display for future viewing, refer to “Saving and Renaming the Results Control Format” on page 8-17.

Archiving Sample Files

When to Archive Sample Files Archive sample files when you feel confident that the channel selections (tracking) used to generate them were correct.

Procedure A sample file is 60 KB to 150 KB in size, depending on the length of the run.

Step	Action
1	To archive sample files, drag the file icon or the run folder containing the files to the floppy disk icon or to an alternative storage device.
2	A 1.4 MB high-density disk holds about 12 files.

Transferring Data to Other Applications

Genotyper Software GeneScan Analysis Software files can be read by Genotyper® software for Windows NT® platform.

Cutting and Pasting Tabular Data

To cut and paste tabular data:

Step	Action								
1	Display the tabular data you want to copy.								
2	Select the rows you want to copy by taking the following action: <table><tr><th>If you want to select...</th><th>Then...</th></tr><tr><td>all tabular data</td><td>Choose Select All (Ctrl+A) from the File menu.</td></tr><tr><td>several consecutive rows</td><td>Shift-click the first and last row in the group you want to select.</td></tr><tr><td>several rows that are not listed next to each other</td><td>Ctrl+click the rows.</td></tr></table>	If you want to select...	Then...	all tabular data	Choose Select All (Ctrl+A) from the File menu.	several consecutive rows	Shift-click the first and last row in the group you want to select.	several rows that are not listed next to each other	Ctrl+click the rows.
If you want to select...	Then...								
all tabular data	Choose Select All (Ctrl+A) from the File menu.								
several consecutive rows	Shift-click the first and last row in the group you want to select.								
several rows that are not listed next to each other	Ctrl+click the rows.								
3	Select Copy (Ctrl+C) from the Edit menu.								
4	Optional: To view the contents of the clipboard before pasting, from the Edit menu, select Show Clipboard .								
5	Open the new application and click where you want to place the information.								
6	Select Paste (Ctrl+V) from the Edit menu.								

Creating a Text File

To create a text file from tabular data:

Step	Action
1	Display the tabular data.
2	Select Export Table from the File menu. The Save As dialog box opens.
3	Choose a name and file location in the dialog box and click Save .

Printing Results

Overview

In This Chapter Topics in this chapter include the following:

Topics	See Page
About Printing	11-1
Printing Run Results Automatically	11-2
Printing Selected Sample Files	11-3

About Printing

Ways You Can Print You can print the results of the run automatically at the end of the run or interactively, as selected sample files or display combinations.

If You Get Unexpected Results You might initially get unexpected results from autoprinting or the Print One command if you switch printers.

Step	Action
1	The first time you print after changing the printer configuration, select Print Setup (Ctrl+J) from the File menu and click OK .
2	Select Print (Ctrl+P) from the File menu and use the Print dialog box.

Printing Run Results Automatically

Introduction You can specify that the results are printed in the 310 and 377 Data Collection software or in the GeneScan® Analysis Software.


From the Data Collection Software When you choose automatic printing from the 310 and 377 Data Collection software, the GeneScan Analysis Software prints a separate page for each designated Sample file, showing electropherograms and tabular data as specified in the Auto-Analysis Defaults (see step 2 on page 2-6).

Choose automatic printing in the 310 and 377 Data Collection software as follows:

On this instrument...	Choose...
ABI PRISM® 310	Auto-Print in the Injection list.
ABI PRISM® 377	Auto-Print in the Run Sheet.

For information about setting up your run, refer to the instrument user manual.


Procedure To print the results automatically as the samples are analyzed:

Step	Action
1	In the Analysis Control window, select the check box labeled Print Results . 
2	Select the samples you want to analyze. For more information, see “Analyzing a Sample File” on page 4-26.
3	Click the Print Setup button to specify the samples and the format. For more information, see “Specifying the Format for Printed Results” on page 3-8.
4	Click OK . The results are printed after the results are analyzed.

Printing Selected Sample Files

Introduction Print selected sample files by using the Results Control window or by choosing the sample file.

Setting Printing Options To set the printing options:

Step	Action
1	Select Preferences from the Settings menu and Results Display from the submenu. The Results Display Preferences dialog box opens.
2	Use the Printing Preferences section, change the electropherogram height and page breaks.  Depending on how you set these options, the format that prints may be different from what is on the screen. For information on printing saved Results Control formats, see “Working with a Previously Saved Display” on page 8-18.

Printing from the Results Control Window Print selected samples after analysis, regardless of whether you choose automatic printing.

To print results for selected sample files after analysis:

Step	Action
1	In the Results Control window, select the dye/samples and format you want to print. Use the same technique as you did to select the format and the dye/samples to display the data. For more information, see “Using the Results Control Window” on page 8-7.

To print results for selected sample files after analysis: *(continued)*


Step	Action	
2	You can take the following action:	
	You can either...	Then...
	Click the Print button.	The Print dialog box opens.
	Select Print (Ctrl+P) from the File menu.	Make any changes to the settings and click OK .
	Select Print One from the File menu.	The sample files are printed. Note The Print dialog box does not appear.

Printing from the File Menu

To print a sample file from the File menu:

Step	Action
1	Select Open from the File menu. The Open Existing dialog box opens. Note You can also double-click the sample file name in the folder containing the files. If the GeneScan® Analysis Software is not running, the software starts and opens the sample file.
2	Click the Sample icon. An Open dialog box opens.
3	In the dialog box, navigate to the folder and select the sample file that you want to open.
4	Click Open . The Sample File window opens. For more information about the Sample File window, see page 4-8.
5	Select one of the five views of the Sample File window, and select Print or Print One from the File menu. Note If you select Print One , then the Print dialog box does not appear.

Creating GeneScan Analysis Modules



Overview

In This Chapter Analysis modules provide the auto-analysis feature with the parameters to use for the GeneScan® Analysis Software. For more information about analysis parameters, refer to Chapter 5, “Working with Analysis Parameters.”

This appendix includes the following topics:

Topic	See Page
About GeneScan Analysis Modules	A-2
Creating GeneScan Analysis Modules	A-5

About GeneScan Analysis Modules

Introduction The GeneScan Analysis Software modules contain the following analysis options:

- ◆ Analysis range to use
- ◆ Statistical method used to fit the standard curve to the size standards data

Analysis Module Format GeneScan analysis modules have the file name format *filename.gsp* and are stored in the Params folder.

The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params

Files Referenced by the Module GeneScan analysis modules reference two types of companion files that contain other analysis parameter information.

Default Settings The following table describes the default settings for the analysis parameters:


File type	Description
Sizecaller standard file	<p>Size standards are specific DNA fragments of known sizes. After defining the peaks of a size standard, the GeneScan® Analysis Software matches this definition to the internal size standard included with the run. The software assigns the defined size values to the appropriate peaks of the internal size standard, and uses this information with the selected sizecalling method to size all unknown fragments.</p> <p>Size standard setting is used only during auto-analysis by 3100 and 3700 Data Collection software.</p> <p>For more information, refer to “What Are Size Standards” on page 7-2.</p>

The following table describes the default settings for the analysis parameters: *(continued)*

File type	Description						
Sizecaller parameter file	<p>There are two sizecaller parameter files. The format is <i>filename.scp</i>, and they are stored in the SizeStandards folder at the following directory location:</p> <p>D:\AppliedBio\Shared\Analysis\Sizecaller\Params</p> <table><tr><th>File</th><th>Use this file...</th></tr><tr><td>ABISizecallerAutoAnalysis.scp</td><td>when the run data is being analyzed automatically (the first time it is analyzed).</td></tr><tr><td>ABISizecallerGSAnalysis.scp</td><td>if the data is re-analyzed.</td></tr></table> <p>Note Do not move or delete these files.</p>	File	Use this file...	ABISizecallerAutoAnalysis.scp	when the run data is being analyzed automatically (the first time it is analyzed).	ABISizecallerGSAnalysis.scp	if the data is re-analyzed.
File	Use this file...						
ABISizecallerAutoAnalysis.scp	when the run data is being analyzed automatically (the first time it is analyzed).						
ABISizecallerGSAnalysis.scp	if the data is re-analyzed.						

**Relationship
Between a Module
and the Files**

The relationship between an analysis module and the two companion files is illustrated below.



**Never Modify the
Files**

Never modify the sizecaller standard file or the sizecaller parameter file. If you modify one of these files, it will no longer work.

Creating GeneScan Analysis Modules



Summary of the Procedure

Follow these steps to create an analysis module:



Step	Action
1	Review the size standards data and select the analysis parameters that remove unwanted noise and peaks outside of the size range of the standards. See step 1 in the “Creating a Size Standard File” procedure below.
2	Create a size standard file (<i>filename.szs</i>) for the reviewed standards. See step 2 on page A-6 to step 10 on page A-9.
3	Review the sample file of a sample to be sized, selecting the analysis parameters that optimize the appearance of the data. See step 1 on page A-9 to step 4 on page A-11.
4	Save the analysis parameters as a new .gsp file, referencing the size standard file just created. See step 5 on page A-11 to step 8 on page A-12.

Creating a Size Standard File



To create a size standard file:

Step	Action
Reviewing and Choosing a Size Standard as a Template	
1	Review the size standard data and optimize the analysis parameters.
2	<p>Select New (Ctrl+N) from the File menu.</p> <p>The Create New dialog box opens.</p> 
3	<p>Click the Size Standard icon.</p> <p>The Select Sample File browser box opens.</p> 
4	<p>Navigate to the Completed folder at the following directory location:</p> <p>D:\AppliedBio\DataExtractor\Completed</p>




To create a size standard file: *(continued)*

Step	Action
5	<p>Select the GeneScan Analysis Software sample file, with the extension .fsa, that you want to use as a template.</p> 
6	<p>Click Open.</p> <p>The Select Dye and Analysis Parameters dialog box opens.</p> 
7	<p>From the Dye pop-up menu, select the dye that was used to label the size standard DNA fragments.</p>

To create a size standard file: *(continued)*



Step	Action																								
8	<p>From the Analysis Parameters pop-up menu, select Analysis Parameters.</p>  <p>This references the current analysis parameter setting rather than a specific analysis parameter file.</p>																								
9	<p>Click OK.</p> <p>The following is an example of the dialog box that opens:</p>  <table border="1"> <thead> <tr> <th>Peak</th><th>Data Points</th><th>Size</th></tr> </thead> <tbody> <tr> <td>1</td><td>1957</td><td>0</td></tr> <tr> <td>2</td><td>1981</td><td>0</td></tr> <tr> <td>3</td><td>2000</td><td>0</td></tr> <tr> <td>4</td><td>2010</td><td>0</td></tr> <tr> <td>5</td><td>2029</td><td>0</td></tr> <tr> <td>6</td><td>2036</td><td>0</td></tr> <tr> <td>7</td><td>2050</td><td>0</td></tr> </tbody> </table> <p>Template File Info: File: Eb10d2~3.fsa Run Date: Tue, Feb 02, 1999 Run Time: 10:41:53 AM Dye: R Parameters: <Analysis Parameters></p>	Peak	Data Points	Size	1	1957	0	2	1981	0	3	2000	0	4	2010	0	5	2029	0	6	2036	0	7	2050	0
Peak	Data Points	Size																							
1	1957	0																							
2	1981	0																							
3	2000	0																							
4	2010	0																							
5	2029	0																							
6	2036	0																							
7	2050	0																							

To create a size standard file: *(continued)*



Step	Action
10	<p>In the Size column, enter the known sizes of the standard's peaks.</p> 
Saving the Size Standard	
1	<p>Select Save from the File menu.</p> <p>The Save this document as dialog box opens.</p>
2	<p>Navigate to and open the Size Standards folder (see below).</p> <p>The folder contains size standards (.szs) files.</p> <p>The path is D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards</p> 
3	<p>Enter a name in the File name text box for the size standards files, and click Save.</p> <p>The dialog box closes and the file is saved to the correct location for auto-analysis to read.</p>
4	<p>Click the Close button () in the newly created <i>Filename.szs</i> dialog box.</p> <p>The dialog box closes.</p>

Creating an Analysis Parameter File

To create an analysis parameter file:

Step	Action
1	<p>In the GeneScan Analysis Software, select New from the File menu. The Create New dialog box opens.</p> 
2	<p>Click the Analysis Parameters icon. An untitled Analysis Parameters dialog box opens.</p> <p>Note To display the orange Peak Amplitude Threshold use the scroll bar under the values and scroll to the right.</p> <p>Note The Size Standard option is only specific to the 3100 and 3700 instruments.</p> 


To create an analysis parameter file: *(continued)*

Step	Action
3	Complete the dialog box using the definitions in the “Sizecaller Algorithm Flowchart” on page 5-4.
4	<p>In the AutoAnalysis Only group box, select from the pop-up menu the size standard that you just created.</p> 
5	<p>Select Save from the File menu.</p> <p>The Save this document as dialog box opens.</p>
6	<p>Navigate to and open the Params folder.</p> <p>The path is: D:\AppliedBio\Shared\Analysis\Sizecaller\Params</p> 

To create an analysis parameter file: *(continued)*

Step	Action
7	Enter a file name for the analysis parameter file in the File name text box, and click Save . This saves the file to the correct location for auto-analysis to read.
8	Click the Close button in the newly created <i>Filename.gsp</i> dialog box.

Sizecalling Methods



Overview

In This Appendix Topics in this appendix includes the following:

Topic	See Page
Least Squares Method	B-2
Cubic Spline Interpolation Method	B-4
Local Southern Method	B-5
Global Southern Method	B-7

Least Squares Method

About This Method Both Least Squares Methods (2nd Order and 3rd Order) use regression analysis to build a best-fit sizecalling curve. This curve compensates for any fragments that may run anomalously. As a result, this method normally results in the least amount of error for all the fragments, including the size standards and the samples.

Depending on whether you choose the 2nd or 3rd Order Least Squares Method in the Analysis Parameters dialog box, the resulting size curve is either a quadratic function or a cubic function. The software uses the known standard fragments and the associated scan number positions to produce a sizing curve based on multiple linear regression.

Least Squares Sizecalling Examples The first figure below shows the 2nd Order Least Squares sizecalling curve, and the second figure shows the 3rd Order Least Squares sizecalling curve.




Figure B-1 2nd Order Least Squares sizecalling curve




Figure B-2 3rd Order Least Squares sizecalling curve


Advantages In nearly all instances in the “Least Squares Sizecalling Examples” above, the mobility of an individual DNA fragment is coincident with the best curve fit of the entire data set. Stated differently, the mobility of most DNA fragments is strictly length dependent. This method automatically compensates for fragments that run anomalously.

GeneScan® Analysis Software calculates a best-fit least squares curve for all samples, regardless of the sizecalling method you choose. The curve is black in the Standard Sizing Curve window.

Cubic Spline Interpolation Method

About This Method

By definition, the Cubic Spline Method forces the sizing curve through all the known points of the selected GeneScan size standard. Although this produces exact results for the values of the standards themselves, it does not compensate for standard fragments that may run anomalously.



Possible Local Sizing Inaccuracy

Mobility of any DNA fragment can be affected by its sequence and by secondary and tertiary structure formation. If any internal size standard fragment has anomalous mobility, the Cubic Spline Method may exhibit local sizing inaccuracy.


For example: Assume that a standard fragment is close in molecular length to an unknown sample fragment. Assume further that the standard fragment runs anomalously. The Cubic Spline Method assigns the official value to this standard fragment, even though it may be slightly incorrect. The size of the unknown fragment is then likely to be calculated incorrectly as well.

Note This method does not determine the amount of sizing accuracy error.

Local Southern Method

About This Method

The Local Southern Method determines the sizes of fragments by using the reciprocal relationship between fragment length and mobility, as described by E. M. Southern (1979).



The Equation

The following equation attempts to describe the reciprocal relationship between the mobility, m , and the length, L_0 , of the standard fragments:

$$L = [c/(m - m_0)] + L_0$$

How This Method Works

This method, which is similar to the Cubic Spline Method, uses the four fragments closest in size to the unknown fragment to determine a best-fit line value. Using this method, only the region of the size ladder near the fragment of unknown length is analyzed.


Note Size estimates may be off if any of the standard fragments run anomalously.

The following table summarizes how the Local Southern Method works:

Step	Action
1	The fitting constants of the curve are calculated for each group of three neighboring points on the standard. A separate curve is created for each set of three points.
2	A curve is then created by using three standard points (two points below and one point above the fragment), and a fragment size is determined.
3	Another curve is created by looking at an additional set of three points (one point below and two points above the fragment) and another value is assigned.
4	The two size values are averaged to determine the unknown fragment length.

Global Southern Method

About This Method This method is similar to the Least Squares Method in that it compensates for standard fragments that may run anomalously. The method creates a best-fit line through all the available points, and then uses values found on that line to calculate the fragment values.



The Equations The following table describes how the equations work:

Equation	Description
$L = [c/(m - m_0)] + L_0$	Attempts to describe the reciprocal relationship between the mobility, m, and the length, L ₀ , of the standard fragments.
$\sum i(L_i - (c/(m_i - m_0) + L_0))^2$	The fitting constants L ₀ , m ₀ , and c are calculated by a least squares fit to minimize the following quantity.

How This Method Works

All points in the standard are weighted equally, and the curve is not constrained to go through any specific point. The software can analyze a large range of fragment sizes with this method.

DNA fragments that are...	Are sized using...
not bracketed within the size standard curve	a second-order least squares curve extrapolation.
bracketed within the size standard curve	the method that was chosen.

For best results, use a standard that brackets all the fragments of interest.

GeneScan Size Standards



Overview

About the Size Standards

The GeneScan Analysis Software comes with several ready-to-use size standard definition files that you can choose from to analyze fragments run on the ABI PRISM® 3700 DNA Analyzer. The size standards are stored in the SizeStandards folder.

The path is:

D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards

See also “About Size Standards” on page 7-2.

Size Standards Included

The following table lists the ready-to-use size standards:

Size Standard	See Page
GS 120.szs	C-2
GS 350 All.szs	C-4
GS 400HD.szs	C-9
GS 500 All.szs	C-11

The table below lists size standards in which some of the fragments have been set to 0.

Note You can easily remove any of these sizes by opening the definition and setting any unwanted sizes to 0.


Size Standard	See Page
GS 350 377.szs	C-7
GS 350-250.szs	C-8
GS 500 377.szs	C-13
GS 500-250.szs	C-14

GeneScan 120 Size Standard

About This Size Standard	You can use the GeneScan-120 LIZ™ size standard to determine fragment lengths between 15 and 120 base-pairs.									
Special Uses	This size standard was designed to provide accurate sizing of short DNA fragments. Therefore, it is particularly useful for SNP analysis. All fragments have been checked for migration that is true to size under a wide variety of run conditions. There are no anomalous fragments.									
How It Is Prepared	All aspects of the preparation of the GeneScan 120 LIZ size standard are proprietary. Each fragment contains a single LIZ fluorophore.									
Fragment Lengths	<p>The following table lists the lengths of the nine fragments that comprise the GeneScan 120 LIZ size standard:</p> <table><tr><td>15</td></tr><tr><td>20</td></tr><tr><td>25</td></tr><tr><td>35</td></tr><tr><td>50</td></tr><tr><td>62</td></tr><tr><td>80</td></tr><tr><td>110</td></tr><tr><td>120</td></tr></table>	15	20	25	35	50	62	80	110	120
15										
20										
25										
35										
50										
62										
80										
110										
120										
Denaturing Electropherogram	The GeneScan 120 LIZ size standard is made of single-stranded DNA fragments. The following figure shows the peak patterns of GeneScan 120 fragments run under denaturing conditions. Fragments were run using the 3700 POP-5 polymer at 60 °C.									

Electropherogram of GeneScan 120 LIZ

The following is an electropherogram of GeneScan 120 LIZ:



GeneScan 350 All Size Standard

About This Size Standard	The GeneScan 350 All size standard contains sizes for all fragments in the GS 350 size standard. This size standard is useful for sizing fragments between 35 and 350 base-pairs. The native fragments are uniformly spaced to provide accurate sizecalling.												
How It Is Prepared	The GeneScan 350 All size standard is prepared by Pst 1 digestion of plasmid DNA, followed by ligation of a TAMRA or ROX-labeled 22-mer oligodeoxynucleotide to the cut ends. A subsequent enzymatic digestion with BstU 1 yields DNA fragments containing a single TAMRA or ROX dye (see “GeneScan 350 Molecular Lengths” below).												
GeneScan 350 Molecular Lengths	<p>The following table lists the GeneScan 350 Denatured Fragment Molecular Lengths (Nucleotides):</p> <table><tr><td>35</td><td>160</td></tr><tr><td>50</td><td>200</td></tr><tr><td>75</td><td>250</td></tr><tr><td>100</td><td>300</td></tr><tr><td>139</td><td>340</td></tr><tr><td>150</td><td>350</td></tr></table>	35	160	50	200	75	250	100	300	139	340	150	350
35	160												
50	200												
75	250												
100	300												
139	340												
150	350												


**Running Under
Denaturing
Conditions**

The following table describes running the GeneScan 350 All standard under denaturing conditions:

Like the GeneScan 2500 and GeneScan 1000 standard...	However, like the GeneScan 250 standard...	Consequently
the GeneScan 350 standard is made of double-stranded DNA fragments.	the GeneScan 350 standard has only one labeled strand.	<p>under denaturing conditions, even if the two strands migrate at different rates, only the labeled strand is detected. Refer to “Electropherogram of GeneScan 350” below.</p> <p>Because of this, split peaks are avoided that result when two strands move through a denaturing polymer at different rates.</p>

**Electropherogram
of GeneScan 350**


The following is an electropherogram of GeneScan 350 run under denaturing conditions:



**Double-Stranded
GeneScan 500
Fragments**

The following figure shows the sizes of double-stranded GeneScan 500 fragments. Use these values to size fragments run under native conditions.

IMPORTANT An asterisk (*) for the 250 and 340 base-pair peaks denotes peaks resulting from abnormal migration of double strands that did not completely separate under denaturing conditions when analyzed on the 3100 and 3700 instruments. Do not use these peaks to size samples. The peaks show smaller values than the actual size of the fragments.



GeneScan 350 377 Size Standard

About This Size Standard

The GeneScan 350 377 size standard contains all GS 350 fragment sizes except the 35 and 50 base-pair sizes. The two smallest fragments (35 and 50 base-pairs) are often lost in the primer peak on gel instruments. The construction of this size standard differs from the GS 500 377.szs. This is because its electropherogram begins only after the 50 base-pair peak (depending on the run conditions), so the two smallest sizes are not present and do not need to be set to 0 (as in the GS 500 377 size standard).

This size standard can be used on any instrument.


GeneScan 350 377 Molecular Lengths

The following table lists the GeneScan 350 denatured fragment molecular lengths (nucleotides):

75	200
100	250
139	300
150	340
160	350

Electropherogram of GeneScan 350 377

The following is an electropherogram of GeneScan 350 377 run under denaturing conditions:



GeneScan 350-250 Size Standard

About This Size Standard

The GeneScan 350-250 size standard contains all GS 350 fragment sizes except the 250 base-pair size, which has been set to 0.


GeneScan 350-250 Molecular Lengths

The following table lists the GeneScan Denatured Fragment Molecular Lengths (Nucleotides):

35	160
50	200
75	0
100	300
139	340
150	350

Electropherogram of GeneScan 350-250

The following is an electropherogram of GeneScan 350-250 run under denaturing conditions:



GeneScan-400HD Size Standard

About This Size Standard Use the GeneScan-400HD (High Density) size standard to determine fragment lengths between 50 and 400 base-pairs.

Special Uses The high density of marker bands in this standard makes it particularly useful for microsatellite analysis. All fragments have been checked for migration that is true to size under a wide variety of run conditions on all ABI PRISM® instruments. There are no anomalous fragments (compared with the 250-bp fragment in GeneScan 350 or 500 on the 3700 Analyzer).

GeneScan-400HD is the recommended size standard for use with the ABI PRISM Linkage Mapping Sets.

How It Is Prepared All aspects of the preparation of the GeneScan-400HD size standard are proprietary. Each fragment contains a single ROX fluorophore.


Fragment Lengths The following table lists the lengths of the 21 fragments that make up the GeneScan-400HD size standard:

50	160	260	360
60	180	280	380
90	190	290	400
100	200	300	
120	220	320	
150	240	340	

Denaturing Electropherogram Although the GeneScan-400HD size standard is made of double-stranded DNA fragments, only one of the strands is labeled. Consequently, even if the two strands migrate at different rates under denaturing conditions you will not need to worry about peak splitting. The following figure shows the peak patterns of GeneScan-400HD fragments run under denaturing conditions. Fragments were run using the 3700 POP-6™ polymer at 50 °C.

**Electropherogram
of
GeneScan-400HD**

The following is an electropherogram of GeneScan-400HD:



GeneScan 500 All Size Standard

About This Size Standard

The GeneScan 500 All size standard contains all fragments in the GS 500 size standard. This size standard is useful for sizing fragments between 35 and 500 base-pairs. The native fragments are uniformly spaced to provide accurate base calling.

How It Is Prepared

The GeneScan 500 All size standard is prepared by Pst 1 digestion of plasmid DNA, followed by ligation of a TAMRA or ROX-labeled 22-mer oligodeoxynucleotide to the cut ends. A subsequent enzymatic digestion with BstU 1 yields DNA fragments containing a single TAMRA or ROX dye (see “GeneScan 350 Molecular Lengths” below).

GeneScan 500 Molecular Lengths

The following table lists the GeneScan 500 denatured fragment molecular lengths (nucleotides):

35	250
50	300
75	340
100	350
139	400
150	450
160	490
200	500


Running Under Denaturing Conditions

Like the GeneScan 2500 and GeneScan 1000 standard, the GeneScan 500 standard is made of double-stranded DNA fragments. However, with the GeneScan 500, only one strand of the double-stranded DNA is labeled, whereas the other two standards have labels on both strands. Consequently, under denaturing conditions, even if the two strands migrate at different rates, only the one labeled strand is detected. Because of this, split peaks are avoided that result when two strands move through a denaturing polymer at different rates.

Refer to “Electropherogram of GeneScan 500” below.


Electropherogram of GeneScan 500

The following is an electropherogram of GeneScan 500 run under denaturing conditions:



Double-Stranded GeneScan 500 Fragments

The following figure shows the sizes of double-stranded GeneScan 500 fragments. Use these values to size fragments run under native conditions.



GeneScan 500 377 Size Standard

About This Size Standard

The GeneScan 500 377 size standard includes all GS 500 fragment sizes except the 35 and 50 base-pairs, which have been set to 0. The two smallest fragments (35 and 50 base-pairs) are often lost in the primer peak on gel instruments.

This size standard can be used on any instrument.


GeneScan 500 377 Molecular Lengths

The following table lists the GeneScan 500 377 denatured fragment molecular lengths (nucleotides):

0	250
0	300
75	340
100	350
139	400
150	450
160	490
200	500

Electropherogram of GeneScan 500 377

The following is an electropherogram of GeneScan 500 377 run under denaturing conditions:



GeneScan 500-250 Size Standard

About This Size Standard

The GeneScan 500-250 size standard has the 250 base-pair peak set to 0, since this peak does not migrate as it should on capillary instruments.

Note The 340 base-pair peak is still present in this file and may not migrate properly. You can check this for the run conditions by changing the 340 base-pair definition to 0, and then use the GeneScan Analysis Software to size the size standard. If you are not happy with GeneScan's size for the 340 base-pair fragment, then it is not migrating properly and you should not use it for your definition.


GeneScan 500-250 Molecular Lengths

The following table lists the GeneScan 500-250 denatured fragment molecular lengths (nucleotides):

35	0
50	300
75	340
100	350
139	400
150	450
160	490
200	500

Electropherogram of GeneScan 500-250

The following is an electropherogram of GeneScan 500-250 run under denaturing conditions:



Troubleshooting the GeneScan Software

D

Introduction

In This Appendix The tables in this section present information about problems you might experience with your GeneScan® Analysis Software runs and suggest possible causes and corrections.

Topics in this appendix include the following:

Topics	See Page
Troubleshooting Projects and Results	D-2
Troubleshooting Gel Data	D-4
Troubleshooting Genotyping Results	D-7
GeneScan Analysis Software Error Messages	D-8

Troubleshooting Projects and Results

Table Description The following table describes the problem, probable cause, and correction for troubleshooting projects and results:

Problem	Probable Cause	Correction
File name is dimmed in Size Standard or Parameters column in Analysis Control window.	The file has been moved from the folder of its original location.	<ul style="list-style-type: none"> a. Move the file back to its original location. b. Reset preferences to specify the new folder location. c. Create or select a new file.
Peaks appear on display but the GeneScan Analysis Software does not detect them (cannot select them in electropherogram display).	<ul style="list-style-type: none"> ◆ Peak Amplitude Threshold set too high. ◆ Minimum Peak Half Width set too high. ◆ Electrophoresis run too quickly resulting in poor resolution. 	<ul style="list-style-type: none"> a. Adjust minimum peak height to include smallest peaks desired and re-analyze. b. Reduce minimum peak half-width setting and re-analyze. c. Repeat electrophoresis at reduced power. <p>For more information, refer to “Peak Detection Options” on page 5-8.</p>
At the position of one strong peak additional colors appear underneath the peak.	Off-scale data not multicomponented correctly.	<ul style="list-style-type: none"> a. Repeat electrophoresis; load less sample. b. Regenerate sample files.
Peaks appearing in a dye color that should not be present.	Bleed-through from other colors because of off-scale data.	Repeat electrophoresis; load less sample.

The following table describes the problem, probable cause, and correction for troubleshooting projects and results: *(continued)*

Problem	Probable Cause	Correction
Peak centers seem to be incorrect in electropherogram.	<ul style="list-style-type: none"> ◆ Resolution of the gel might be inadequate (ABI PRISM® 377). ◆ Signal-to-noise ratio might be too low. 	Repeat electrophoresis at lower power.
Software cannot display the sizing curve for a sample.	<ul style="list-style-type: none"> ◆ Sample's in-lane size standard does not match defined size standard. ◆ Sample file was not sized called. 	<ul style="list-style-type: none"> ◆ Re-analyze the sample file with a different size standard or create a new one. <p>In the Analysis Control window:</p> <ul style="list-style-type: none"> a. Ctrl+click dye/sample that represents Size Standard for the sample file. b. Select a size standard definition file and re-analyze.
Peaks disappear in the electropherogram.	Included the primer peak in the analysis.	Re-analyze sample file without primer peak.

Troubleshooting Gel Data

Table Description

The following table lists the problem, probable cause, and correction for troubleshooting of gel data:

Problem	Probable Cause	Correction
At the position of one strong peak, additional colors appear underneath the peak.	<ul style="list-style-type: none"> ◆ Off-scale data not multicomponented correctly. ◆ Poor / incorrect matrix. ◆ Gel Image not multicomponented. 	<ul style="list-style-type: none"> a. Repeat electrophoresis; load less sample. b. Attach a new gel matrix, regenerate the gel image. c. Regenerate the gel image with multicomponenting selected.
Peaks appearing in a dye color that should not be present.	Bleed-through from other colors because of off-scale data.	Repeat electrophoresis; load less sample.
TAMRA-labeled size standard appears yellow on the gel display.	Collected using: <ul style="list-style-type: none"> ◆ Filter set A (ABI 373), or ◆ Virtual filter A (ABI PRISM 377). 	Repeat electrophoresis with: <ul style="list-style-type: none"> ◆ Filter set B (ABI 373), or ◆ Virtual filter C (ABI PRISM 377).
TET-labeled products not seen on gel display.	Collected using: <ul style="list-style-type: none"> ◆ Filter set A (ABI 373), or ◆ Virtual filter A (ABI PRISM 377). 	Repeat electrophoresis with: <ul style="list-style-type: none"> ◆ Filter set B (ABI 373), or ◆ Virtual filter C (ABI PRISM 377).

The following table lists the problem, probable cause, and correction for troubleshooting of gel data: *(continued)*

Problem	Probable Cause	Correction
Signal showing up in neighboring lanes.	Leaking wells of gel.	<ul style="list-style-type: none"> ◆ Consider using a square-tooth comb instead of a shark-tooth comb. ◆ If using 96 lanes, then rerun gel using protocol in the <i>ABI PRISM 377 DNA Sequencer 96-Lane Upgrade User's Manual</i> (P/N 4305423).
	Signal intensity very high and signal is being detected in neighboring lanes due to closeness of spacing.	Move tracker lane position from center of band to edge of the band away from strong signal and extract as usual.
		Use 1- or 2-lane averaging to extract lanes.
HEX-labeled products appear green on gel display.	Collected using: <ul style="list-style-type: none"> ◆ Filter set A (ABI 373), or ◆ Virtual filter A (ABI PRISM 377). 	Repeat electrophoresis with: <ul style="list-style-type: none"> ◆ Filter set B (ABI 373), or ◆ Virtual filter C (ABI PRISM 377).
Collection time was sufficient, but only a small portion of gel displayed.	<ul style="list-style-type: none"> ◆ Gel Image Processing preferences did not include enough scans to display entire gel. ◆ Electrophoresis power too low. 	a. Regenerate gel image with new scan range. b. Adjust to correct settings; repeat electrophoresis.

The following table lists the problem, probable cause, and correction for troubleshooting of gel data: *(continued)*

Problem	Probable Cause	Correction
Improper tracking results.	Bad matrix.	Attach new matrix.
	Sample Sheet not filled out properly.	Fill out Sample Sheet properly.
	Comb types set improperly.	a. Fix and type in gel preferences. b. Retrack gel.
	Peak height or red signal too low.	Rerun gel with more size standard.

Troubleshooting Genotyping Results

Table Description The following table describes the problem, probable cause, and correction for troubleshooting the Genotyper® software results:

Problem	Probable Cause	Correction
Allele peaks seen in correct molecular weight range, with additional peaks seen outside this range.	<ul style="list-style-type: none"> ◆ Bleed-through from other colors because of off-scale data. ◆ Primers not fully optimized. 	<ul style="list-style-type: none"> a. Repeat electrophoresis; load/inject less sample. b. Check optimization.
With allele peaks of high intensity, the GeneScan Analysis Software calls many small peaks.	<ul style="list-style-type: none"> ◆ Background above minimum peak height. ◆ Too much PCR product loaded. 	<ul style="list-style-type: none"> a. Adjust minimum peak height threshold; re-analyze. b. Repeat electrophoresis; load/inject less sample. <p>For more information, refer to “Peak Detection Options” on page 5-8.</p>
A homozygous individual shows a dip at the top of an allele peak which may be called as two separate peaks.	Truncated single peak because of off-scale data can appear as two peaks.	Repeat electrophoresis; load/inject less sample.
Warning message: “Could not complete ‘Run Macro’ command because the labeled peak could not be found.”	The first allele peak for one or more loci in the allelic ladder is lower than the preset minimum peak height specification in the categories list.	<ul style="list-style-type: none"> a. Adjust minimum peak height; re-analyze. b. Repeat electrophoresis; load/inject less sample. <p>For more information, refer to “Peak Detection Options” on page 5-8.</p>

GeneScan Analysis Software Error Messages

Introduction This section includes two tables:

- ◆ Analysis Log Error Messages
- ◆ Error Messages When Defining Size Standards

Analysis Log Error Messages The following table describes the error messages you might encounter in the GeneScan Analysis Software Log:

Analysis Log Error Message	Comment/Correction	Refer To
The Analysis Range parameter does not include enough data points. Check your analysis parameters.	Make sure the Analysis Range in your analysis parameters contains at least 250 data points.	“Setting Analysis Parameters” on page 5-5.
The Range of Data Points parameter to analyze is too large. Check your analysis parameters.	Specify a smaller range in the analysis parameters.	“Setting Analysis Parameters” on page 5-5.
The analysis Parameters could not be accessed. Check your Analysis Parameters Setting.	Make sure the Analysis Parameters file specified in the Analysis Control window is valid and accessible.	“Setting Analysis Parameters” on page 5-5.

Error Messages When Defining Size Standards

The following table describes the error messages you might encounter while defining size standards:

Error Message	Comment/Correction	Refer To
The affected sample file is not available. Locate the sample file and try again.	If the sample file name is dim in the Analysis Control window, the GeneScan Analysis Software has not located the sample file. You can instruct the program to search for the sample file.	"Finding Missing Sample Files" on page 2-14.
A Dye Standard is not selected for the affected sample file. Select a Dye Standard and try again.	Select the dye/sample that represents the standard by Ctrl+clicking the appropriate dye/sample field.	
The affected sample file does not have a valid Analysis Parameters Selection. Select new analysis parameters and try again.	Select either: ♦ The default program parameters (<Analysis Parameter >), or ♦ A valid analysis parameters file in the Analysis Control window.	"Using Analysis Parameter Files" on page 5-13.
No peaks were found within the Analysis Range. Check your analysis parameters.	♦ Make sure the Peak Amplitude Threshold setting allows for detection of the peaks in your sample. ♦ If peaks in your data are narrow, make sure the Minimum Peak Half Width is a small number.	"Sizecaller Algorithm Flowchart" on page 5-4.

GeneScan Analysis Software Files



Table of Files The following table lists the files that the GeneScan Analysis Software reads, writes, and, in most cases, creates. The software does not create gel files, and creates Sample files only through lane extraction from a gel.

Table E-1 GeneScan Analysis Software files

File Type	Created by	Location	Where is it used
.fsa	3700 Data Extractor 3100 Data Extractor Gel Processor (377) 310 Data Collection	In a folder often with the title "Run Folder"	GeneScan software can be used to display, view and edit any .fsa file.
.gel	377 Data Collection	In a folder often with the title "Run Folder"	The gel file contains the raw data collected during the instrument run. The gel file is tracked and data is extracted into sample files by Gel Processor.
filename.gsp	Shipped with GeneScan software. Users can create custom files.	D:\AppliedBio\Shared\Analysis\Sizecaller\Params	These files specify certain ranges and methods used during data analysis. Users can create custom .gsp files.
filename.szs	Shipped with GeneScan software. Users can create custom files.	D:\AppliedBio\Shared\Analysis\Sizecaller\SizeStandards	These files are used to identify peak sizes for specified size standards run under certain conditions. Users can define the standards after running them on the instruments.

Table E-1 GeneScan Analysis Software files *(continued)*

File Type	Created by	Location	Where is it used
filename.prj	GeneScan Software	User can specify location. This file is often stored in a folder with the sample files in of the corresponding project	Projects contain references to Sample files. Sample files of a single project can be from one or multiple runs. Projects allow a group of data to be organized, displayed and analyzed together.
filename.mxt	GeneScan Software, using the Matrix making feature		Matrix files must be created for each instrument. The matrix standards are collected on 310 and 377 instruments. These files are used to create .mxt files using GeneScan Analysis software.
Sample log.log	GeneScan Software	D:\AppliedBio\GeneScan\Bin	Contains a running record of analysis performed by the software
Analysis log.log	GeneScan Software	D:\AppliedBio\GeneScan\Bin	Contains a running record of analysis performed by the software

Technical Support

F

Contacting Technical Support

You can contact Applied Biosystems for technical support by telephone or fax, by e-mail, or through the Internet. You can order Applied Biosystems user documents, MSDSs, certificates of analysis, and other related documents 24 hours a day. In addition, you can download documents in PDF format from the Applied Biosystems Web site (please see the section "To Obtain Documents on Demand" following the telephone information below).

To Contact Technical Support by E-mail

Contact technical support by e-mail for help in the following product areas:

Product Area	E-mail address
Genetic Analysis (DNA Sequencing)	galab@appliedbiosystems.com
Sequence Detection Systems and PCR	pcrlab@appliedbiosystems.com
Protein Sequencing, Peptide and DNA Synthesis	corelab@appliedbiosystems.com
Biochromatography, PerSeptive DNA, PNA and Peptide Synthesis systems, CytoFluor [®] , FMat [™] , Voyager [™] , and Mariner [™] Mass Spectrometers	tsupport@appliedbiosystems.com
Applied Biosystems/MDS Sciex	api3-support@sciex.com
Chemiluminescence (Tropix)	tropix@appliedbiosystems.com

**Hours for
Telephone
Technical Support**

In the United States and Canada, technical support is available at the following times:

Product	Hours
Chemiluminescence	8:30 a.m. to 5:30 p.m. Eastern Time
Framingham support	8:00 a.m. to 6:00 p.m. Eastern Time
All Other Products	5:30 a.m. to 5:00 p.m. Pacific Time

**To Contact
Technical Support
by Telephone or
Fax**

In North America

To contact Applied Biosystems Technical Support, use the telephone or fax numbers given below. (To open a service call for other support needs, or in case of an emergency, dial **1-800-831-6844** and press **1.**)

Product or Product Area	Telephone Dial...	Fax Dial...
ABI PRISM® 3700 DNA Analyzer	1-800-831-6844 , then press 8	1-650-638-5981
DNA Synthesis	1-800-831-6844 , then press 21	1-650-638-5981
Fluorescent DNA Sequencing	1-800-831-6844 , then press 22	1-650-638-5981
Fluorescent Fragment Analysis (includes GeneScan® applications)	1-800-831-6844 , then press 23	1-650-638-5981
Integrated Thermal Cyclers (ABI PRISM® 877 and Catalyst 800 instruments)	1-800-831-6844 , then press 24	1-650-638-5981
ABI PRISM® 3100 Genetic Analyzer	1-800-831-6844 , then press 26	1-650-638-5981
Bioinformatics (includes BioLIMS®, BioMerge®, and SQL GT™ applications)	1-800-831-6844 , then press 25	1-505-982-7690
Peptide Synthesis (433 and 43X Systems)	1-800-831-6844 , then press 31	1-650-638-5981
Protein Sequencing (Procise® Protein Sequencing Systems)	1-800-831-6844 , then press 32	1-650-638-5981

Product or Product Area	Telephone Dial...	Fax Dial...
PCR and Sequence Detection	1-800-762-4001 , then press 1 for PCR, 2 for the 7700 or 5700, 6 for the 6700 or dial 1-800-831-6844 , then press 5	1-240-453-4613
Voyager™ MALDI-TOF Biospectrometry and Mariner™ ESI-TOF Mass Spectrometry Workstations	1-800-899-5858 , then press 13	1-508-383-7855
Biochromatography (BioCAD® Workstations and Poros® Perfusion Chromatography Products)	1-800-899-5858 , then press 14	1-508-383-7855
Expedite™ Nucleic acid Synthesis Systems	1-800-899-5858 , then press 15	1-508-383-7855
Peptide Synthesis (Pioneer™ and 9050 Plus Peptide Synthesizers)	1-800-899-5858 , then press 15	1-508-383-7855
PNA Custom and Synthesis	1-800-899-5858 , then press 15	1-508-383-7855
FMAT™ 8100 HTS System and CytoFluor® 4000 Fluorescence Plate Reader	1-800-899-5858 , then press 16	1-508-383-7855
Chemiluminescence (Tropix)	1-800-542-2369 (U.S. only), or 1-781-271-0045	1-781-275-8581
Applied Biosystems/MDS Sciex	1-800-952-4716	1-650-638-6223

Outside North America

Region	Telephone Dial...	Fax Dial...
Africa and the Middle East		
Africa (English Speaking) and West Asia (Fairlands, South Africa)	27 11 478 0411	27 11 478 0349
South Africa (Johannesburg)	27 11 478 0411	27 11 478 0349
Middle Eastern Countries and North Africa (Monza, Italia)	39 (0)39 8389 481	39 (0)39 8389 493
Eastern Asia, China, Oceania		
Australia (Scoresby, Victoria)	61 3 9730 8600	61 3 9730 8799
China (Beijing)	86 10 64106608	86 10 64106617
Hong Kong	852 2756 6928	852 2756 6968
Korea (Seoul)	82 2 593 6470/6471	82 2 593 6472
Malaysia (Petaling Jaya)	60 3 758 8268	60 3 754 9043
Singapore	65 896 2168	65 896 2147
Taiwan (Taipei Hsien)	886 2 2358 2838	886 2 2358 2839
Thailand (Bangkok)	66 2 719 6405	66 2 319 9788
Europe		
Austria (Wien)	43 (0)1 867 35 75 0	43 (0)1 867 35 75 11
Belgium	32 (0)2 712 5555	32 (0)2 712 5516
Czech Republic and Slovakia (Praha)	420 2 61 222 164	420 2 61 222 168
Denmark (Naerum)	45 45 58 60 00	45 45 58 60 01
Finland (Espoo)	358 (0)9 251 24 250	358 (0)9 251 24 243
France (Paris)	33 (0)1 69 59 85 85	33 (0)1 69 59 85 00
Germany (Weiterstadt)	49 (0) 6150 101 0	49 (0) 6150 101 101
Hungary (Budapest)	36 (0)1 270 8398	36 (0)1 270 8288
Italy (Milano)	39 (0)39 83891	39 (0)39 838 9492
Norway (Oslo)	47 23 12 06 05	47 23 12 05 75
Poland, Lithuania, Latvia, and Estonia (Warszawa)	48 (22) 866 40 10	48 (22) 866 40 20

Region	Telephone Dial...	Fax Dial...
Portugal (Lisboa)	351 (0)22 605 33 14	351 (0)22 605 33 15
Russia (Moskva)	7 095 935 8888	7 095 564 8787
South East Europe (Zagreb, Croatia)	385 1 34 91 927	385 1 34 91 840
Spain (Tres Cantos)	34 (0)91 806 1210	34 (0)91 806 1206
Sweden (Stockholm)	46 (0)8 619 4400	46 (0)8 619 4401
Switzerland (Rotkreuz)	41 (0)41 799 7777	41 (0)41 790 0676
The Netherlands (Nieuwerkerk a/d IJssel)	31 (0)180 331400	31 (0)180 331409
United Kingdom (Warrington, Cheshire)	44 (0)1925 825650	44 (0)1925 282502
All other countries not listed (Warrington, UK)	44 (0)1925 282481	44 (0)1925 282509
Japan		
Japan (Hacchobori, Chuo-Ku, Tokyo)	81 3 5566 6100	81 3 5566 6501
Latin America		
Del.A. Obregon, Mexico	305-670-4350	305-670-4349

**To Reach
Technical Support
Through the
Internet**

We strongly encourage you to visit our Web site for answers to frequently asked questions and for more information about our products. You can also order technical documents or an index of available documents and have them faxed or e-mailed to you through our site. The Applied Biosystems Web site address is

<http://www.appliedbiosystems.com/techsupp>

To submit technical questions from North America or Europe:

Step	Action
1	Access the Applied Biosystems Technical Support Web site.
2	Under the Troubleshooting heading, click Support Request Forms , then select the relevant support region for the product area of interest.
3	Enter the requested information and your question in the displayed form, then click Ask Us RIGHT NOW (blue button with yellow text).
4	Enter the required information in the next form (if you have not already done so), then click Ask Us RIGHT NOW . You will receive an e-mail reply to your question from one of our technical experts within 24 to 48 hours.

**To Obtain
Documents on
Demand**

Free, 24-hour access to Applied Biosystems technical documents, including MSDSs, is available by fax or e-mail or by download from our Web site.

To order documents...	Then...
by index number	a. Access the Applied Biosystems Technical Support Web site at http://www.appliedbiosystems.com/techsupp b. Click the Index link for the document type you want, then find the document you want and record the index number. c. Use the index number when requesting documents following the procedures below.

To order documents...	Then...
by phone for fax delivery	<p>a. From the U.S. or Canada, call 1-800-487-6809, or from outside the U.S. and Canada, call 1-858-712-0317.</p> <p>b. Follow the voice instructions to order the documents you want.</p> <p>Note There is a limit of five documents per request.</p>
through the Internet for fax or e-mail delivery	<p>a. Access the Applied Biosystems Technical Support Web site at http://www.appliedbiosystems.com/techsupp</p> <p>b. Under Resource Libraries, click the type of document you want.</p> <p>c. Enter or select the requested information in the displayed form, then click Search.</p> <p>d. In the displayed search results, select a check box for the method of delivery for each document that matches your criteria, then click Deliver Selected Documents Now (or click the PDF icon for the document to download it immediately).</p> <p>e. Fill in the information form (if you have not previously done so), then click Deliver Selected Documents Now to submit your order.</p> <p>Note There is a limit of five documents per request for fax deliver but no limit on the number of documents you can order for e-mail delivery.</p>

License and Warranty



Applera Corporation Software License and Limited Product Warranty

Warranty Statement	CAREFULLY READ THE FOLLOWING TERMS AND CONDITIONS, WHICH APPLY TO THE SOFTWARE ENCLOSED (THE "SOFTWARE"). YOUR OPENING OF THIS PACKAGE INDICATES YOUR ACCEPTANCE OF THEM. IF YOU DO NOT ACCEPT THEM, PROMPTLY RETURN THE COMPLETE PACKAGE, AND YOUR MONEY WILL BE RETURNED. THE LAW PROVIDES FOR CIVIL AND CRIMINAL PENALTIES FOR ANYONE IN VIOLATION OF THE LAWS OF COPYRIGHT.
-------------------------------	---

Copyright	The SOFTWARE, including its structure, organization, code, user interface, and associated documentation, is a proprietary product of Applied Biosystems group of Applera Corporation, and is protected by international laws of copyright. Title to the SOFTWARE, and to any portion of the SOFTWARE, shall at all times remain with Applera Corporation and its subsidiaries.
------------------	--

License	<p>1. You may use the SOFTWARE on a single computer (or in a single network, if your software is designated as a network version). You may transfer the SOFTWARE to another single computer (or network, if a network version), but the SOFTWARE may never be installed on more than one computer (or more than one network, if a network version) at any one time.</p> <p>2. Registration codes for multiple users are available as follows: Access to the number of copies of SOFTWARE is provided by a registration code. For example, if you have a registration code that enables you to three copies of SOFTWARE simultaneously, you cannot install the SOFTWARE in more than three separate computers.</p>
----------------	---

3. You may make one copy of the SOFTWARE for backup purposes.
4. You may not transfer this SOFTWARE license to a third party.

Restrictions

1. You may not copy, transfer, rent, modify, use, or merge the SOFTWARE, or the associated documentation, in whole or in part, except as expressly permitted in this Agreement.
2. You may not reverse assemble, decompile, or otherwise reverse engineer the SOFTWARE.

Limited Warranty

For a period of 90 days after purchase of the SOFTWARE, Applied Biosystems warrants that the SOFTWARE will function substantially as described in the documentation delivered with the SOFTWARE. If you believe that you have discovered an error which causes substantial deviation from that documentation, and after written notification from the purchaser, and Applied Biosystems is able to reliably reproduce that error at Applied Biosystems' facility, then Applied Biosystems will have the option to either (i) correct the error in a subsequent release of the SOFTWARE, which shall be supplied to you free of charge, or (ii) accept a return of the SOFTWARE from the purchaser, and refund your purchase price of the SOFTWARE. Applied Biosystems does not warrant that the SOFTWARE will meet your requirements, or will be error-free, or will conform exactly to the documentation.

Limitation of Liability

EXCEPT AS SPECIFICALLY STATED IN THIS AGREEMENT, THE SOFTWARE IS PROVIDED AND LICENSED "AS IS". THIS WARRANTY IS GIVEN IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. APPLIED BIOSYSTEM'S SOLE LIABILITY FOR BREACH OF WARRANTY SHALL BE LIMITED TO A REFUND OF THE PURCHASE PRICE FOR SUCH PRODUCT. IN NO EVENT WILL APPLIED BIOSYSTEMS BE LIABLE FOR ANY DAMAGES, INCLUDING INCIDENTAL OR CONSEQUENTIAL DAMAGES, EVEN IF APPLIED BIOSYSTEMS HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Term

You may terminate this Agreement by destroying all copies of the SOFTWARE and documentation. Applied Biosystems may terminate this Agreement if you fail to comply with all of its terms, in which case

you agree to return to Applied Biosystems all copies of the SOFTWARE and associated documentation.

Miscellaneous

1. Failure to enforce any of the terms and conditions of this Agreement by either party shall not be deemed a waiver of any rights and privileges under this Agreement.

2. In case any one or more of the provisions of this Agreement for any reason shall be held to be invalid, illegal, or unenforceable in any respect, such invalidity, illegality, or unenforceability shall not affect any other provisions of this Agreement, and this Agreement shall be construed as if such invalid, illegal, or unenforceable provisions had never been contained herein.

3. This Agreement shall be construed and governed by the laws of the State of California.

4. This Agreement constitute the entire agreement between Applied Biosystems and you concerning the SOFTWARE.

5. This Software License Agreement supersedes any other license agreement which may be contained in this software package.

Glossary

This glossary defines special terminology used in the *GeneScan Analysis Software User's Manual*. The terms are listed in alphabetical order. Many terms are defined in the text of the manual, so if you do not find a term here, check the index to see if you can locate it in the manual.

align-by-size curve	<p>Curve created by the GeneScan® Analysis Software for aligning data by size.</p> <p>The software calculates a best-fit, least-squares curve for all samples. This is a third-order curve when you use the Third Order Least Squares sizecalling method; for all other sizecalling methods it is a second-order curve.</p> <p>This curve is black in the Standard Sizing Curve window, although when the sizing curve and this curve match, they overlap so you see only the sizing curve.</p> <p>See also sizecalling curve and size standard spline interpolation curve.</p>
Align By Size/Data Point	<p>Shows the horizontal scale of the electropherograms by fragment size or by data point.</p> <p>You can display data by size only if you ran an internal size standard with your samples and sized called the data.</p>
analysis parameters	<p>Options that specify certain ranges and methods used during analysis using the GeneScan Analysis Software.</p> <p>The software has default analysis parameters that are stored in the project itself.</p> <p>These parameters apply globally, unless you create your own parameters files for use with specific protocols.</p>
baselining	<p>Adjusting the baselines of detected dye colors to the same level for a better comparison of relative signal intensity.</p>

Sequence Collector	<p>Sequence Collector (formerly BioLIMS™) is a genetic information management system that provides a relational database for storage and retrieval of DNA sequence and fragment data.</p> <p>In addition to the database itself, Sequence Collector contains a set of software applications and tools for maintaining and interacting with the database.</p> <p>The Sequence Collector database resides on a UNIX workstation and uses a Sybase® or an Oracle® database server.</p> <p>The client applications run on Windows NT®-based computers and/or on UNIX workstations.</p>
data point	<p>The 3700 Data Collection program samples data as it passes by the detector.</p> <p>Each “sampling” is stored as a data point.</p>
dye color indicator	<p>Left color box in the Results Control window and the legend of the Results Display.</p> <p>In the Results Display, click this box to move the associated electropherogram to the front.</p> <p>In the Results Control and the Results Display windows, double-click this box to change the dye scale, or Ctrl+double-click it to reset the dye scale to the default.</p>
dye/sample	<p>Individual sample labeled with a single dye within a sample file.</p> <p>Sample files normally contain up to four dye/samples, depending on how many labeled samples you included in each lane or injection of your 3700 Data Collection program run.</p>
electropherogram	<p>Four-color picture of a sequence, showing peaks that represent the bases. The term is used interchangeably with chromatogram in this manual.</p>
grid	<p>Spreadsheet-like display used for entering data in tabular format.</p> <p>The Analysis Control and Results Control windows display grids for entering sample information.</p>

internal size standard	<p>Also called internal lane or injection size standard, DNA fragment of known size that you include with your run. On the ABI PRISM 3700 DNA Analyzer you include the size standard with each injection.</p> <p>Running an internal lane standard results in particularly accurate and precise molecular length determination because the internal lane standard and the unknown fragments undergo exactly the same electrophoretic forces.</p> <p>The software can then compensate for band-shift artifacts caused by variations in the run.</p>
legend	<p>Informational text that appears beneath electropherogram panels in the results displays.</p> <p>You can show or hide legends, and use the color boxes displayed in them to bring specified electropherograms to the front of the panel, or customize the colors.</p>
matrix file/multicomponent matrix	<p>File used to adjust for the spectral overlap between the fluorescent dyes used on the ABI PRISM® instruments. A mathematical matrix of the spectral overlaps is created and the inverse matrix is used to correct the data during analysis. Matrix files are stored in the AppliedBio folder inside the Macintosh® computer System folder, or in the matrix folder in the GeneScan Analysis Software folder. The values of the matrix are stored in the gel file (ABI 373 and ABI PRISM 377) and in the Sample files.</p> <p>For more information, see Chapter 6, "Making a Matrix File."</p>
overlaid	<p>Displayed together so they overlap.</p> <p>In the GeneScan Analysis Software Results Display window, all electropherograms in a single panel are overlaid.</p> <p>You can bring a specific one to the front by clicking the color box that represents it in the legend.</p>
plot color indicator	<p>Right color box in the Results Control window and the legend of the Results Display.</p> <p>In the Results Display window, click this box to move the associated electropherogram to the front.</p> <p>In the Results Control and the Results Display windows, double-click this box to change the plot color, or Ctrl+double-click it to reset the plot color to the default.</p>

preferences	<p>Defaults you can set so that certain parameters are automatically applied when you are working with a project.</p> <p>The GeneScan Analysis Software remembers preferences and applies them globally to all new projects.</p>
project	<p>File containing links to a set of sample files that you want to analyze and display together.</p> <p>A project can contain sample files from multiple runs. Adding a sample file to a project creates a reference to the file. It does not copy the file into the project.</p>
project options	<p>Formatting information you can set for the current project. Project options are remembered by the project when you open it again.</p>
sample files	<p>Computer files that contain raw and analyzed data.</p> <p>Sample files are created directly by the 3700 DNA Analyzer. Sample files contain data such as peak locations, sizecalling values, and a record of analysis settings.</p>
sizecalling curve	<p>Curve created by the GeneScan Analysis Software for sizecalling.</p> <p>The software calculates this curve based on the sizecalling method you specify for data analysis.</p> <p>This curve is read in the Standard Sizing Curve window. When it matches the align-by-size curve, the two overlap so you see only this curve.</p> <p>See also “align-by-size curve” and “size standard spline interpolation curve.”</p>
size standard	<p>Specific DNA fragments of known sizes.</p> <p>After you define the peaks of a size standard, the GeneScan Analysis Software matches this definition to the internal lane or injection standard that you include with your run.</p> <p>The software assigns the defined size values to the appropriate peaks of the internal lane or injection standard, and uses this information with the selected sizecalling method to size all unknown fragments.</p>

size standard spline
interpolation curve

Curve created by the GeneScan Analysis Software for aligning data by size.

The software creates this curve if you use the Local Southern or Cubic Spline Interpolation sizecalling method and the size standard data does not match the best-fit curve, which is normally used for aligning the data by size.

This curve is blue in the Sizing Curve window. See also align-by-size curve and sizecalling curve.

tiled

Displayed so they do not overlap.

The GeneScan Analysis Software displays tiled electropherogram panels in the Results Display.

If you display more than one electropherogram in each panel, all electropherograms in the panel are overlaid.

Index

A

ABI 373

- generating sample files
procedure 6–15 to 6–16

ABI Prism 310

- automatic analysis diagram 2–2, 2–3

ABI PRISM 310

- loading and running dye standards 6–9 to 6–11

ABI PRISM 377

- generating sample files 6–15 to 6–16
- loading and running dye standards 6–12 to 6–14

Analysis Control window

- using to analyze project files 3–2 to 3–5
- using to analyze sample file 3–6 to 3–15

Analysis Log, using 9–36 to 9–37

analysis modules

- about A–2 to A–4
- creating
 - analysis parameter file A–10 to A–12
 - size standard file A–6 to A–9

analysis parameter file A–10 to A–12

analysis parameters

- about 5–2
- setting 5–5 to 5–12
 - Analysis Range options 5–7
 - Auto Analysis Only option 5–12
 - Baselining options 5–12
 - Data Processing options 5–7
 - default settings 5–5, A–2
 - displaying analysis parameters 5–6
 - Peak Detection options 5–8, 5–9 to 5–10
 - Sizecall Range options 5–10
 - Sizecalling Method options 5–10 to 5–11
- using 5–13 to 5–17
 - changing existing parameters 5–17
 - creating custom parameters 5–16
 - deleting custom parameters 5–17
 - different parameters 5–15
 - displaying default parameters 5–16

analysis parameters (continued)

- using (continued)
 - the same parameters 5–13 to 5–14

analysis results

- changing how displayed and printed 8–13 to 8–15

Results Display window

- about 8–4 to 8–6
- using 8–7 to 8–12

steps to evaluating

- evaluating results 8–2
- updating the results 8–16
- using the Analysis Log 9–36 to 9–37
- ways to display 8–3
- See Also* Sample Results View

analyzing

- project files 3–2 to 3–5
- sample files 3–6 to 3–15

AppliedBio folder, files installed in 7–6

archiving 10–4

Auto Analysis, setting analysis parameters 5–12

automatic analysis, setting up 2–4 to 2–6

B

Baseline Window Size, setting analysis parameters 5–12

C

CD-ROM drive 1–10

colors

- assigning to electropherograms 9–23
- defining in electropherograms 9–24 to 9–27

copying data to other applications 10–5

Cubic Spline Method B–4

customer support. *See* technical support F–1

D

Data Collection software

- automatically printing results from 11–2
- deleting, custom analysis parameters 5–17
- disk space recommended 1–10

Documents on Demand F-6
dye color, changing in electropherograms 9-14
dye scale, changing 9-28 to 9-30

E

electropherograms
 about 9-2 to 9-3
 changing dye scale 9-28 to 9-30
 defining custom colors 9-24 to 9-27
 displaying electropherogram data 9-9 to 9-11
 displaying electropherograms and tabular data 9-4 to 9-8
 using to verify peak detection 9-38
 using to verify results 9-31 to 9-32
 verifying size calculations 9-33 to 9-35
 working with electropherogram data 9-12 to 9-23
 assigning standard colors 9-23
 changing horizontal scale 9-21
 changing the dye color 9-14
 changing vertical scale 9-22 to 9-23
 displaying peak positions 9-14
 displaying x-and y-axis positions 9-13
 highlighting peaks 9-15
 making active window 9-13
 scrolling the display 9-16
 showing data by fragment size 9-20
 about 9-20
 showing off-scale data 9-18
 using legends to change the display 9-15
 zooming in and out 9-17
e-mail, address for technical support F-1
EPT Data View 4-24 to 4-25
 displaying the view 4-24
 example 4-25
 what it displays 4-24
error messages D-8 to D-9

F

files, software files, table of E-1 to E-2
folders
 AppliedBio folder, files installed in 7-6
 defining locations 3-16

G

gel files
 troubleshooting gel data D-4 to D-5
GeneScan
 about the size standards C-1
 analysis modules
 about A-2 to A-4
 creating A-6 to A-12
 creating analysis parameter file A-10 to A-12
 See Also analysis parameters
 archiving 10-4
 copying, data to other applications 10-5
 registering 1-8
 related manuals 1-7
 saving
 projects 10-3
 Results Displays 10-3
 sample files 10-3
 why save 10-2
 size standard
 GeneScan 350 All C-4 to C-7
 GeneScan 400HD C-9
 GeneScan 500 All C-11 to C-12
 GeneScan 500 to 250 C-14
 software files, table of E-1 to E-2
GeneScan 500 All size standard C-11 to C-12
GeneScan-400 HD size standard C-9 to C-10
Genotyper software, troubleshooting software
 results D-7
Global Southern Method B-7 to B-8

H

hard drive partitions 1-11
hardware and software requirements to run
 program 1-10
help. *See* technical support F-2
horizontal scale, changing in
 electropherograms 9-21

I

Internet address
 Documents on Demand F-6

L

Least Squares sizecalling method B-2 to B-3

legends, using to change display in
 electropherograms 9–15
Local Southern Method B–5 to B–6

M

Macintosh, converting sample files to the
 NT 4–3
macintosh, converting sample files to the
 Windows NT platform 4–4
manual
 related manuals 1–7
matrix files
 about 6–2 to 6–6
 assigning to sample files 6–23 to 6–24
 bad matrix files, causes 6–26 to 6–27
 choosing a scan range 6–17 to 6–19
 defined Glossary–3
 evaluating the matrix file 6–25
 generating files (ABI 373 and ABI PRISM
 377) 6–15 to 6–16
 installing new file
 sample file 4–27
 loading and running dye standards
 ABI PRISM 310 6–9 to 6–11
 ABI PRISM 377 6–12 to 6–14
 new matrix file, generating 6–20 to 6–21
 process of creating new file 6–7 to 6–8
 saving and naming 6–22
monitor, recommended 1–10
multicomponent matrix defined Glossary–3

O

off-scale data
 showing in electropherograms 9–18
Operating System requirement 1–10

P

partitions, hard drive 1–11
PCR Single-stranded Conformation
 Polymorphism (SSCP)
 Minimum Peak Half Width setting 5–8
peak detection
 setting analysis parameters 5–8, 5–9 to
 5–10
 verifying 9–38

peak positions, displaying in
 electropherograms 9–14
peaks, highlighting in electropherograms 9–15
printing
 about printing 11–1
 automatically printing run results 11–2
 printer recommended 1–10
 sample files 11–3 to 11–4
projects
 analyzing project files 3–2 to 3–5
 creating
 finding missing sample files 2–14 to
 2–15
 process diagrams, using the ABI Prism
 310 2–2, 2–3
 using to manage sample files 2–6 to
 2–7
 working with 2–8 to 2–13
 creating a new project 2–9 to 2–13
 defining folder locations 3–16
 importance of saving 10–2
 opening an existing project 2–8
 removing samples from a project 2–13
 saving, procedure 10–3
 troubleshooting D–2 to D–3

R

Raw Data View 4–22 to 4–23
 displaying the view 4–22
 example 4–23
 what it displays 4–22
 what to evaluate 4–23
registering the software 1–8
Registration Number 1–8
requirements, hardware and software 1–10
Results Control window
 saving and renaming 8–17 to 8–19
 important considerations 8–17
 renaming current display 8–19
 saving the format 8–17
 using previously saved formats 8–18
 working with previously saved
 displays 8–18
Results Display window
 about 8–4 to 8–6
 using 8–7 to 8–12

Results Displays

- importance of saving 10-2
- saving, procedure 10-3

S

Sample File window

- about 4-8
- EPT Data View 4-24 to 4-25
- Raw Data View 4-22 to 4-23
- Sample Info View 4-11 to 4-19
- Sample Results View 4-9 to 4-11
- Size Curve View 4-20 to 4-21

sample files

- about 4-2 to 4-3
- analyzing 4-26 to 4-27
 - installing new matrix file 4-27
 - procedure 4-26
 - using Analysis Control window 3-6 to 3-15
- archiving 10-4
- converting Macintosh files to the NT 4-3
- converting Macintosh files to the Windows NT platform 4-4
- file size 1-10
- finding missing sample files 2-14 to 2-15
- matrix file, assigning to sample files 6-23 to 6-24
- opening, procedure 4-7
- printing 11-3 to 11-4
- removing from a project 2-13

Sample File window

- about 4-8
- EPT Data view 4-24 to 4-25
- Raw Data View 4-22 to 4-23
- Sample Info View 4-11 to 4-19
- Sample Results View 4-9 to 4-11
- Size Curve View 4-20 to 4-21
- saving, procedure 10-3
- unlocking sample files 2-12
- using project to manage files 2-6 to 2-7
- why save 10-2

Sample Info View 4-11 to 4-19

- description of view 4-11 to 4-19

Sample Results View 4-9 to 4-11

saving

- matrix file 6-22
- projects 10-3

saving (continued)

- Results Control format 8-17 to 8-19
- Results Displays 10-3
- sample files 10-3
- why save 10-2

Show Clipboard command 10-5

Size Curve View 4-20 to 4-21

size standard file A-6 to A-9

size standards

- about size standards 7-2
- about the size standards C-1
- defining, how to 7-3 to 7-8
- file, creating for analysis module A-6 to A-9

GeneScan 350 377 C-7

GeneScan 350 All C-4 to C-7

GeneScan 400HD C-9

GeneScan 500 377 C-13

GeneScan 500 All C-11 to C-12

GeneScan 500 to 250 C-14

using 7-9 to 7-14

Sizecall Range, setting analysis

parameters 5-10

sizecalling methods

Cubic Spine Method B-4

Global Southern Method B-7 to B-8

Least Squares Method B-2 to B-3

Local Southern Method B-5 to B-6

setting analysis parameters 5-10 to 5-11

software

files, table of E-1, E-1 to E-2

software, registering 1-8

SSCP, Minimum Peak Half Width setting 5-8

T

tabular data displays

- about 9-2 to 9-3
- displaying electropherogram and tabular data 9-4 to 9-8
- using to verify peak detection 9-38
- using to verify results 9-31 to 9-32

technical support F-1 to F-7

e-mail address F-1

Internet address F-6

regional sales offices F-4 to F-5

telephone/fax (North America) F-2

troubleshooting

- bad matrix files, causes 6–26 to 6–27
- error messages D–8 to D–9
- gel data D–4 to D–5
- Genotyping software results D–7
- projects and results D–2 to D–3

U

- unlocking sample files 2–12

V

- vertical scale, changing 9–22 to 9–23
- viewing scale, changing for
 - electropherograms 9–17

W

- WWW address
 - Applied Biosystems F–6
 - Documents on Demand F–6

X

- x-axis positions, displaying in
 - electropherograms 9–13

Y

- y-axis positions, displaying in
 - electropherograms 9–13

Z

- zooming, electropherograms 9–17



850 Lincoln Centre Drive
Foster City, California 94404-1128
USA

Headquarters

850 Lincoln Centre Drive
Foster City, CA 94404 USA
Phone: +1 650.638.5800
Toll Free (In North America): +1 800.345.5224
Fax: +1 650.638.5884

Worldwide Sales and Support

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches into 150 countries on six continents. For sales office locations and technical support, please call our local office or refer to our web site at www.appliedbiosystems.com or to the Technical Support and Training appendix in this document.

www.appliedbiosystems.com



Applera Corporation is committed to providing the world's leading technology and information for life scientists. Applera Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Printed in the USA, 01/2001
Part Number 4308923 Rev. B

an **Applera** business